## **INgezim® ASFV CROM Ag 2.0** R.11.AS2.K.42

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# INgezim® ASFV CROM Ag 2.0

R.11.AS2.K.42

INgezim® ASFV CROM Ag 2.0 is a double antibody lateral flow assay based on the use of a monoclonal antibody (MAb) specific of VP72 of African Swine Fever Virus (ASFV).

## KIT FEATURES


### APPLICATION

Detection of African Swine Fever Virus (antigen) in swine blood samples.

### TECHNICAL BASIS

This diagnostic device consists of a plastic cassette with three windows that protects an immunochromatographic strip located inside. On the nitrocellulose membrane, an anti-protein control antibody (control line) and an anti-VP72 antibody (test line) are fixed at the level of the C/T window.

If a sample is positive, the antigen present in the sample will be recognized by the recombinant antibody bonded to the black latex beads. As a result, a black coloration appears in the region of the test line. The latex/antibody/antigen complexes migrate across the membrane and, upon reaching the region of the test line, the virus is again recognized by the antibody immobilized in that area. The presence of the control line indicates that the test is valid and that immunochromatography has been performed correctly.

## ASSAY VALIDATION


### ANALYTICAL SENSITIVITY

Analytical sensitivity was evaluated by serial dilution of the VP72 protein in a negative matrix. The obtained results indicated that the assay is able to detect the sample of 0.16µg/ml, which is equivalent to 3.5 ng of protein.

### DIAGNOSTIC SENSITIVITY


#### Experimental samples

A total of 125 experimental bloods characterized as positive by PCR were tested at the European reference laboratory for ASF, CISA-INIA, CSIC. The results were evaluated in terms of days post-infection and compared with the results obtained with the previous version of the commercial assay INgezim® ASFV CROM Ag. Both assays showed comparable sensitivity to the PCR technique. It is 87% for samples collected at days 7-10 p.i. and 47% at days 11-20 p.i.

#### Field Samples

106 positive wild boar samples and 126 positive pig samples collected during surveillance campaigns in different reference laboratories were analyzed, classified by the PCR technique certified in each Laboratory. The results are shown in the graph. The difference in percentages is since in pigs a symptomatological control is carried out and the samples were collected from live animals. However, in wild boar the samples were taken from carcasses and, therefore, they are damaged samples in which the quality of the antigen can be compromised.

### DIAGNOSTIC SPECIFICITY

165 negative bloods collected from ASF-free Spanish farms, 150 bloods from American pigs collected by the USDA (USA) that were characterised as negative by PCR and problematic with the commercial assay INgezim® ASFV CROM Ag (previous version) and 193 wild boar bloods collected during surveillance campaigns in different reference laboratories and classified as negative were analyzed. The general specificity of the assay for pigs was 99.4%, improving the parameters of the previous version of the commercial assay. The specificity in wild boar was 97.4%.

### KIT COMPOSITION

- Immunochromatographic devices
- Dropper with chromatography buffer
- 20 µL disposable pipette

Registration No. 10671-RD  
EXPIRY DATE: 24 MONTHS  
Stored at 4°C-25°C

**GOLD STANDARD DIAGNOSTICS MADRID, S.A.**  
C/Hermanos García Noblejas 39, 8º  
28037 MADRID (SPAIN)  
Tel: (+34)91 3680501  
<http://www.goldstandarddiagnostics.com/>