INgezim® ASFV-R

R.11.ASF.K1

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INgezim® ASFV-R

R.11.ASF.K1

INgezim® ASFV-R is an immunoenzymatic assay based on indirect ELISA technique, which uses a monoclonal antibody (Mab) specific to porcine IgG and recombinant proteins cp312 and p30 of African Swine Fever Virus (ASFV).


KIT FEATURES

APPLICATION

Detection of specific antibodies to cp312 and p30 proteins of ASFV in serum, blood (fresh or on paper) and spleen exudate samples, from swine and wild boar.

TECHNICAL BASE

1. Plates are supplied coated with antigen (ASFV recombinant protein). Samples are added to the wells and incubated.
2. If the samples contain specific antibodies to cp312 and p30 proteins, they will bind to the antigen.
3. When the conjugate (monoclonal antibody specific to porcine IgG, conjugated with peroxidase, AcM-PO) is added, it will bind to the antibodies bound to the antigen. This binding is revealed by colorimetric reaction after addition of the substrate.

RESULTS INTERPRETATION

The assay establishes two cut offs, which will classify the samples as **Positive** or **Negative**, depending on the value of the optical density of the sample in the assay, considering a range of optical densities close to the cut offs as **Doubtful** results.

ASSAY VALIDATION

ANALYTICAL SENSITIVITY

- **EU Reference Serum (CISA-INIA).** 3 batches of the weakly positive EU reference serum for ASFV (ASF-CL27, CL-28 and CL-34) were analysed. The obtained results indicated that the assay is able of detecting this serum as positive.
- **Detection of specific antibodies of different genotypes.** A study with sera from animals infected with different isolates of genotypes (I, II and X) was carried out. The assay can detect antibodies against different ASFV genotypes.
- **Precocity of the test.** INgezim® ASFV-R detects specific antibodies from day 7 post-infection in 15% of cases, between days 15 and 21 p.i. in 50% of the cases and from the 22nd day p.i. in 100% of the cases.

ANALYTICAL SPECIFICITY


There is not cross-reaction with antibodies specific to other agents.

DIAGNOSTIC SPECIFICITY

622 swine and wild boar sera were analysed, all of them from ASFV free areas. The obtained results indicated 100% specificity with samples from free areas.

DIAGNOSTIC SENSITIVITY

49 samples from experimentally intranasal inoculated animals with BeninDMGF strain, 39 sera from swine experimentally infected with different isolates from I, II and X genotypes, and 18 field samples from ASFV positive swine were analysed. The obtained results indicated 94.3% sensitivity.

CORRESPONDENCE WITH IMMUNOBLOTTING (IB)

86 field samples were analysed, which had been previously classified by IB. The obtained results indicated 100% sensitivity and specificity in swine samples and 100% sensitivity and 92% specificity in wild boar samples.

KIT COMPOSITION

- 96 well microtitration plates
- Vials with Positive Control
- Vials with Negative Control
- Vials with Conjugate
- Bottles with Washing Solution
- Bottles with Diluent
- Bottles with substrate (ABTS)
- Bottles with Stop Solution

Spanish registration nº 11128RD
EXPIRATION: 18 MONTHS
Stored at 2°C-8°C

GOLD STANDARD DIAGNOSTICS MADRID, S.A.
C/Hermanos García Noblejas 39, 8º
28037 MADRID (SPAIN)
Tel: (+34)91 3680501
<http://www.goldstandarddiagnostics.com/>