INgezim® EHDV Compac

R.12.EHD.K3

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INgezim® EHDV Compac

R.12.EHD.K3

INgezim® EHDV Compac is an immunoenzymatic assay based on the blocking ELISA technique, which uses a recombinant VP7 protein of the Epizootic Hemorrhagic Disease (EHD) in combination with an antibody specific to the same protein.


KIT FEATURES

APPLICATION

Detection of antibodies specific to the Epizootic Hemorrhagic Disease Virus (EHDV) protein in serum, plasma or blood from cows, goats, sheep or deer.

TECHNICAL BASIS

1. The plates are supplied coated with the VP7 recombinant protein of EHDV. Samples are added to the wells and incubated.
2. If the samples contain antibodies specific to the VP7 recombinant protein, they will bind to the antigen.
3. When the conjugate (Monoclonal Antibody specific of VP7 recombinant protein, peroxidase-labeled, mAb-PO) is added, it will bind to the protein if there are no specific antibodies in the sample (negative). If these antibodies were present (positive sample), they would block the coated antigen preventing the binding of anti-VP7 antibodies.
4. This binding is revealed by colorimetric reaction after the addition of substrate.

RESULTS INTERPRETATION

The assay establishes two cut-offs, which will classify samples as **positive** or **negative**, depending on the value of the optical density of the sample in the test.

ASSAY VALIDATION

ANALYTICAL SENSITIVITY

Serial factor 2 dilutions of a positive serum and a negative serum were analyzed, and the detection limit was compared to the detection limit of the commercial test approved for this use: ID Screen® EHDV Competition. Both tests were able to detect the serum as positive up to a dilution of 1/32.

ANALYTICAL SPECIFICITY

To determine whether there is a cross-reaction with antibodies developed against other agents, animal sera positive for antibodies against bluetongue virus and tuberculosis and negative for EHDV were analyzed by date of collection (samples from Spain before 2021). The obtained results indicated that there is no cross-reaction with the specific antibodies of these agents.

DIAGNOSTIC SENSITIVITY

254 cow sera collected from farms with compatible symptoms, 9 goat/sheep sera collected during routine sampling on farms, 94 sera collected in 2023 during wild animal surveillance campaigns, at the University of Córdoba, classified positive by the commercial product: ID Screen® EHDV Competition based on the same protein were analyzed. The obtained results indicated a relative sensitivity or PPA (Positivity Percentage Agreement) compared to the ID Screen® EHDV competition of 100%.

DIAGNOSTIC SPECIFICITY

751 cow sera (collected before the virus entered Spain), 671 goat and sheep sera (collected before the virus entered Spain or at the same time as the positive samples and characterised as negative by ID Screen® EHDV Competition), as well as 276 deer sera (collected before the virus entered Spain or at the same time as the positive samples and characterised as negative by ID Screen® EHDV Competition), were analyzed. The obtained results indicated a relative specificity or NPA (Negativity Percentage Agreement) of 99.8%.

USE OF WHOLE BLOOD AND PLASMA AS A MATRIX

To determine if blood and plasma matrices could be used in the test, paired blood, plasma and serum samples from positive animals were analyzed. The results showed that the test can also be used in whole blood and plasma matrix.

KIT COMPOSITION

- Microtitration plates
- Vials with Positive Control
- Vials with Negative Control
- Vials with Conjugate
- Bottles with Washing Solution
- Vials with Diluent
- Bottles with Substrate (TMB)
- Bottles with Stop Solution

Spanish registration n°:
EXPIRATION: 24 MONTHS
Stored at 2°C-8°C

GOLD STANDARD DIAGNOSTICS MADRID, S.A.
C/Hermanos García Noblejas 39, 8º
28037 MADRID (SPAIN)
Tel.: (+34)91 3680501
<http://www.goldstandarddiagnostics.com/>