INgezim® Leishmania 2.0

R.15.LS2.K1

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INgezim® Leishmania 2.0

R.15.LS2.K1

INgezim® Leishmania 2.0 is based on an indirect enzymatic immunoassay (indirect ELISA), using a recombinant protein of a *L. infantum* kinesin in combination with an A/G protein. It is a DIVA assay for differentiation of infected or vaccinated animals.

KIT FEATURES

APPLICATION

Detection of specific antibodies to *Leishmania* in canine serum or blood. This is a DIVA assay for the differentiation of infected or vaccinated animals. The assay can be performed either manually or with automated equipment.

TECHNICAL BASIS OF THE KIT

1. Plates are supplied coated with recombinant protein of *L. infantum*. The samples are added to the wells and incubated.
2. If the samples contain *Leishmania* specific antibodies, these will bind to the antigen.
3. When the specific Peroxidase substrate is added, a colorimetric reaction will be detected in the presence of the conjugate.

INTERPRETATION OF THE RESULTS


The assay establishes two cut offs: positive and negative.

Positive samples: OD > the positive cut off.

Negative Samples: OD is ≤ the negative cut off.

Doubtful Samples: OD is between the two cut offs.

Titration: The titer of the sample will be the last dilution showing and OD higher than the positive cut off.

ASSAY VALIDATION

DIAGNOSTIC SENSITIVITY

152 canine sera classified as positive by INgezim® LEISHMACROM were analyzed. The animals used were 69 sera from Spain, 15 sera from Brazil and 68 sera from Germany, of which 50 were analyzed by Indirect Immunofluorescence Assay (IFA). The obtained results indicated a relative sensitivity or PPA with respect to INgezim® LEISHMACROM of 93% (141/152).

The diagnostic sensitivity with respect to IFA was 85% in this study using the 50 IFA-positive samples.

DIAGNOSTIC SPECIFICITY

A total of 276 canine sera from classified as negative by the commercial product INgezim® LEISHMACROM were analyzed. The animals used were 212 negative sera from Spain, 21 sera from Brazil, 43 sera from Germany, of which 24 were analyzed by IFA.

The obtained results indicated a relative specificity or NPA with respect to INgezim® LEISHMACROM of 99.6%.

The diagnostic specificity for IFA was 95.8% in this study using the 24 IFA-negative samples.

ANALYTICAL SPECIFICITY

Sera from canine Babesiosis and Leishmaniosis negative animals were analyzed by IFA. The obtained results indicated that there is no cross-reaction with specific antibodies to this agent.

At the same time, other studies were carried out with animals infected by other agents, in which it was demonstrated that the kit does not show cross-reaction with canine *Babesiosis*, *Giardia duodenalis*, *Dirofilaria repens*, *Toxocara canis*, *Ehrlichia canis*, *Ancylostoma caninum* or *Anaplasmosis*.

DOSE-EFFECT COHERENCE

Serial factor 2 dilutions of a positive serum and a negative serum were analyzed.


The obtained values indicated that the assay responds correctly to the dilution of a sample by varying proportionally the obtained OD 450nm.

ANALYTICAL SENSITIVITY

The same serum was analyzed as for the determination of dose-effect coherence previously characterized by IFA with titer 1/3200. The assay is able to detect as positive the 1/6400 dilution of this serum.

USE OF AUTOMATIC EQUIPMENT

The assay can be performed either manually or on automatic equipment. In order to confirm that the test works correctly, a validation was carried out with automatic equipment.

COMPOSITION OF THE KIT

- 96 wells microtitration plates
- Vials with Positive Control
- Vials with Negative Control
- Vials with conjugate
- Bottle with washing solution
- Bottle with diluent
- Bottle with substrate (TMB)
- Bottle with stop solution

Registration number in Spain:
EXPIRATION: 24 MONTHS
Stored at 2°C-8°C

GOLD STANDARD DIAGNOSTICS MADRID, S.A.
C/Hermanos García Noblejas 39, 8º
28037 MADRID (SPAIN)
Tel: (+34)91 3680501
<http://www.goldstandarddiagnostics.com/>