
INgezim® EHDV Compac

Prod Ref: 12.EHD.K.3

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KIT COMPOSITION

REAGENT	2 Plates		5 Plates	
	Uni.	Vol.	Uni.	Vol.
Microtitration plates	2		5	
Vials Positive Control	1	1 mL	2	1 mL
Vials Negative Control	1	1 mL	2	1 mL
Vials with Peroxidasa Conjugate 100x concentrated	1	350 µL	2	350 µL
Bottles with Diluent (DE01-01), ready to use	1	60 mL	1	125mL
Bottles with Washing Solution 25x concentrated	1	125mL	1	125mL
Bottles with Substrate (TMB), ready to use	1	30 mL	1	60 mL
Bottles with Stop Solution, ready to use	1	60 mL	1	60 mL

OTHER MATERIALS AND REAGENTS NEEDED NOT PROVIDED WITH THE KIT

- Distilled or deionized water.
- Micropipettes.
- Disposable micropipette tips.
- Plate washing device.
- Test tubes.
- ELISA Reader (450 nm filter).

TECHNICAL BASIS

This kit is based on the blocking enzyme-linked immunosorbent assay (Blocking ELISA) technique, which is briefly described below:

The wells of the plate are coated with the target antigen, VP7 of the Epizootic hemorrhagic disease virus (EHDV). The serum to be tested is added into the wells and, after a washing step, a peroxidase-labeled monoclonal antibody specific for the VP7 protein is added to the wells. If the test serum contains antibodies to EHDV, these will bind to the protein immobilized on the plate, preventing binding of the monoclonal antibody. If the test serum does not contain such antibodies, the peroxidase-labeled monoclonal antibody will bind to the antigen immobilized on the plate. After removal of any unbound material from the plate by sequential washing steps, the presence of bound monoclonal antibody is detected by the addition of the chromogenic substrate, which turns into a measurable blue color in the presence of peroxidase. The appearance of color in the well indicates the absence of virus-specific antibodies in the tested sample, while the absence of color indicates the presence of antibodies in the sample.

PRECAUTIONS AND WARNINGS FOR USERS

1. Read the instructions for use carefully.
2. Bring all reagents to room temperature prior to use.
3. Do not mix instructions or reagents from different kits and avoid any contamination of the reagents.
4. Do not use components after expiration dates and do not mix components from different lots.
5. Do not eat, drink, or smoke while specimens or kit reagents are being handled and do not pipette by mouth.
6. Use a new tip for each sample to be tested.
7. The substrate is extremely sensitive to light and contamination, so it is recommended to remove only the volume to be used from the bottle and to never return the excess substrate to the bottle.
8. Stop solution is a strong acid solution. In case of accidental contact with skin, wash gently with water.
- 9.

STORAGE OF COMPONENTS

Keep all the reagents between +2°C and +8°C.

INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipette suitable for dispensing 300 µL in each well. After the incubation periods, the washing steps must be done following the instructions:

- Throw out the content of the plate by tightly overturning the plate to avoid the possible mixture of different wells content. As a precaution, the wells should be emptied over a cuvette containing 0.1M NaOH solution, as the test samples may contain infectious agents.
- Dispense a volume of 300 µL of washing solution in each well.
- Softly mix the plate, avoiding the exchange of liquid between wells.
- Turn over the plate abruptly to empty the wells.
- Repeat the process as many times as indicated on the instructions of the kit. After the washing step, verify that the next reagent to be added to the plate is ready to use. Do not keep the plate dry.
- After the last washing step, shake the plate upside down on an absorbent filter paper to eliminate all washing buffer.

PREPARATION OF SAMPLES

Serum, blood, or plasma samples:

Serum, blood, or plasma samples collected from cattle, goat, sheep, or deer can be used in the kit (do not use highly hemolyzed or contaminated samples, they can affect assay's result). Samples must be assayed at a 1/2 dilution in the supplied diluent (DE01-01). This dilution could be done directly in the plate by adding 50 µL of diluent and 50 µL of sample into each well.

Dried blood on filter paper:

Add 2-3 drops of blood onto a filter paper, for example TFN specimen collection card (Ref. 8.460.0004.A-N, Ahlstrom Munksjö) and allow to dry overnight at room temperature (20°C - 25°C). Samples collected by this protocol are stable at 4°C for longer periods of time than fresh blood samples.

Add 250 µL/well of DE01-01 diluent to a low or medium binding capacity microtiter plate. Take two pieces of the filter paper using a punch and add them into the same well to be wetted with the diluent. Incubate overnight (16-24 h) at 37°C ± 1°C.

Add 100 µL per well of the extracted sample to perform the assay. DO NOT DILUTE.

The use of dried blood samples on filter paper may show slightly lower sensitivity than the use of the other matrices, especially in weakly positive samples.

PREPARATION OF REAGENTS

Washing solution:

Dilute one part of the concentrate washing solution provided in the kit with 24 parts of distilled or deionized water (i.e., 40 mL of concentrate plus 960 mL of dH₂O). When ready, this solution remains stable between +2°C and +8°C until the expiration date described on the label of the concentrated solution.

Control sera:

The control sera must be treated as test samples, diluted 1/2 in diluent prior to be used. This dilution can be done directly in the plate by adding 50 µL of diluent plus 50 µL of control serum.

Conjugate. To be prepared just prior to use:

Prepare a 1/100 dilution in the diluent included in the kit (DE01-01):

- For one strip, 8 wells, it is recommended to dilute 10 µL of the 100x-concentrated conjugate in diluent (DE01-01) to a final volume of 1 mL.
- For a complete plate, it is recommended to dilute 110 µL of the 100x-concentrated conjugate in diluent (DE01-01) to a final volume of 11 mL.

Shake gently before using.

It is recommended to only prepare the needed volume per each assay, since the leftover reagent must be discarded.

TEST PROCEDURE

1. Before performing the assay, let all kit components reach room temperature (20°C to 25°C).
2. If dilution is directly prepared in the plate, add first 50 µL of the supplied diluent to each well. Then, add 50 µL of control sera and of the samples to be tested. Softly shake the plate for an appropriate homogenization of the samples, careful to avoid liquid transfer between wells. It is recommended to include control sera in duplicate for each assay. Seal the plate and incubate for 1 h ± 5 min at 25 ± 1°C.
3. Empty the wells into a container with NaOH 0.1M solution and wash 3 times as previously described.
4. Add 100 µL of specific conjugate to each well. Seal the plate and incubate for 30 ± 5 min at 25 ± 1°C.
5. Wash 3 times as previously described.
6. Add 100 µL of substrate to each well. It is recommended to use multichannel pipettes for increased speed and homogeneity. Incubate the plate for 10 ± 1 min at room temperature, (20 - 25°C).
7. Add 100 µL of stop solution to each well following the same order used for substrate dispensation.
8. Read the Absorbance at 450 nm (Abs450) of each well using an ELISA reader within 5 min after the addition of the stop solution.

READING AND RESULT INTERPRETATION

VALIDATION OF THE TEST:

The test results should be considered valid only if the Abs450 of the Negative Control is, at least, 5 times higher than the Abs450 of the Positive Control, and if the Abs450 of the negative is higher to 1.

$$1) \frac{Abs_{Negative\ control}}{Abs_{positive\ control}} \geq 5$$

$$2) Abs_{Negative\ control} \geq 1$$

RESULTS INTERPRETATION:

When samples are tested in duplicate, Abs450 of the sample should be calculated as the arithmetic mean of values in both wells.

To interpretate assay's results, calculate the following percentage:

$$Competition \% = \left(\frac{Abs_{Negative\ Control} - Abs_{Sample}}{Abs_{Negative\ Control} - Abs_{Positive\ Control}} \right) \times 100$$

*Where Abs means Absorbance at 450 nm.

For deer:

- Samples with a S/N higher or equal to 40 % would be considered **positive**.
- Samples with a S/N lower than 40 % would be considered **negative**.
-

For cattle, goat, or sheep:

- Samples with a S/N higher or equal to 60 % would be considered **positive**.
- Samples with a S/N lower than 60 % would be considered **negative**.

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