INgezim® LEISHMANIA 2.0

Prod. Ref: 15.LS2.K1

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KIT COMPOSITION

<i>REAGENT</i>	<i>1 Plates</i>		<i>10 Plates</i>	
	Uni.	Vol.	Uni.	Vol.
Microtitration plates	1		10	
Vials containing Negative Control Sera, ready to use	1	3,5 mL	1	10 mL
Vials containing Positive Control Sera, ready to use	1	3,5 mL	1	10 mL
Vials containing concentrated peroxidase conjugate	1	1,5 mL	1	13 mL
Bottles with Washing Solution 10x concentrated	1	100 mL	8	250 mL
Bottles with Diluent (DE03-01), ready to use	1	125 mL	5	250 mL
Bottles with Substrate (TMB), ready to use	1	15 mL	3	60 mL
Bottles with Stop Solution, ready to use	1	15 mL	2	125 mL

OTHER MATERIALS AND REAGENTS NEEDED NOT PROVIDED WITH THE KIT

Distilled or deionized water.
Micropipettes from 5 to 200 µL.
Disposable micropipette tips.
Washing plates device.
Test tubes from 50 to 250 mL.
ELISA Reader (450 nm filter)

TECHNICAL BASIS

This kit is based on an indirect enzymatic immunoassay (Indirect ELISA). The technique is briefly described below: The recombinant antigen, specific of Leishmania antibodies, is fixed on polystyrene plates. When a sample contains Leishmania-specific antibodies, they will bind to the recombinant antigen adsorbed on the plate. After a washing step to eliminate all non-fixed material, the presence of dog antibodies can be detected using a specific peroxidase conjugate. After the addition of the substrate a colorimetric reaction will appear which can be measured by a spectrophotometer at 450nm. Please note that our kit uses recombinant antigen obtained by genetic engineering. This method warranties the total absence of infectivity in any of the components of the kit.

The assay, due to the use of an antigen that is not present in any of the vaccines currently on the market, can be considered a DIVA assay able of differentiating between vaccinated and infected animals. Under the current situation, there is no risk of vaccines interfering with the results. However, this situation could change in the future, so it is recommended to consider this possible scenario.

PRECAUTIONS AND WARNINGS FOR USERS

- Read the instructions for use carefully.
- Bring all reagents to room temperature at least during 2 hours prior to use. VERY IMPORTANT!
- Do not mix reagents or use instructions from different kits.
- Avoid any contamination of the reagents of the kit.
- Do not use components after expiry dates and do not mix components from different batches.
- There should be no eating, drinking or smoking when handling reagents or samples.
- Do not pipette by mouth.
- Use a new tip for each serum sample.
- Systematically include a positive control and a negative control each time the assay is run.
- Stop solution is a strong acid. Handle with care. In case of eye or skin contact, wash immediately with plenty of water.
- Substrate must be handled with care, as it is very sensitive to light and contamination. Remove only the volume to be used, with sterile tip or by decantation, and never return the leftover to the bottle.

STORAGE OF COMPONENTS

Store all plates and reagents between +2°C and +8°C, to keep them stable until expiry date.

INFORMATION ABOUT THE WASHING STEPS

The washing steps can be done using an automatic plate washer or a multichannel pipetting device suitable for dispensing 300 μ l on each well.

After the incubation periods, the washing steps must be done following the instructions below:

- Throw away the content of the plate by briskly turning the plate over in order to avoid the possible exchange of the contents from one well to another.
- Dispense a volume of 300 μ L of washing solution on each well.
- Shake the plate gently, avoiding cross-contamination between wells.
- Briskly turn the plate over to empty the wells.
- Repeat the process as many times as indicated in the kit's instructions.
- Prior to emptying the content of the plate in the last washing step, verify that the next reagent to be added is ready to use. Do not let the plate dry longer than strictly necessary.
- After the last washing step, tap the plate upside down on absorbent filter paper to remove any remaining washing solution.

PREPARATION OF SAMPLES

Sera samples must be tested at 1/100 dilution in diluent DE03-01 (i.e.: 5 μ L of serum + 495 μ L of serum diluent DE03-01). For titration we recommend assaying two-fold dilution from the 1/100.

PREPARATION OF REAGENTS

- Washing solution: Dilute one part of the concentrated washing solution provided in the kit with 9 parts of distilled or deionized water (100 mL of concentrated solution and 900 mL of water). Once prepared this solution remains stable when stored between +2°C and +8°C.
 - Control sera (+) and (-): Both controls are ready to use. Do not dilute.
 - Serum diluent: This reagent is ready to use. Do not dilute.
 - Conjugate: This reagent is concentrated. Dilute in supplied diluent depending on the followed procedure. Homogenize well the solution prior to use. Prepare only the needed quantity for one use, as leftover must be discarded.
- o Manual procedure: dilute 1/10:
- The necessary and sufficient quantity for a plate is 1100 μ L of conjugate in 9,900 mL of diluent.
 - The necessary and sufficient quantity for a 8 well strip is 100 μ L of conjugate in 0,900 mL of diluent.
- o Automated procedure: dilute 1/70:
- The necessary and sufficient quantity for a plate is 158 μ L of conjugate in 10,842 mL of diluent.
 - The necessary and sufficient quantity for a 8 well strip is 15 μ L of conjugate in 0,985 mL of diluent.

TEST PROCEDURE

MANUAL PROCEDURE

- Prior to starting the test, bring all reagents to room temperature (at least 2 hours).
- Add 100 μ L of positive control and 100 μ L of the negative control to each well (ready to use).
- Add 100 μ L of DILUTED sample into appropriate wells. Testing samples and controls in duplicates is recommended. Seal the plate and incubate for 10 minutes (\pm 2 min) at room temperature (18-25°C).
- Wash 4 times following the procedure previously described.
- Add 100 μ L of the 1/10 DILUTED conjugate to each well. Seal the plate and incubate for 10 minutes (\pm 2 min) at room temperature (18-25°C).
- Wash 4 times following the described procedure.
- Add 100 μ L of the substrate solution to each well. Keep the plate at room temperature for 5 minutes not exposed to light. In order to speed up this process, it is advisable to use a multichannel pipette.
- Add 100 μ L of the stop solution to each well. This reagent must be added following the same order in which the substrate was added.
- Read the OD of each well at 450 nm within 5 minutes after the addition of the stop solution.

AUTOMATED PROCEDURE

- Prior to starting the test, bring all reagents to room temperature (at least 2 hours). If performing the test on ELISA automatic systems, distribution and identification of samples and controls (duplicates recommended) should be established following manufacturer's instructions carefully. Select the required number of microtiter strips or wells and insert them into the holder.
- Add 100 µL of positive control and 100 µL of the negative control to each well (ready to use).
- Add 100 µL of DILUTED sample into appropriate wells.
- Incubate for 60 minutes (\pm 2 min) at 26°C.
- Wash 4 times following the procedure previously described.
- Add 100 µL of the 1/70 DILUTED conjugate to each well. Incubate for 30 minutes (\pm 2 min) at room temperature (26°C).
- Wash 4 times following the described procedure.
- Add 100 µL of the substrate solution to each well. Keep the plate at room temperature (26°C) for 15 minutes not exposed to light.
- Add 100 µL of the stop solution to each well.
- Read the OD of each well at 450 nm within 5 minutes after the addition of the stop solution.

READING AND RESULT INTERPRETATION

Measure and record the O.D. at 450nm for samples and controls. Determine the arithmetic mean (X) if test was performed in duplicates.

A. VALIDATION CRITERIA

Positive Control X (O.D.) \geq 0,80

Negative Control X (O.D.) \leq Negative cut off

B. INTERPRETATION

Negative cut off = O.D. positive control x 0,30

Positive cut off = O.D. positive control x 0,35

Samples will be considered:

- negative if its O.D. is equal or lower than the negative cut off.

- positive if its O.D. is equal or higher than the positive cut off.

- doubtful if its O.D. is between both cut off. In these cases, it is suggested to repeat the test after 3-4 weeks; if the O.D. remains lower than the positive cut off, sample will be considered as negative.

If titration of problem sample, title will be the last dilution showing an O.D. higher than the positive cut off.

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www.goldstandarddiagnostics.com

CERTIFIED

MANAGEMENT SYSTEM
QUALITY - ISO 9001
ENVIRONMENTAL - ISO 14001

9191.INGE

9175.INGE

Distribuido en
Distributed in

por
by