INgezim® PPA Compac 2.0

Prod Ref: 11.ASF.K3

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KIT COMPOSITION

<i>REAGENT</i>	<i>2 Plates</i>		<i>5 Plates</i>	
	Uni.	Vol.	Uni.	Vol.
Microtitration plates	2		5	
Vials with Positive Control	1	1mL	2	1mL
Vials with Negative Control	1	1mL	2	1mL
Bottles with Peroxidase Conjugate, ready to use	1	30mL	2	30mL
Bottles with Washing Solution 25x concentrated	1	125mL	1	125mL
Bottles with Diluent (DE31-01), ready to use	1	60mL	1	60mL
Bottles with Substrate (TMB), ready to use	1	30mL	1	60mL
Bottles with Stop Solution, ready to use	1	60mL	1	60mL

OTHER MATERIALS AND REAGENTS NEEDED NOT PROVIDED WITH THE KIT

Distilled or deionized water.
Micropipettes.
Disposable micropipette tips.
Washing plates device.
Test tubes.
ELISA Reader (450 nm filter).

TECHNICAL BASIS

This kit is based on the blocking enzymatic immunoassay technique (**Blocking ELISA**), briefly described below:

The plate wells were previously sensitized with the target antigen, the vp72 of the African swine fever virus, the major structural protein of the virus with high antigenic capability. The sera to be tested are added to the wells and, after a washing step, a specific monoclonal antibody to the vp72 protein labelled with peroxidase is added. When a serum sample contains specific antibodies against the virus, they will bind to the antigen absorbed on the plate, blocking the binding of the monoclonal antibody. If the serum sample does not contain specific antibodies, the monoclonal antibody labelled with peroxidase will be the one that binds to the coated antigen.

Non-fixed material is eliminated from the plate by several washing steps. Afterwards, the presence of bound monoclonal antibody is visualized by the addition of the chromogenic substrate, which acquires a measurable blue color in the presence of peroxidase. In this way, the appearance of color in the well indicates the absence of specific antibodies to the virus in the assayed serum, while the absence of color indicates the presence of specific antibodies.

PRECAUTIONS AND WARNINGS FOR USERS

1. Read the instructions for use carefully.
2. Bring all reagents to room temperature prior to use.
3. Do not mix instructions or reagents from different kits and avoid any contamination of the reagents.
4. Do not use components after expiration dates and do not mix components from different lots.
5. Do not eat, drink, or smoke while specimens or kit reagents are being handled and do not pipette by mouth.
6. Use a new tip for each serum sample.
7. The substrate is extremely sensitive to light and contamination, so it is recommended to remove only the volume to be used from the bottle and to never return the excess substrate to the bottle.
8. Stop solution is a strong acid solution. In case of accidental contact with skin, wash gently with water.

STORAGE OF COMPONENTS

Keep all the reagents between +2°C and +8°C.

Control sera are stable for one month. In case they are not going to be used in this period, we recommend storing them in aliquots at -20°C.

INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µL on each well. After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by tightly overturning the plate to avoid the possible mixture of different wells content. As a precaution, the wells should be emptied over a cuvette containing 0.1M NaOH solution, as the test samples may contain infectious agents.
- Dispense a volume of 300 µL of washing solution on each well.
- Softly mix the plate, avoiding the contamination between wells.
- Turn over the plate abruptly to empty the wells.
- Repeat the process as many times as indicated on the instructions of the kit. After the washing step, verify that the next reagent to be added to the plate is ready to use. Do not keep the plate dry.
- After the last washing step, shake it upside down on an absorbent filter paper to eliminate all washing buffer.

PREPARATION OF SAMPLES

Samples to be used in the assay must be porcine sera. **Do not use highly hemolyzed or contaminated samples. You might obtain false positive results.** Samples has to be assayed at a 1/2 dilution in the supplied serum diluent (DE31).

This dilution could be done directly in the plate by adding 50 µL of serum diluent and 50 µL of serum into each well.

PREPARATION OF REAGENTS

Washing solution:

Dilute one part of the concentrate washing solution provided in the kit with 24 parts of distilled or deionized water (i.e., 40 mL of concentrate plus 960 mL of dH₂O). When ready, this solution remains stable between +2°C and +8°C until the expire date described at the label of the concentrated solution.

Control sera:

The control sera has to be treated as test samples, diluted 1/2 in diluent prior to be used. This dilution can be done directly in the plate by adding 50 µL of diluent plus 50 µL of control sera. Once opened, control sera are stable for one month between +2°C and +8°C. In case they are not going to be used in this period, we recommend storing them in aliquots at -20°C.

TEST PROCEDURE

1. Before performing the assay, let all kit components reach room temperature (20°C to 25°C).
2. If dilution is directly prepared in the plate, add first 50 µL of the supplied serum diluent to each well. Then, add 50 µL of control sera and of the samples to be tested. Softly shake the plate for an appropriate homogenization of the samples, careful to avoid liquid transfer between wells. It is recommended to include control sera in duplicate for each assay. Seal the plate and **incubate for 1 h ± 5 min at 36 ± 1°C or overnight* (16-20 hours) at +4°C to + 8°C.**
3. Empty the wells into a container with NaOH 0.1M solution and wash 4 times as previously described.
4. Add 100 µL of specific conjugate to each well. Seal the plate and **incubate for 30 minutes at 36 ± 1°C.**
5. Wash 5 times as previously described.
6. Add 100 µL of substrate to each well. It is recommended to use multichannel pipettes for increased speed and homogeneity. Keep the plate for **15 min at room temperature (20°C to 25°C).**
9. Add 100 µL of stop solution to each well following the same order used for substrate dispensation.
7. Read the optical density at 450 nm of each well using an ELISA reader within 5 min after the addition of stop solution.

* The overnight incubation shows a higher sensitivity than the 1-hour protocol in samples with low antibody titer.

READING AND RESULTS INTERPRETATION

VALIDATION OF THE TEST:

Test results should be considered valid only if the OD of the Negative Control is, at least ,4 times higher than the OD of the Positive Control:

$$\frac{\text{Neg control DO (NC)}}{\text{Pos control DO (PC)}} \geq 4$$

RESULTS INTERPRETATION:

When samples are tested in duplicate, OD of the sample should be calculated as the arithmetic mean of OD values in both wells.

Cut off calculation for optical density:

The positive and negative cut off values are calculated as follows:

- Positive Cut Off = NC - [(NC - PC) x 0.5]
- Negative Cut Off = NC - [(NC - PC) x 0.4]

Where: NC = OD of the Negative Control Serum

PC = OD of the Positive Control Sera.

- Serum samples with an OD lower than or equal to the Positive Cut Off, are considered **POSITIVE** to ASFV antibodies.
- Serum samples with an OD higher than or equal to the Negative Cut Off, are considered **NEGATIVE** to ASFV antibodies.
- Serum samples with OD values between both cut offs are considered **DOUBTFUL**. We recommend re-testing these animals one more time or applying a different technique to check this serum (IFI, Indirect ELISA, etc.).

Cut off calculation for blocking percentage:

For calculating the blocking percentage (x %) of a sample you should proceed as follows:

$$X \% = \frac{\text{NC} - \text{Sample OD}}{\text{NC} - \text{PC}} \times 100$$

For blocking percentage, the following cut off should be applied:

- POSITIVE samples, X% ≥ 50%
- NEGATIVE samples, X% ≤ 40%
- DOUBTFUL samples, 40% < X% < 50%. We recommend re-testing these animals one more time or applying a different technique to check this serum (IFI, Indirect ELISA, etc.).

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