INgezim® PRRS Blue

Prod Ref: 11.PR3.K1

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KIT COMPOSITION

<i>REAGENT</i>	<i>2 Plates</i>		<i>5 Plates</i>	
	Uni.	Vol.	Uni.	Vol.
Microtitration plates	2		5	
Vials with Peroxidase Conjugate, ready to use	1	30mL	2	30mL
Bottles with Washing Solution 25x concentrated	1	125mL	1	125mL
Bottles with Substrate (TMB), ready to use	1	30mL	1	60mL
Bottles with Stop Solution, ready to use	1	60mL	1	60mL
Bottles with Diluent (DE35-01), ready to use	1	60mL	2	60mL
Vials containing Negative Control Sera, ready to use	1	4mL	2	4mL
Vials containing Positive Control Sera, ready to use	1	4mL	2	4mL

OTHER MATERIALS AND REAGENTS NEEDED NOT PROVIDED WITH THE KIT

Distilled or deionized water.
Micropipettes from 5 to 200 µL.
Disposable micropipette tips.
Washing plates device.
Test tubes from 50 to 250 mL.
ELISA Reader (650 nm filter)

TECHNICAL BASIS

This kit is based on an indirect enzymatic immunoassay (**Indirect ELISA**). The technique is briefly described below:

The recombinant antigen, specific for antibodies of both European and American strains of PRRS virus, is fixed on polystyrene plates. When a serum sample contains specific antibodies against the virus, they will bind to the recombinant antigen adsorbed on the plate. After a washing step to eliminate all non-fixed material, the presence of swine antibodies can be detected using a specific peroxidase conjugate. After the addition of the substrate a colorimetric reaction will appear which can be measured by a spectrophotometer at 650 nm.

Please note that our kit uses purified antigen obtained by the expression of the ORF7 from American and European strains of PRRS virus. This method warrants the total absence of infectivity in any of the components of the kit.

PRECAUTIONS AND WARNINGS FOR USERS

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature prior to use.
3. Do not mix instructions or reagents from different kits and avoid any contamination of the reagents.
4. Do not use components after expiration dates and do not mix components from different lots.
5. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled and do not pipette by mouth.
6. Use a new tip for each serum sample.
7. The substrate is extremely sensitive to light and contamination, so it is recommended to remove only the volume to be used and never return the excess substrate to the bottle.
8. The stop solution can sometimes precipitate when stored at 4°C. This must be checked, and precipitates must be dissolved before use.

STORAGE OF COMPONENTS

Keep all the reagents between +2°C and +8°C.

INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µL on each well. After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turnover of the plate to avoid the possible mixture of the content from one well to another. As a precaution, the wells should be emptied over a cuvette containing 0.1M NaOH solution, as the test samples may contain infectious agents.
- Dispense a volume of 300 µL of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit. After the washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry.
- After the last step of washing shake the plate turned over an absorbent filter paper.

PREPARATION OF SAMPLES

Sera samples must be tested at 1/40 dilution in diluent DE35-01 (i.e.: 5 µL of serum + 195 µL of serum diluent DE35-01).

PREPARATION OF REAGENTS

Washing solution:

Dilute one part of the concentrated washing solution provided in the kit with 24 parts of distilled or deionized water (40 mL of concentrated solution and 960 mL of water). Once prepared this solution remains stable when stored at +4°C.

Control sera (+) and (-):

Both controls are ready to use. **Do not dilute.**

Conjugate:

This reagent is ready to use. **Do not dilute.**

Serum diluent:

This reagent is ready to use. **Do not dilute.**

TEST PROCEDURE

1. All reagents must be allowed to come to room temperature before use.
2. Add 100 µL of samples diluted as it is specified in the previous instructions to the wells of the plate. Add 100 µL of positive and negative control serum. We recommend running samples and controls in duplicate for confirmatory purposes. Seal the plate and incubate for 30 min ± 2 at room temperature (20°C to 25°C).
3. Wash 3 times.
4. Add 100 µL of conjugate to each well. Seal the plate and incubate for 30 min ± 2 at room temperature (20°C to 25°C).
5. Wash 5 times.
6. Add 100 µL of substrate solution to each well. Keep the plate at room temperature (20°C to 25°C) for 15 min. In order to speed up this process, it is advisable to use a multichannel pipette. Count the time since the first well had been filled.
7. Add 100 µL of stop solution to each well following the same order in which the substrate was added.
8. Read the OD at 650 nm within 5 min after the addition of stop solution.

READING AND RESULT INTERPRETATION

Measure and record the **A (650)** for samples and controls. Determine the arithmetic mean ($\bar{X}$).

VALIDITY CRITERIA:

- $CP \bar{X} - CN \bar{X} \geq 0.30$
- $CN \bar{X} \leq 0.30$

S/P RATIO CALCULATE:

The presence or absence of antibody to PRRSV is determined by calculating the S/P ratio for each sample

$$S/P = \frac{\text{Sample A (650)} - NC \bar{X}}{PC \bar{X} - NC \bar{X}}$$

Negative S/P values must be considered with value 0.

INTERPRETATION:

- Samples with a S/P **higher or equal** than 0.4 should be considered positive to PRRS antibodies.
- Samples with a S/P **lower** than 0.4 should be considered negative to PRRS antibodies.

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Designed & manufactured in Spain by:

GOLD STANDARD DIAGNOSTICS MADRID SA
C/Hermanos García Noblejas, 41 2ª planta

28037- MADRID (SPAIN)

Tlf: +34 91368.05.01/04

Fax: +34 91 408.75.98

E-mail: info.spain@eu.goldstandarddiagnostics.com

www.goldstandarddiagnostics.com

Distribuido en
Distributed in

por
by