

INgezim ADV gE PLUS

R.11.GEP.K3

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
INgezim ADV gE PLUS

R.11.GEP.K3

INgezim ADV gE PLUS is a immunoenzymatic assay based on a blocking ELISA assay technique, which uses a monoclonal antibody (MAb) specific to Aujeszky Disease Virus (ADV) gE protein.

TECHNICAL BASIS OF THE KIT

- Plates are coated with inactivated ADV antigen. Serum samples are added and incubated.
- If the samples contain specific antibodies to ADV, they will bind to the antigen. (**Two options for incubation are possible: over night at 4°C or 2 hours at RT**)
- When the MAb-PO (ready to use) specific to ADV gE protein is added, it will bind to the protein only if there are no antibodies in the sample blocking the antigen (negative animals). In case there were antibodies blocking the antigen (infected animals), the conjugate would not be able to bind it. The binding is detected by the development of a colorimetric reaction after the addition of the substrate.

APPLICATION

Detection of specific antibodies to ADV gE protein in porcine sera samples. The assay allows the specific differentiation between antibodies produced by infection or vaccination. The assay has been developed to avoid specificity problems due to unspecific antibodies present in sera from breeding animals, without losing sensitivity in the test.

RESULT INTERPRETATION

Two cut off are used for the results interpretation: positive and negative. The samples will be considered as **Positive** when its OD value is equal to or lower than the positive cut off. The samples will be considered as **Negative** when its OD value is equal to or higher than the negative cut off. The samples with an OD value between both cut offs will be considered as **Doubtful**.

SENSITIVITY AND SPECIFICITY

1. *Using reference sera from O.I.E.*

The international reference serum ADV-GI (1/8) and the equivalent one are detected as required by the O.I.E.

2. *Specificity using sera from multi-vaccinated animals and from animals infected with different pathogens other than ADV.*

A set of 495 sera from multi-vaccinated animals with different deleted vaccines (GESKIPUR, NOVIPORVAC, SUVAXYN, AUGIVAC GP, NIA 4 and BEGONIA), and 40 sera from experimentally infected animals (with ASFV, Erysipelas, TGEV and PRRS) were analyzed with INgezim ADV gE. The results obtained were negative in 100% of the sera samples.

3. *Specificity using sera from breeding animals.*

509 sera from breeding animals obtained from 11 different problematic herds were analyzed. The results obtained showed **99.8% specificity and 98.9% sensitivity when these kinds of sera are analyzed.**

4. *Specificity using sera from experimentally vaccinated animals.*

Sera of 18 experimentally vaccinated animals with a deleted vaccine were analyzed. In the same study, sera of 8 non vaccinated animals were checked as controls. The results obtained indicated that **the assay was 100% specific as no positive results were obtained.**

5. *External validation*

This product has been included in the "SEROAUJ 2015", organized by the OIE Reference Laboratory (ANSES, France). In this report it is described that INgezim ADV Ge PLUS is able to detect a 1/128 dilution of the serum S96AUJ02 which indicates that it is able to detect also as positive the dilution 1/64 which is equivalent to dilution 1/8 of the OIE Reference Serum ADV-1.

INgezim ADV gE PLUS showed 100% correspondence with the expected results in this study

COMPOSITION OF THE KIT

- Microtitration plates of 96 wells.
- Vials with Positive Control
- Vials with Negative Control
- Vials with Peroxidase Conjugate
- Bottle with Washing Solution
- Bottle with serum diluent
- Bottle with stop solution
- Bottle with substrate (TMB) ready to use

PRODUCT MANUFACTURED BY INGENASA
REGISTRATION NUMBER 1126 RD

SHELF LIFE: 15 months
Stored 2°C-8°C

Ed. 020217