

## INgezim BRSV Compac

R.12.BRS.K3

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

## INgezim BRSV Compac

R.12.BRS.K3

**INgezim BRSV Compac** is an enzymatic assay based on a competition ELISA technique, which uses a monoclonal antibody (MAb) specific for Bovine Respiratory Syncytial Virus (BRSV), and a semi-purified protein extract of BRSV.

### TECHNICAL BASIS OF THE KIT

1. Plates are coated with BRSV antigen. SerUM samples and conjugate are added at the same time and incubated.
2. If the samples contain specific antibodies to BRSV (positive samples), they will compete with the conjugate for the union to the antigen and the conjugate will not be able to bind it. If there are not antibodies in the sample (negative animals), the conjugate will bind to the proteins coated. The binding is detected by the development of a colorimetric reaction after the addition of the substrate.

### APPLICATION

Detection of specific antibodies to BRSV, in bovine sera samples.


### INTERPRETATION OF THE RESULTS

One cut off is used for the results interpretation. The samples will be considered as **Positive** when their OD value is lower than the cut off. The samples will be considered as **Negative** when their OD value is higher than the cut off.

### VALIDATION

#### Diagnostic specificity.

To evaluate the diagnostic specificity a panel of 554 sera from BRSV free areas was analyzed. Results obtained indicated **99% diagnostic specificity**.

#### Comparison with a Reference ELISA

In order to evaluate the assay, a panel of 135 sera previously catalogued as positive by PrioCHECK® BRSV Ab were used. Results obtained indicated that the sensitivity of INgezim BRSV Compac in comparison with the Reference ELISA was 99%.

		PrioCHECK® BRSV Ab
		POS
INgezim® BRSV COMPAC	POS	134
	NEG	1

### COMPOSITION OF THE KIT

- Microtitration plates of 96 wells
- Vials with Positive Control
- Vials with Negative Control
- Vials with Peroxidase Conjugate
- Bottle with Washing Solution
- Bottle with diluent
- Bottle with stop solution
- Bottle with substrate (TMB) ready to use

REGISTRATION NUMBER 3046 RD  
PRODUCT MANUFACTURED BY INGENASA

SHELF LIFE: **18 months**  
Stored at 2°C-8°C

Ed.020217