

INgezim BTV COMPAC 2.0

R.12.BTV.K3

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

INgezim BTV COMPAC 2.0

R.12.BTV.K3

INgezim BTV Compac 2.0 is an enzymatic assay based on a blocking ELISA technique, which uses monoclonal antibody (Mab) specific for Blue Tongue Virus (BTV) VP7 protein, and a recombinant antigen.

TECHNICAL BASIS OF THE KIT

- Plates are coated with BTV antigen (recombinant VP7 protein). Samples (milk or serum) are added and incubated (Two options for incubation).
- If the samples contain specific antibodies to BTV, they will bind to the antigen.
- When Mab-PO specific of VP7 is added, only if there are no antibodies in the sample blocking the antigen (negative animals), it will bind to the protein. In case the sample contains antibodies blocking the antigen (infected animals), the conjugate will not be able to bind to it. The binding is detected by the development of a colorimetric reaction after the addition of the substrate.

APPLICATION

- Detection of specific antibodies to BTV VP7 protein.
- Use of milk and serum samples
- Two options for sample incubations depending on the sensitivity required and the kind of samples.
- Being a blocking assay, the detection of antibodies is independent of the specie, so it can be used for sheep, goats and cattle livestock.

INTERPRETATION OF THE RESULTS

Two cut offs are used for the results interpretation: Samples are considered **Positive** if their OD value at 450nm is equal to or lower than the positive cut off. Samples are considered **Negative** if their OD value at 450 nm is equal to or higher than the negative cut off. Samples with an OD value within the range of both cut offs is considered **Doubtful**.

VALIDATION

SERUM SAMPLES

1. Specificity

1262 sera (92 sheep, 92 goats y 1078 cattle) from BTV free Spanish herds were analysed. The results indicate that the specificity is **100%**.

2. Diagnostic Sensitivity (Experimental Studies)

Animals infected with the 24 serotypes of BTV were studied. The results obtained demonstrate that the assay is able to detect antibodies specific of the 24 serotypes.

3. Diagnostic Sensitivity (Field Studies)

A set of 172 samples previously classified as positive for BTV, were analysed. All of them showed positive results, except one. Hence, the sensitivity of the assay in this study was 99.4%.

MILK SAMPLES

In order to evaluate the performance of the assay using milk samples, a set of those samples were analysed.

- 443 individual samples previously classified as negative.
- 139 individual samples previously classified as positive.

The results indicate a specificity of 100% and sensitivity higher than 94%.

COMPOSITION OF THE KIT

- Microtitration plates of 96 wells
- Vial with BTV Positive Control.
- Vial with BTV Negative Control.
- Vial with Conjugate specific for BTV.
- Bottles with Washing Solution concentrated.
- Bottles with diluent ready to use.
- Bottle with substrate (TMB) ready to use.
- Bottle with stop solution.

REGISTRATION NUMBER 1125-RD

PRODUCT MANUFACTURED BY INGENASA

Shelf life: **18 months**
Store at 2°C-8°C

Ed.020217