

INgezim Neospora caninum

R.12.NC.K1

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

INgezim Neospora caninum

R.12.NC.K1

INgezim Neospora caninum is an immunoenzymatic assay based on an indirect ELISA technique to detect specific antibodies to *Neospora caninum* in bovine samples.

TECHNICAL BASIS OF THE KIT

1. Plates are coated with inactivated *Neosporacanium* antigen. Samples are added to the wells and incubated.
2. If the sample contains antibodies to *Neospora caninum*, they will bind to the antigen.
3. After adding a Polyclonal-HRPO specific for bovine IgG, it will bind to the IgGs of the sample previously bound to the antigen. This binding is detected by the development of the colorimetric reaction after the addition of the substrate.

APPLICATION

Detection of specific antibodies to *Neospora caninum* in bovine sera samples.

INTERPRETATION OF THE RESULTS

The assay uses two different cut offs. Samples will be considered **Positive** if their OD value is higher than or equal to the positive cut off; **Negative** if their OD value is lower than the cut off; **Doubtful** if their OD are within the range of both cut offs.

VALIDATION

1. Correlation with Canadian Gold Standard ELISA

		Reference ELISA			
		POSITIVE	NEGATIVE	DOUBTFUL	TOTAL
INGEZIM NEOSPORA	POSITIVE	135	0	0	135
	NEGATIVE	1	875	1	876
	DOUBTFUL	17	8	0	25
	TOTAL	153	883	0	1036

- o Considering doubtful sera as positive, the relative sensitivity is 99,4% (95% CI=096,1-100%) and the relative specificity, 100% (95% CI= 99,7-100%)
- o Considering doubtful sera as negative, the relative sensitivity is 88,2% (95% CI=81,8-92,7%) and the relative specificity, 99,1% (95% CI= 98,1-99,6%)

COMPOSITION OF THE KIT

- 96 wells microtitration plates
- Vials with positive control
- Vials with negative control
- Vials with conjugate
- Bottle with washing solution
- Bottle with diluent
- Bottle with substrate
- Bottle with stop solution

REGISTRATION NUMBER 0800-RD

PRODUCT MANUFACTURED BY INGENASA

IT-73840 IT-73780 ISO 14001:2015 9191.INGE ISO 9001:2015 9175.ING2

Shelf life: **18 months**
Store at 2°C-8°C

Ed.020217