
## INgezim AHSV Compac PLUS

R.14.AHS.K3

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

## INgezim AHSV Compac PLUS

R.14.AHS.K3

**INgezim AHSV COMPAC PLUS** is an enzymatic assay based on a blocking ELISA technique, which uses a monoclonal antibody (MAb) specific to African Horse sickness Virus (AHSV) VP7 protein, and a recombinant antigen (VP7 protein of AHSV).

### TECHNICAL BASIS OF THE KIT

1. Plates are coated with recombinant antigen (VP7 protein of AHSV). Serum samples are added and incubated.
2. If the samples contain specific antibodies to AHSV VP7 protein, they will bind to the antigen.
3. When a MAb-PO specific of AHSV VP7 protein is added, only if there are no antibodies in the sample blocking the antigen (negative animals), it will bind to the protein. In case the sample contains antibodies blocking the antigen (infected or vaccinated animals), the conjugate will not be able to bind to it. The binding is detected by the development of a colorimetric reaction after the addition of the substrate.

### APPLICATION

Detection of specific antibodies to AHSV VP7 protein in equine serum samples

### INTERPRETATION OF THE RESULTS

Two cut offs are used for the results interpretation. Samples will be considered **Positive**, if their OD value is equal to or lower than the positive cut off; **Negative**, if their OD value is equal to or higher than the negative cut off and **Doubtful** if their OD value is within the range of both cut offs.


### VALIDATION

#### DIAGNOSTIC SENSITIVITY & SPECIFICITY

##### HORSES

9 sera from field horses infected by AHSV and 152 sera from AHSV vaccinated horses, previously classified as positive by seroneutralization (SN), were analysed. The results obtained indicated 97.5% sensitivity in relation to SN being the doubtful sera those from vaccinated animals.

A set of 698 sera of AHSV negative horses (514 from free zones and 184 negative by SN) were analysed. The results obtained indicated 99.5% specificity.

■ Negative; ■ vaccinated; ■ infected

##### DONKEYS

There were also analysed 16 sera of donkeys classified as positive and 6 as negative, giving all of them the expected results.

#### ANALYTICAL SENSIBILITY

A set of sera of animals experimentally infected with different serotypes of AHSV, were analysed (9 sera of Guinea Pig infected with the 9 serotypes of AHSV, 8 goat infected with 8 serotypes of AHSV except serotype 4, 1 rabbit infected with serotype 4 of AHSV). Results obtained indicated that the assay is able to detect antibodies specific of the 9 AHSV serotypes. Moreover, 65 sera of camels negative to Ab specific of AHSV (although they don't suffer the disease, AHSV Ab specific can be found) were analysed, being the specificity 100%.

#### ANALYTICAL SPECIFICITY

A set of sera of animals experimentally infected with different agents near AHSV, were analysed (6 sera of Guinea Pig infected with serotypes 1, 2, 13, 14, 22 and 23 of BTV; 4 sheep infected with serotypes 19, 5, 10 and 20 of BTV; 1 goat infected with Akabane virus and 3 cows infected with serotypes 1, 2 and 3 of EHDV). Results obtained indicated that the assay does not cross-react with antibodies to any of the related agents.

#### OIE REFERENCE SERA

The following OIE standard sera from CISA have been used: negative horse, positive and weak positive horses, positive and weak positive mules, positive and weak positive donkeys. The results obtained were the expected.

#### COMPOSITION OF THE KIT

- Microtitration plates of 96 wells
- Vials with Positive Control
- Vials with Negative Control
- Vials with Peroxidase Conjugate
- Bottle with Washing Solution
- Bottle with Diluent
- Bottle with Stop Solution
- Bottle with Substrate (ABTS)
- Bottle with Substrate Buffer

REGISTRATION NUMBER 1390 RD  
PRODUCT MANUFACTURED BY INGENASA

SHELF LIFE: **18 months**  
Stored at 2°C-8°C

Ed. 250117