

INgezim West Nile IgM

R.14.WNV.K2

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

INgezim WEST NILE IgM

R.14.WNV.K2

INgezim WEST NILE IgM is a capture immunoenzymatic assay which uses a monoclonal antibody specific to equine IgM and a monoclonal antibody specific to WNV E protein (Domain III).

TECHNICAL BASIS OF THE KIT

1. Plates are coated with a monoclonal antibody (MAb) specific of equine IgM. Samples are added and incubated.
2. If the sample contains specific IgMs to WNV, they will be captured by the α -equine IgM MAb which is coating the plate.
3. Once the WNV antigen is added, the E protein will be captured by the specific IgM from the sample. A negative antigen is also used to determine unspecific reactions.
4. The WNV specific MAb-HRPO added later will bind to the antigen in case it has been captured by the specific IgM of the sample (infected or vaccinated animals). This binding is detected by adding a specific substrate.

APPLICATION

Detection of specific IgM antibodies to West Nile Virus E protein

INTERPRETATION OF THE RESULTS

The assay establishes two cut offs. Samples will be considered **Positive**, if their Positivity index (IP= ODPositive Ag- ODNegative Ag) is higher than the positive cut off; **Negative**, if their IP value is lower than the negative cut off. IP in the range between the positive and negative cut offs will be considered doubtful.

VALIDATION

VALIDATION WITH REFERENCE SERA

- Two reference Sera from APHIS (USDA) were used:
- Equine serum infected with strain N.Y.99, positive by seroneutralization.
 - Equine serum infected with strain N.Y.99, IgM positive.

Both sera were detected as positive.

ANALYTICAL SENSITIVITY

Different extractions made to vaccinated horses were analyzed. Animals were bled at days 0, 14, 21, 28, 42, 56, 70 and 85 p.v. The assay was able to detect specific IgM from **day 14 p.v.**


RING TRIAL (Proficiency test ANSES 2013)

The conclusions of the Final Report regarding sensitivity were: The Ingezim WN IgM kit appeared to be more sensitive than IDEXX, IDVet or in-house IgM capture ELISAs with:

- A serum sampled from the WNV lineage 1-infected pony on day 35 post-infection and diluted at 1/6. Found positive with the Ingezim WN IgM kit and negative with the 3 other tests.
- A serum sampled from the WNV lineage 2-infected pony and diluted at 1/24 and 1/48. Found positive with the Ingezim WN IgM kit, doubtful with the IDEXX IgM WNV Ab kit and negative with the IDScreen WN IgM Capture or in-house ELISAs.

DIAGNOSTIC SPECIFICITY AND SENSITIVITY

- **Study 1** 14 equine sera infected by WNV, 1 SN positive reference serum, 1 WNV specific IgM positive serum, 3 equine WNV vaccinated sera, positive by other commercial assays, 404 equine sera from WNV free areas and 28 field WNV infected sera from 3 outbreaks in Spain (2010-2012) and positive by another commercial ELISA. **The sensitivity and specificity of the assay were higher than 99.5%.**

- **Study 2:** 17 equine sera from a Linage 2 outbreak occurred in Italy in 2013 and previously classified as positive. INgezim WNV IgM detected 16 of them as positive.

COMPOSITION OF THE KIT

- Microtitration plates of 96 wells
- Vials with positive antigen
- Vials with negative antigen
- Vials with Positive Control
- Vials with Negative Control
- Vials with Peroxidase Conjugate
- Bottle with Washing Solution
- Bottle with diluent
- Bottle with stop solution
- Bottle with substrate (TMB) ready to use

PRODUCT MANUFACTURED BY INGENASA
REGISTRATION NUMBER 2863 RD

IT-73840 ISO 14001:2015 ISO 9001:2015
IT-73780 9191.INGE 9175.ING2

SHELF LIFE: **12 months**
Stored at 2°C-8°C

Ed.250117