

## INgezim PARVO CANINO

R.15.CPV.K1

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

## INgezim PARVO CANINO

R.15.CPV.K1

**Ingezim PARVO CANINO** is based on an indirect enzymatic immunoassay (indirect ELISA), which uses a monoclonal antibody (MAb), specific to canine immunoglobulins (Igs) and a canine Parvovirus (CPV) VP2 recombinant protein as antigen.

### TECHNICAL BASIS OF THE KIT

1. Plates are supplied coated with CPV antigen. Serum samples are added and incubated.
2. If the samples contain specific antibodies to CPV, they will bind to the antigen.
3. When a MAb-PO specific to canine immunoglobulins is added, it will bind to the Igs previously bound to the antigen. The binding is detected by the development of a colorimetric reaction after the addition of the substrate.

### APPLICATION

Detection and/or titration of specific antibodies to CPV in canine sera samples.

### INTERPRETATION OF THE RESULTS

One cut off is used for the results interpretation:

**Screening:** IP values higher than the cut off are considered **Positive**; IP values lower than the cut off are considered **Negative**.


**Titration:** The titre of each sample will be the last dilution at which the IP value is higher than the cut off.

### VALIDACIÓN


1. **Correlation with inhibition haemagglutination test (IHA). Field sera.** 100 canine field sera were analysed in this study. The results obtained indicated 95% sensitivity and a specificity percentage higher than 99%.
2. **Analytical sensitivity in relation to IHA:**
  - a. **Study 1:** 20 sera of animals experimentally vaccinated and challenged were analyzed. 19 extractions were taken and analyzed by both assays. The results obtained indicated that ELISA is able to detect specific antibodies of CPV at day 10p.v. in 35% of cases and IHA in 25%.
  - b. **Study 2:** 10 sera of experimentally vaccinated dogs were analyzed. 20 extractions were taken and analyzed by both assays. The results obtained indicated that the ELISA is able to detect specific antibodies of CPV at day 10p.v. in 20% of cases and IHA 0%.

### 4. ELISA & IHA titres correlation.

- a. **Field studies:** The graph indicates the correlation between titres of both techniques in 61 animals analyzed.

- b. **Experimental studies:** The graph indicates the titre correlation in a study made with 20 dogs experimentally vaccinated and challenged at day 73 post vaccination.

### COMPOSITION OF THE KIT

- 96 wells microtitration plates
- Vials with Positive Control
- Vials with Negative Control
- Vials with conjugate
- Bottle with washing solution
- Bottle with diluent
- Bottle with substrate
- Bottle with stop solution

PRODUCT MANUFACTURED BY INGENASA

Shelf life: **18 months**  
Store at: 2°C-8°C

Ed. 02017