

INgezim INFLUENZA A

R.10.FLU.K3

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

INgezim INFLUENZA A

R.10.FLU.K3

INgezim INFLUENZA A is an enzymatic assay based on a blocking ELISA technique, which uses a monoclonal antibody (MAb) specific to Influenza A virus

TECHNICAL BASIS OF THE KIT

- Plates are coated with inactivated Influenza A antigen. Sera samples are added and incubated.
- If the samples contain specific antibodies to Influenza A, only if there are no antibodies in the sample blocking the antigen (negative animals), it will bind to the protein. In case the sample contains antibodies blocking the antigen (infected animals), the conjugate will not be able to bind to it. The binding is detected by the development of a colorimetric reaction after the addition of the substrate.

APPLICATION

Detection of specific antibodies to the nucleoprotein of Influenza A virus, in avian, porcine and equine serum samples. The assay is useful for the detection of antibodies to different strains of Influenza A virus.

INTERPRETATION OF THE RESULTS

Two cut off are used for the results interpretation: positive and negative. Samples will be considered **Positive** if their OD value is equal to or lower than the positive cut off. Samples will be considered **Negative** if their OD value is equal to or higher than the negative cut off. Samples with an OD value between both values will be considered **Doubtful**.

VALIDATION

POULTRY

Sensitivity

Sera from experimentally vaccinated animals (turkeys and chickens) and vaccinated birds from a zoo were analysed. Different vaccines and different doses were used (H7N4, H5N9, H5N9+H7N4, H5N9+H7N1), showing all of them positive results, indicating that the assay is able to detect antibodies to the nucleoprotein of different strains of Influenza A virus.

Specificity

613 samples from different sources (turkey, ostrich, hen, chicken, partridge, vulture, flamingo, kite, etc) and catalogued as negative were analysed and the specificity was 98%.


SWINE

Sensitivity and specificity

Sera from 303 experimentally vaccinated animals and 75 control animals were analysed. Different vaccines were used (H1N1, H1N2, H3N2) Results obtained indicated that the assay is able to detect antibodies to the nucleoprotein of different strains of Influenza A virus showing 88% sensitivity and 97% specificity.

Correlation with IHA

A set of positive and negative sera classified by IHA were analysed by INgezim Influenza A and the results obtained indicated a very accurate correlation in relation to the reference technique


HORSES

Correspondence with IHA

142 equine samples with different titers of IHA have been analyzed by ELISA INgezim® Influenza A and by ID Screen® Influenza A. Results obtained are shown in the table below:

	IHA	
	Sen	Spe
INgezim®	97%	92%
ID Screen®	98%	73%

The assay has participated in the "2016 PTS porcine influenza detection in serum (GD)" and in the "2016 PTS avian influenza antibody detection in chicken serum (GD)"

COMPOSITION OF THE KIT

- Microtitration plates of 96 wells
- Vials with Positive Control
- Vials with Negative Control
- Vials with Peroxidase Conjugate
- Bottle with Washing Solution
- Bottle with Diluent
- Bottle with Stop Solution
- Bottle with Substrate (TMB)

SPANISH REGISTRATION NUMBER
FOR POULTRY 1058 RD
PRODUCT MANUFACTURED BY INGENASA

SHELF LIFE: **15 months**
Stored at 2°C-8°C

Ed. 100117