

INGEZIM PPA Compac

R.11.PPA.K3

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
INGEZIM PPA Compac

R.11.PPA.K3

INGEZIM PPA COMPAC is an enzymatic assay based on a blocking ELISA assay technique, which uses a monoclonal antibody (MAb) specific to African Swine Fever Virus (ASFV) VP72 protein.

TECHNICAL BASIS OF THE KIT

1. Plates are coated with inactivated ASFV antigen (VP72 protein). Sera samples are added and incubated.
2. If the samples contain specific antibodies to ASFV, they will bind the antigen.
3. When a MAb-PO specific of VP72 is added, only if there are no antibodies in the sample blocking the antigen (negative animals), it will bind to the protein. In case the sample contains antibodies blocking the antigen (infected animals), the conjugate will not be able to bind to it. The binding is detected by the development of a colorimetric reaction after the addition of the substrate.

APPLICATION

Detection of specific antibodies to African Swine Fever Virus, in porcine sera samples.

INTERPRETATION OF THE RESULTS

Two cut offs are used for the results interpretation. Samples will be considered **Positive** if their OD value is equal to or lower than the positive cut off. Samples will be considered **Negative** if their OD value is equal to or higher than the negative cut off. Samples with OD value between both cut offs will be considered **Doubtful**.

DIAGNOSTIC SENSITIVITY

208 sera from endemic countries in West and East Africa and from old outbreaks in Spain were analyzed. These sera had previously been classified by the OIE Reference Indirect ELISA as positive samples. The results obtained indicated a 99% sensitivity respect the OIE Technique.

ANALYTICAL SENSITIVITY

Samples from pigs experimentally infected with different isolates of ASFV were analyzed. The isolates used for the experimental infection were: E70/E75 p72 genotype I (Spain); p72 genotype I (Sardinia); Ken05.Tk1 p72 genotype X (Kenya); p72 genotype IX (Kenya); p72 genotype II (Armenia); Malau 83, genotype VIII and Rep. Dominicana 78, genotype I.

The results obtained indicated that the assay is able to detect antibodies specific of all of these isolates. These antibodies are detected between days 10 and 25 post infection, depending on the isolate used.

DIAGNOSTIC SPECIFICITY

1. Testing sera from endemic areas in Africa
167 African pig sera from ASFV endemic areas were analysed. These sera had been classified as negative by the OIE reference Indirect ELISA. The results obtained showed 100% specificity.
2. Testing sera from ASFV non infected areas
1043 sera from ASFV non-infected areas were analysed. These sera had been classified as negative by the OIE reference Indirect ELISA. The results obtained showed 100% specificity.
3. Testing problematic sera
A set of 23 samples with negative results by Immunoblotting and doubtful results by the Indirect ELISA was analysed. The results obtained by INGEZIM PPA COMPAC showed a correspondence of 100% with the Immunoblotting technique

COMPOSITION OF THE KIT

- Microtitration plates of 96 wells
- Vials with Positive Control
- Vials with Negative Control
- Vials with Peroxidase Conjugate
- Bottle with Washing Solution
- Bottle with diluent
- Bottle with stop solution
- Bottle with substrate (TMB) ready to use

REGISTRATION NUMBER 335 RD
PRODUCT MANUFACTURED BY INGENASA

SHELF LIFE: **18 months**
Stored at 2°C-8°C

Ed.020217