



eurofins

Technologies

GENE*Spin*

**Kit for isolation of high-quality DNA
from food and feed samples**

Cat. no. 5224400605

For 50 DNA extractions

Table of Contents:

1.	GENESPIN - INTRODUCTION	3
2.	KIT CONTENTS OF GENESPIN	4
3.	MATERIALS NOT CONTAINED IN THE KIT	5
4.	EQUIPMENT NEEDED	5
5.	PRODUCT DESCRIPTION	6
6.	THE BASIC PRINCIPLE	7
7.	STORAGE AND HOMOGENIZATION OF SAMPLES	9
8.	ELUTION PROCEDURES	10
9.	PREPARATION AND STORAGE OF WORKING SOLUTIONS	11
10.	PROTOCOL FOR FOOD SAMPLES	12
10.1	Important hints and remarks	16
10.2	Special hints for difficult matrices	17
11.	TROUBLESHOOTING	18
12.	SAFETY INSTRUCTIONS	20
13.	PRODUCT WARRANTIES AND SATISFACTION GUARANTEE	24
14.	IMPORTANT NOTES	25
15.	TECHNICAL SERVICE	26

Kits, their components and instructions
for use are subject to alterations.

GENESpin

**Kit for DNA isolation
from food and feed**

Cat. no. 5224400605

1. GENESpin - Introduction

GENESpin provides an optimal lysis and DNA extraction system for nearly all types of food samples. Resulting eluates can be used for all kinds of subsequent detection methods, e.g. real-time and conventional PCR (for very difficult sample types causing PCR inhibition, additional use of our DNA Cleaning Columns, cat. no. 5224700310, is recommended).

Eurofins GeneScan Technologies GmbH is a producer of test kits for food and feed analysis offering DNA extraction kits and test kits for genetically modified organisms (GMOs), food pathogens, plant and animal species.

GENESpin silica membrane spin technology allows fast and effective purification of nucleic acids from various matrices.

The silica membranes are optimized for high DNA recoveries and low binding efficiencies for impurities. For further information regarding DNA purification, please feel free to contact us.

2. Kit contents of **GENESpin**

GENESpin Lysis Buffer, 100 mL

GENESpin Binding Buffer, 30 mL

GENESpin Wash Buffer 1, 30 mL

GENESpin Wash Buffer 2, 12 mL

GENESpin Elution Buffer, 13 mL

GENESpin Proteinase Buffer, 1.8 mL

GENESpin Proteinase K, lyophilized, 6 mg

50 GENESpin Columns

200 Collection tubes (2 mL)

Attention:

GENESpin Binding Buffer and GENESpin Wash Buffer 1 contain guanidinium hydrochloride and/or detergents!

Wear gloves and goggles!

Guanidinium hydrochloride can form highly reactive compounds when combined with bleach (sodium hypochlorite). DO NOT mix bleach or acidic solutions with sample preparation waste!

GENESpin

Kit for DNA isolation from food and feed

Cat. no. 5224400605

3. Materials not contained in the kit

- 96 – 100% ethanol
- 1.5 mL microcentrifuge tubes
- Pipet tips

4. Equipment needed

- Water bath or incubator for 65°C and 70°C
- Heating block for 65°C
- Pipettors
- Microcentrifuge
- Vortex mixer
- Equipment for sample disruption/
preparation/homogenization

5. Product description

- GENESpin is designed for the isolation of genomic DNA from food and feed samples of plant and animal origin.
- GENESpin allows processing of up to 200 mg material (larger amounts are possible through upscaled use of lysis buffer). Typical yields for GENESpin are in the range of 0.1-10 µg DNA, however, are sample-dependent.
- The eluted DNA is ready for use in subsequent reactions like PCR or real-time PCR. (For very demanding matrices, additionally use DNA Cleaning Columns, # 5224700310.)

Sample material:	200 mg (or upscaled)
Fragment size:	300 bp – approx. 50 kbp
Binding capacity:	30 µg
Typical DNA yield:	0.110 µg
A260/A280:	1.6 – 1.9
Elution volume:	100 µL
Preparation time:	30 min/6 preps

GENESpin

Kit for DNA isolation from food and feed

Cat. no. 5224400605

6. The basic principle

DNA is extracted from homogenized food samples with a lysis buffer containing chaotropic salts, denaturing agents and detergents. The standard isolation procedure ensures lysis with *GENESpin* Lysis Buffer. Lysis mixtures are cleared by centrifugation or filtration in order to remove contaminations and residual cellular debris. The clear supernatant is mixed with the *GENESpin* Binding Buffer and ethanol to create conditions for optimal binding to the *GENESpin* silica membrane, which was selected for this purpose due to its unique DNA binding properties.

After washing with two different buffers for efficient removal of potential PCR inhibitors, DNA is eluted in low salt buffer or water and is ready-to-use for subsequent analysis or amplification.

Food samples are very heterogeneous and contain many different compounds like fat, cocoa or polysaccharides which can lead to

suboptimal extraction or subsequent analysis of DNA.

GENE*Spin* guarantees good recovery rates also for small genomic DNA fragments < 1 kb from processed, complex food matrices (e.g. ketchup or spice) with very low DNA content and degraded DNA.

In consequence, it is recommended to design primers for amplification of short DNA fragments (80-150 bp) for subsequent PCR.

GENESpin

Kit for DNA isolation from food and feed

Cat. no. 5224400605

7. Storage and homogenization of samples

The lysis procedure is most effective with well homogenized samples with small particle size. Suitable methods include grinding with pestle and mortar in the presence of liquid nitrogen or with steel beads, as well as any type of commercial homogenizer, bead mill etc. After homogenization and treatment of the sample with GENESpin Lysis Buffer, the lysate can be cleared by centrifugation or filtering.

Methods for sample homogenization/ grinding

Mortar and pestle in the presence of liquid nitrogen

- Commercial homogenizers, e.g. bead mills
- VA steel beads (7 mm diameter, 4-5, beads in 15 ml tube): Chill the tube in liquid nitrogen. Vortex about 30 s. Repeat procedure until the entire sample is ground to powder. Remove the beads gently with a magnet. Keep the material frozen throughout the whole procedure. Do not use nitrogen inside the tube to avoid sticking of sample material to the beads.

8. Elution procedures

It is possible to adapt elution method and volume of GENESpin Elution Buffer to the subsequent application as follows:

- **Complete yields:** 90-100% of bound nucleic acids can be eluted by performing two elution steps with 2 x 100 μ L. Combine eluates and measure yield.
- **Highly concentrated eluates:** With minimal elution volumes (25 - 50 μ L), 60 - 80% of bound nucleic acids can be eluted, resulting in highly concentrated eluates. GENESpin Elution Buffer can be replaced by TE buffer or water (adjust pH to 8 - 8.5).

GENESpin

**Kit for DNA isolation
from food and feed**

Cat. no. 5224400605

9. Preparation and storage of working solutions

Store kit at room temperature (18 - 25°C) before first use. Storage at 4°C is possible but may cause precipitation of salts in buffers. If precipitations occur, warm buffer to 25 - 37°C to dissolve the precipitate before use.

Prepare the following reagents before first use:

- **GENESpin Wash Buffer 2:** Add 48 mL of 96-100% ethanol to GENESpin Wash Buffer 2, mark the label of the bottle to indicate that ethanol was added. Store at room temperature for up to one year.
- **GENESpin Proteinase K:** Before the first use of the kit, add 600 µL of GENESpin Proteinase Buffer to dissolve the lyophilized GENESpin Proteinase K. Proteinase K solution is stable for 6 months at - 20°C.

10. Protocol for food samples

1 Homogenize samples

Homogenize about 0.2 g material with a commercial homogenizer.

homogenize samples


2 Lyse cells

Transfer the resulting powder to a 2 mL collection tube. Preheat GENESpin Lysis Buffer to 65°C immediately before use and add 550 µL GENESpin Lysis Buffer (65°C). Mix carefully (15 s), add 10 µL GENESpin Proteinase K and mix again (2-3 s).

If the lysis buffer volume is not large enough to dissolve the sample completely, add more buffer (and Proteinase K proportionally) until the sample is totally resuspended.

Incubate at 65°C for 30 min.

Optional: add 10 µL RNase A (20 mg/ml) per 550 µL lysis buffer, mix well, incubate for 30 min at room temperature.



+ 550 µL
GENESpin Lysis
Buffer
(65°C)

+ 10 µL Prot. K



65°C 30 min

*(optional:
+ 10 µL RNase
A per 550 µL
GENESpin Lysis
Buffer,
RT 30 min)*

GENESpin

Kit for DNA isolation from food and feed

Cat. no. 5224400605


Afterwards, centrifuge the mix for 10 min ($> 10.000x g$) to pellet contaminants and cell debris.



*10 min
> 10.000 x g*

3 Adjust DNA binding conditions

Transfer the clear supernatant into a new centrifuge tube capable of holding at least 3 sample volumes. Add 1 volume GENESpin Binding Buffer and 1 volume ethanol (e.g. for 300 μ L sample: add 300 μ L binding buffer + 300 μ L ethanol). Vortex the mixture for 30 sec.



*+ 1 vol.
GENESpin
Binding Buffer
and
+ 1 vol. ethanol*

vortex 30 s

4 Bind DNA

For each sample, place a GENESpin Column into a new 2 mL collection tube and pipette 700 μ L mixture onto the column. Centrifuge for 1 min at $11.000 x g$. Discard flow-through. Repeat the procedure with remaining sample.



*load
sample
(700 μ L)*



*1 min
11.000 x g*

5 Wash and dry

1st washing step

Pipette 400 μ L GENESpin Wash Buffer 1 onto the GENESpin Column. Centrifuge for 1 min at 11.000 x g. Discard flow-through.



+ 400 μ L
GENESpin
Wash Buffer 1



1 min
11.000 x g

2nd washing step

Pipette 700 μ L GENESpin Wash Buffer 2 onto the GENESpin Column. Centrifuge for 1 min at 11.000 x g. Discard flow-through.



+ 700 μ L
GENESpin
Wash Buffer 2



1 min
11.000 x g

3rd washing step

Pipette 200 μ L GENESpin Wash Buffer 2 onto the GENESpin Column. Centrifuge for 2 min at 11.000 x g in order to remove Wash Buffer 2 completely (Residual ethanol from Wash Buffer 2 may inhibit enzymatic reactions). Discard flow-through.



+ 200 μ L
GENESpin
Wash Buffer 2



2 min
11.000 x g

GENESpin
Kit for DNA isolation
from food and feed
Cat. no. 5224400605

6 Elute DNA

Place the GENESpin Column in a **new** 1.5 mL centrifuge tube.



new 1.5 mL tube

Pipette 100 μ L GENESpin Elution Buffer (preheated to 70°C) onto the membrane.

*+ 100 μ L
GENESpin
Elution Buffer
(70°C)*

Incubate for 5 min at room temperature (18 - 25°C).



5 min RT

Centrifuge for 1 min at 11.000 x g to elute the DNA.



*1 min
11.000 x g*

10.1 Important hints and remarks

Considering the generally rather low DNA content in processed food, this protocol should be started with at least 200 mg of sample material.

Lysis buffer was tested for extraction of DNA from various types of samples including food of plant and animal origin.

RNase A use

RNase A (not included in the kit) addition is recommended for RNA-rich samples (i.e. unprocessed samples of plant or animal origin). Add 10 μL of 20 mg/mL stock solution (or equivalent amount) per 550 μL lysis buffer in step 2 or perform a RNase A digestion in the eluate before further use.

A vacuum manifold can optionally be used for acceleration of washing steps. Loading and elution steps should be done with centrifugation as described in the protocol.

GENESpin

Kit for DNA isolation from food and feed

Cat. no. 5224400605

10.2 Special hints for difficult matrices

- Ketchup, sauce and similar fluid samples (0.2 g equivalents) can be mixed with GENESpin Lysis Buffer, 500-1000 μ L each, and incubated with GENESpin Proteinase K as described in the protocol.
- Powdered hygroscopic samples more GENESpin Lysis Buffer than indicated can be used until the lysate is at least semifluid and can be pipetted. Extraction can be improved by pre-incubation of sample with GENESpin Lysis Buffer for 1-2 h.

If according to local regulations, defined amounts of sample have to be analyzed, higher amounts of sample (e.g. 1 or 2 g) can be used with upscaled lysis buffer volumes. We recommend to use only one 300 μ L aliquot of the clear supernatant per GENESpin Column. Otherwise, prepare 2 aliquots and load them step by step onto the GENESpin Column. GENESpin Lysis Buffer can be ordered separately (Cat no. 5224701901).

11. Troubleshooting

Problem	Possible cause
DNA yield is low	<ul style="list-style-type: none"> • Grinding/homogenization of food material was not sufficient • Extraction of DNA from sample during lysis was not sufficient. • Sample contains too much RNA • Suboptimal elution
DNA is degraded	<ul style="list-style-type: none"> • Sample contamination with DNase • Centrifugation speed was too high
DNA purity is poor	<ul style="list-style-type: none"> • DNA contaminants are not fully removed

GENESpin
Kit for DNA isolation
from food and feed
Cat. no. 5224400605

Suggestions



- For most matrices grinding with commercial bead mills/ mixers/homogenizers or steel beads is recommended.
- For higher yields of DNA, the lysis incubation time can be prolonged (up to overnight).
- Add 10-20 μ L of RNase A solution to the Lysis Buffer before heat incubation. If this is not successful, add the enzyme to the cleared lysate and incubate for 30 min at 37°C.
- The DNA can be either eluted in higher volumes (up to 300 μ L) or by repeating the elution step up to three times. Remember that the *GENESpin* Elution Buffer must be preheated to 70°C.
- Check if the pH of your elution buffer is in the range of 8.0 - 8.5. Preferentially use the supplied *GENESpin* Elution Buffer.
- Check working area and pipettes and clean, if necessary.
- Centrifuge at speed indicated in the protocol. Higher speed and prolonged vortexing can lead to shearing of the DNA.
- Repeat washing step with *GENESpin* Wash Buffer 1.

12. Safety instructions

The following components of the *GENESpin* kit contain hazardous contents.

Wear gloves and goggles and follow the safety instructions given in this section.



GHS classification

Component	H-phrases	P-phrases
GENESpin Proteinase K: lyophilized	H315, H319, H334, H335 Danger 	P261, P280, P302+P352, P304+P340 P305+P351+ P338 P342+P311
GENESpin Binding Buffer Guanidine hydrochloride 36-50 %	H302+H332 H315 H319 Warning 	P260, P280 P301+P312 P302+P352 P304+P340 P501

GENESpin

**Kit for DNA isolation
from food and feed**

Cat. no. 5224400605

Component	H-phrases	P-phrases
GENESpin Wash Buffer 1: Guanidine hydrochloride 24-36 %, Ethanol 35-55 %	H226, H302+H332 H315 H319 Warning  	P210 P241 P260 P280 P301+P312 P501

Hazard phrases

H226	Flammable liquid and vapour.
H302+H332	Harmful if swallowed or if inhaled.
H315	Causes skin irritation.
H319	Causes serious eye irritation.
H334	May cause allergy or asthma symptoms or breathing difficulties if inhaled.
H335	May cause respiratory irritation.

Precaution phrases

- | | |
|----------|--|
| P210 | Keep away from heat, sparks, open flames, hot surfaces. No smoking. |
| P241 | Use explosion-proof electrical, ventilating, lighting equipment. |
| P260 | Do not breathe spray, vapours. |
| P261 | Avoid breathing dust. |
| P280 | Wear eye protection, protective gloves, protective clothing. |
| P301+312 | IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. |
| P302+352 | IF ON SKIN: Wash with plenty of water, soap. |
| P304+340 | IF INHALED: Remove person to fresh air and keep comfortable for breathing. |

GENESpin
Kit for DNA isolation
from food and feed
Cat. no. 5224400605

P305+351+338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P342+311	If experiencing respiratory symptoms: Call a POISON CENTER or doctor/physician.
P501	Dispose of contents/container to an approved waste disposal plan.

For further information see Material Safety Data Sheet.

13. Product Warranties and Satisfaction Guarantee

Eurofins GeneScan Technologies GmbH (“GeneScan”) warrants the products manufactured by it will be free of defects in materials and workmanship when used in accordance with the applicable instructions for a period equal to or shorter of one year from the date of shipment of the product(s) or the expiration date marked on the product packaging under the storage conditions, recommended in the instructions and/or on the package. Application protocols published by GeneScan are intended to be only guidelines for the buyers of the products. Each buyer is expected to validate the applicability of each application protocol to his individual application. GeneScan makes no other warranty, expressed or implied. There is no warranty of merchantability or fitness for a particular purpose.

GeneScan’s sole obligation with the respect to the foregoing warranties shall be, at its option, to either replace or to refund the purchase price of the product(s) or part thereof that proves

GENESpin

Kit for DNA isolation from food and feed

Cat. no. 5224400605

defective in materials or workmanship within the warranty period, provided the customer notifies GeneScan promptly of any such defect. GeneScan shall not be liable for any direct, indirect or consequential damages resulting from economic loss or property damages sustained by buyer or any customer from the use of the product(s). A copy of Eurofins GeneScan Technologies GmbH terms and conditions can be obtained on request, and is also provided with our product/price lists.

14. IMPORTANT NOTES

Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

15. TECHNICAL SERVICE

If you have any questions or experience any difficulties regarding this kit or Eurofins GeneScan Technologies products in general, please do not hesitate to contact us.

Eurofins GeneScan Technologies customers are also a major source of information regarding advanced or specialised use of our products. This information is helpful to other scientists as well as to the researchers at Eurofins GeneScan Technologies. We therefore encourage you to contact us if you have any suggestions concerning product performance or new applications and techniques.

For technical assistance and more information please contact the Eurofins GeneScan Technologies customer service department or your local distributor.

GENESpin

**Kit for DNA isolation
from food and feed**

Cat. no. 5224400605

Eurofins GeneScan Technologies GmbH

Engesser Str. 4, 79108 Freiburg, Germany

Phone: + 49-(0)761-5038-200

Fax: + 49-(0)761-5038-111

kits@eurofins.com

www.eurofins.com/kits

GENESpin

Short Protocol

Homogenize	Homogenize 200 mg material	
Lyse cells	+ 550 μ L Lysis Buffer (65°C) + 10 μ L Proteinase K 65°C, 30 min > 10,000 x g, 10 min	
Adjust DNA binding conditions	Clear supernatant (1 volume) + 1 vol. Binding Buffer + 1 vol. ethanol	
Bind DNA	Load 700 μ L sample stepwise (max. loading capacity 750 μ L) 11,000 x g, 1 min	
Wash and dry	1 st wash	+ 400 μ L Wash Buffer 1 11,000 x g, 1 min
	2 nd wash	+ 700 μ L Wash Buffer 2 11,000 x g, 1 min
	3 rd wash	+ 200 μ L Wash Buffer 2 11,000 x g, 2 min
Elute DNA	+ 100 μ L Elution Buffer (70°C) RT, 5 min 11,000 x g, 1 min	