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Aflatoxin M1 supplementary standards

Additional Aflatoxin M1 standard solutions (Cat.nr. HU0040112)

The **Aflatoxin M1 supplementary standards** are two standards solutions containing known amounts of Aflatoxins M1 to be used in combination with **I'screen AFLA M1 (Cat. nr. HU0040002)** or **I'screen AFLA M1 milk (Cat. nr. HU0040001/HU0040021/HU0040041)** kits, in order to step up the number of points in the calibration curve near the European regulatory level for milk at 50 ng/l.

1. PROVIDED REAGENTS

Std A-35 ng/l: 1 plastic vial containing 1.5 ml of a 35 ng/l Aflatoxin M₁ standard solution.

Std B- 65 ng/l: 1 plastic vial containing 1.5 ml of a 65 ng/l Aflatoxin M₁ standard solution.

2. MATERIALS REQUIRED BUT NOT PROVIDED

I'screen AFLA M1 (Cat. nr. HU0040002) or **I'screen AFLA M1 milk (Cat. nr. HU0040001/ HU0040021/HU0040041)** ELISA Kit

3. WARNING AND PRECAUTIONS FOR THE USERS

- This product is for *in vitro* diagnostic use only.
- The solutions contain reagents that may be identified as dangerous substance by the Regulation (EC) N° 1272/2008. For *Material Safety Data Sheet* please contact the technical assistance.
- Handle the reagents with caution, avoiding contact with skin, eyes and mucous membranes.

4. HANDLING AND STORAGE INSTRUCTIONS

- The product is related to each ELISA kit lot. **ATTENTION:** Do not use these calibrators with a different lot of ELISA kit.
- Store the product at +2/+8 °C and never freeze.
- Bring all reagents to room temperature before use (2 hours)
- Do not use components after the expiration date.
- Do not intermix components from different kit lots.
- Do not use photocopies of the instruction booklet. Follow the original instruction booklet that is included with the product.

5. WORKING SOLUTIONS PREPARATION

The solutions are ready to use (**shake gently prior to use**).

6. PROCEDURE FOR USE

Use the product when running the related ELISA kit, following the procedure as reported in chapter 8 "Assay procedure" of the instruction for use of the ELISA kit.


- At point 2 of the procedure, dispense 100 ul of the standards solutions provided within the kit, numbered from 1 to 7, and of the supplementary standard solution A and B. **ATTENTION:** Dispense the Aflatoxin M1 standard solution in the following order:
std 1-2-3-4-A-5-B-6-7
- Proceed as described in chapter 8, without any further modification.

7. CALCULATION OF RESULTS

Calculate the results as reported in chapter 9 of the ELISA kit instruction for use.

Please note: For results calculation, Excel spreadsheets are available. Please contact the technical assistance.

8. CALIBRATION CURVE EXAMPLE

9. EVALUATION OF RESULTS

After processing the results, it is necessary to verify the assay performance. The verification is performed by comparison of obtained data with kit specifications reported in chapter 12 of the instruction for use of **I'screen AFLA M1** or **I'screen AFLA M1 milk** kits.

If the values are outside the specifications given, then the results of the test are not assured, therefore the aflatoxin M1 concentration levels in the samples may not be valid.

In these cases it is advised to check the expiry date of the kit, the wavelength at which the reading was performed, as well as the procedure followed.

If operation errors are not identified as cause, contact our technical assistance.

WARNING: Kit replacement will only be possible in case of return. The kit must be stored in its integral version at +2/+8°C.

10. LIABILITY

Samples evaluated as positive using the kit have to be re-tested with a confirmation method.

Eurofins Technologies Hungary shall not be liable for any damages to the customer caused by the improper use of the kit and for any action undertaken as a consequence of results.

Eurofins Technologies Hungary shall not be liable for the unsafe use of the kit out of the current European safety regulations.