INGEZIM GLUTEN HIDROLIZADO

R.30.GLH.K2

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
INGEZIM GLUTEN HIDROLIZADO

R.30.GLH.K2

Ingezim Gluten Hidrolizado is a direct immunoassay for identify and quantify hydrolyzed prolamins. It uses the R5 Monoclonal Antibody. **Limit of Detection: 0,25 ppm.**

TECHNICAL BASIS OF THE KIT


1. Chemically activated 48 wells plate.
2. If the sample contains the antigen, the prolamins extracted from the sample are fixed to the plate wells, even if they are hydrolyzed.
3. After a washing step, when the R5 MAb-PO is added, it will bind to the fixed prolamins of the sample.
4. Binding is detected by a colorimetric reaction after the addition of the substrate. This reaction is directly proportional to the amount of gluten in the well.

APPLICATION

The assay has been designed for the detection of Gluten in hydrolyzed samples (beer, whiskey, syrup...)

SENSITIVITY

The assay is able to detect **0.25ppm** of GLUTEN, by using a linear regression curve. It leads to more precise results.


SPECIFICITY

The specificity of the MAb used in the assay was determined by the immunoblotting technique. The MAb is able to detect specifically the prolamins of gluten present in WHEAT, RYE and BARLEY (gliadin, secalin and hordein respectively).

COMPARISON WITH THE DAS ELISA

When comparing several beer samples, the recovery it's too much higher by using the DIRECT ELISA than using the DAS (Double Antibody Sandwich) ELISA.

COMPOSITION OF THE KIT

- Microtritation plate of 48 well
- European Gliadin Standard (100X concentrated)
- Vial with Anti prolamins Mab Conjugate
- Bottle with Washing Solution
- Bottle with diluent
- Bottle with substrate (TMB)
- Bottle with stop solution
- Bottle with Extraction Solution

PRODUCT MANUFACTURED BY INGENASA

SHELF LIFE: 10 months
Stored at 4°C

02/02/2018