

INgezim M.HYO COMPAC

R.11.MHYO.K3

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
INgezim M.HYO COMPAC

R.11.MHYO.K3

Ingezim M. HYO Compac is based on a blocking ELISA assay technique, which uses a monoclonal antibody (MAb) specific to *M. hyopneumoniae*.

TECHNICAL BASIS OF THE KIT

1. Plates are coated with *M. hyopneumoniae* antigen. (Serum samples are added and incubated).
2. If the samples contain antibodies specific to *M. hyopneumoniae*, they will bind to the antigen.
3. When a MAb-PO specific for *M. hyopneumoniae* is added, only if there are no antibodies in the sample blocking the antigen (negative animals), it will bind to the protein. In case the sample contains antibodies blocking the antigen (infected animals), the conjugate will not be able to bind to it. The binding is detected by the development of a colorimetric reaction after the addition of the substrate..

APPLICATION

Detection and /or titration of specific antibodies of *M. hyopneumoniae* in porcine sera samples.

RESULT INTERPRETATION

There are two Cut Off Values for the results interpretation: positive and negative. Samples with an OD lower than positive Cut Off Value must be considered **Positive** and samples with OD higher than negative Cut Off Value must be considered **Negative**. Samples with OD value between both Cuts Offs must be considered **Doubtful**.

VALIDATION


SENSITIVITY AND SPECIFICITY

1. Field trials (swine sera from different sources and with different number of births)

The serological analysis for this disease is complicated due to the great variability in the individual response to the infection and to the fact that there are different vaccination plans and products commercialized , which do not allow to predict the results. Furthermore, according to many authors, most of the available monodose vaccines are unable to produce a major seroconversion in immunized animals until several weeks after inoculation (8-12 weeks). With these considerations and taking into account the high prevalence of *M. hyopneumoniae* in intensive porcine farms, it appears to be a high homology of results within the different batches of animals analyzed and that it is possible to find clear differences between groups when INGEZIM M.HYO Compac is used.

2. Experimental studies.

Using INGEZIM M.HYO Compac, samples from piglets taken at different dates after vaccination and vaccinated with different formulas (experimental and commercial ones) were analyzed. The following figure shows the OD values obtained in an experimental study of vaccination and challenge (day 30) made with one of the vaccines used. Different extractions were taken from animals from different sources.

— Negative cut off — Positive cut off

PRODUCT MANUFACTURED BY INGENASA

COMPOSITION OF THE KIT

- Microtitration plates of 96 well divided into 12 strips of 8 wells each.
- Vials with Positive Control serum ready to use
- Vials with Negative Control serum ready to use
- Vials with Peroxidase Conjugate ready to use
- Bottle with Washing Solution 25x concentrated
- Bottle with Diluent ready to use
- Bottle with Substrate (TMB)
- Bottle with Stop Solution

IT-73840 ISO 14001:2015 ISO 9001:2015
IT-73780 9191.INGE 9175.ING2

SHELF LIFE: **12 months**
Stored 2°C-8°C

Ed.020217