

INgezim ADV gB

R.11.GB.K3

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

INgezim ADV gB

R.11.GB.K3

INgezim ADV gB is an enzymatic assay based on a blocking ELISA assay technique, which uses a monoclonal antibody (MAb) specific to Aujeszky Disease Virus (ADV) gB protein.

**TECHNICAL BASIS OF THE KIT**

1. Plates are coated with inactivated ADV antigen. Serum samples are added and incubated.
2. If the samples contain specific antibodies to ADV, they will bind the antigen.
3. When a MAb-PO specific to ADV gB protein is added, it will bind to the protein only if there are no antibodies in the sample blocking the antigen (negative animals). In case there are antibodies blocking the antigen (infected animals), the conjugate will not be able to bind it. The binding is detected by the development of a colorimetric reaction after the addition of the substrate.

APPLICATION

Detection of specific antibodies to ADV's gB protein which is present in both, vaccinated and infected animals

RESULT INTERPRETATION

One Cut Off is established for the results interpretation. Samples with an OD lower than the Cut Off value are Positive and samples with OD higher than the Cut Off value are Negative.

VALIDATION**SENSITIVITY AND SPECIFICITY****1. Using reference sera from the O.I.E.**

Aujeszky's disease ADV-1 equivalent serum is detected, as it is required by the UE Directive. A set of reference sera (G and Q) catalogued as positive or negative by an O.I.E. Reference Laboratory (Maisons Alfort, France) were analyzed. The correspondence with the expected results was **100%**.

2. Sensitivity using experimentally vaccinated sera

Experimentally vaccinated animals were used and the antibodies detection profile was obtained. The assay can detect antibodies specific to gB from day 6 post vaccination.

3. Diagnostic sensitivity and specificity

1740 Spanish field sera previously catalogued by other assays were analyzed. The assay shows 99.5% sensitivity and 99% specificity.

4. Correspondence with the reference technique SN (seroneutralization). Experimental studies.

210 sera from 18 experimentally vaccinated animals (with a deleted vaccine) and 8 negative controls were analyzed. The samples were taken at different days starting at day 0, during 4 months. The results obtained indicated that the correspondence between INgezim ADV gB and Seroneutralization was **95%**.

5. External validation

The product has been included in the Serology Ring Test (ELISA) for Aujeszky disease report of 2011 organized by ANSES. In this report it is described that INgezim ADV gB is able to detect dilution 1/64 of serum S96AUJ02 which is equivalent to dilution 1/8 of OIE Reference Serum ADV-1 (as it is declared in the same report). INgezim ADV gB obtained a correspondence of 100% with the expected results in this study.

COMPOSITION OF THE KIT

- Microtitration plates of 96 wells.
- Vials with Positive Control
- Vials with Negative Control
- Vials with Peroxidase Conjugate
- Bottle with Washing Solution
- Bottle with serum diluent
- Bottle with conjugate diluent
- Bottle with stop solution
- Bottle with substrate (TMB) ready to use

PRODUCT MANUFACTURED BY INGENASA
REGISTRATION NUMBER 809 RD

SHELF LIFE: **18 months**
Stored 2°C-8°C

Ed. 020217