
INgezim CIRCO IgG

R.11.PCV.K1

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

Correlación lineal ($R^2=0,71$);
Correlación polinómica ($R^2=0,91$)

REFERENCIAS

(A)Segalés et al. 6th International Symposium on Emerging and Re-emerging Pig Diseases Barcelona, Spain, 2011); (B)Pileri et al. 22nd International Pig Veterinary Society Congress. Jeju, Korea (2012); (B)Segalés et al. AVEDILA 2013, Madrid Spain

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

Ed.020217

INgezim CIRCO IgG R.11.PCV.K1

INgezim PCV indirect is based on an indirect ELISA technique, which uses a monoclonal antibody (MAb) specific to porcine immunoglobulines (Igs), and a recombinant antigen (PCV2 VP2 protein).

TECHNICAL BASIS OF THE KIT

- Plates are coated with PCV2 antigen (recombinant protein VP2). Serum samples are added and incubated.
- If the samples contain specific antibodies to PCV, they will bind to the antigen.
- When a MAb-PO specific of porcine immunoglobulines is added, it will bind to the Igs bound to the antigen. The binding is detected by the development of a colorimetric reaction after the addition of the substrate.

APLICACION

Detection and /or titration of specific antibodies to PCV in porcine sera samples

INTERPRETATION OF THE RESULTS

There are two Cut Off values for the results interpretation. Samples with an OD higher than the positive Cut Off are **positive** and samples with an OD lower, are **negative**.


The titer of the sample is the last dilution showing an OD higher than the positive cut off.

One well titration: It is possible to obtain an approximated titre of the sera based on its OD value, by using a formula which corresponds with an exponential curve.

VALIDATION OF THE ASSAY

1. Experimental Study 1

A set of sera from three experimentally infected piglets previously vaccinated were used. The samples consisted of different extractions made at days 0, 16 and 46 post vaccination and 14 and 22 post infection. Piglets 40 and 70 showed presence of maternal antibodies at day 0.


2. Experimental study 2.

Sera from 12 experimentally infected animals and 8 non infected animals were analyzed. Samples consisted of extractions taken at days 14 and 21 post infection. The results indicated that the assay is able to detect antibodies at least at day 21 post infection.



3. Studies with piglets.

Studies with piglets from vaccinated sows were performed. Extractions between 3 and 24 weeks old were taken. The results were compared with the ones obtained using INGEZIM CIRCOVIRUS IgG/IgM.

4. Correlation with IPMA

Study A: This study was made with 645 sera of animals of different ages. Study B: This study was made with 1248 sera of animals of different ages. Log₂ PCV2 IPMA titers correspondence with ELISA OD (A) or S/P (B) values.

REFERENCIAS

(A)Segalés et al. 6th International Symposium on Emerging and Re-emerging Pig Diseases Barcelona, Spain, 2011), (B)Pileri et al. 22nd International Pig Veterinary Congress. Jeju, Korea (2012); (B)Segalés et al. AVEDILA 2013, Madrid Spain

COMPOSITION OF THE KIT

- Microtitration plates of 96 wells divided into 12 strips of 8 wells each.
- Vials with Positive Control
- Vials with Negative Control
- Vials with Peroxidase Conjugate
- Bottle with Washing Solution concentrated
- Bottle with Diluent
- Bottle with Substrate TMB
- Bottle with Stop Solution

PRODUCT MANUFACTURED BY INGENASA

SHELF LIFE: **18 months**
Stored at 2°C-8°C

Ed.020217