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

NÚMERO DE REGISTRO 526 RD

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
CADUCIDAD: **18 meses**
Conservado a 2°C-8°C

INgezim IBR is an immunoenzymatic assay based on an indirect ELISA technique, which uses a monoclonal antibody (MAb) specific of bovine immunoglobulins and an inactivated antigen.

TECHNICAL BASIS OF THE KIT

- Plates are coated with inactivated IBR antigen. Samples are added to the wells and incubated.
- If the sample contains antibodies specific to IBR, they will bind to the antigen.
- When a MAb-PO specific for bovine immunoglobulins is added, it will bind to the Igs of the sample previously bound to the antigen. This binding is detected by the development of a colorimetric reaction after the substrate addition.

APPLICATION

Detection of specific antibodies to Infectious Bovine Rhinotracheitis Virus (IBR) in bovine serum, plasma, milk and whey samples.

INTERPRETATION OF THE RESULTS

Two cut off are used for the results interpretation. Samples will be considered **Positive** if their OD value is higher than the positive cut off; **Negative** if their OD value is lower than the negative cut off and **Doubtful** if the OD values are between both cut offs.

VALIDATION

1. Use of O.I.E. reference sera

In order to determine the validity of the assay, 3 O.I.E. Reference Sera were used: EU1 (strong positive), EU2 (weak positive) and EU3 (negative). The correspondence with the expected results was **100%**.

2. Use of previously clasified sera.

Field sera from different sources, positive and negative by blocking ELISA (INgezim® IBR compac 2.0) and indirect ELISA (Svanovir® IBR Ab) were analyzed. Table below shows correspondence between assays:

	INgezim® IBR Compac 2.0	Svanovir® IBR Ab
INgezim® IBR 2.0	93,7%	94,3%
INgezim® IBR Compac 2.0		93,7%

3. Use of FLI reference sera (Friedrich-Loeffler-Institut, Riems, Alemania)

6 ser aware used: Weak positive by seroneutralization (R1C, R2C y R3D) and negtaive by SN (R31B y R32C). Correspondence with the expected results was **100%**.

4. Use of milk samples

90 Field milk samples from different sources were analyzed. Samples consisted on positive and negative ones by a blocking ELISA (INgezim® IBR compact 2.0) and indirect ELISA (Svanovir® IBR Ab). Table below shows the correspondence between assays:

	INgezim® IBR Compac 2.0	Svanovir® IBR Ab
INgezim® IBR 2.0	94,4%	95,5%
INgezim® IBR Compac 2.0		95,5%

COMPOSITION OF THE KIT

- Microtitration plates of 96 wells coated with IBR antigen.
- Vial with Positive Control.
- Vial with Negative Control.
- Vial with Conjugate for milk.
- Vial with Conjugate for serum
- Bottles with Washing Solution concentrated.
- Bottles with diluent for milk
- Bottles with diluent for serum
- Bottle with substrate (TMB) ready to use.
- Bottle with stop solution.

SPANISH REGISTRATION N 526-RD
PRODUCT MANUFACTURED BY INGENASA

Shelf life: **18 months**
Store at 2°C-8°C