



Validation Report
SENSISpec Pistachio ELISA
(Cat. No. HU0030021)

Version 4

CONTENT

1. SCOPE	2
2. PRECISION	2
2.1. REPEATABILITY	2
2.2. REPRODUCIBILITY	3
2.3. REPRODUCIBILITY OF EXTRACTION	4
3. ANALYTICAL SENSITIVITY	5
4. RECOVERY	7
4.1. FOOD SAMPLES	7
4.1.1. HEAT-TREATED MATRICES	11
4.2. RINSE WATER /CIP	11
4.3. SWABS	13
5. LINEARITY	13
6. CROSS-REACTIVITY	15
7. ROBUSTNESS	16
7.1. VARIATION OF EXTRACTION TEMPERATURE	16
7.2. VARIATION OF EXTRACTION TIME	17
7.3. DRIFT	17



Validation Report

SENSISpec Pistachio ELISA

(Cat. No. HU0030021)

Version 4

1. Scope

The **SENSISpec Pistachio ELISA** is designed for the sensitive determination of pistachio in food, based on a polyclonal antibody. The present report describes the validation process and summarizes its results.

2. Precision

2.1. Repeatability

Repeatability (intra-assay variation) was determined by testing three controls of various concentration levels in 24fold replicates, either in sample diluent and cookie matrix. After identification of possible outliers the following results were obtained.

Table 1: Repeatability of the SENSISpec Pistachio ELISA in sample diluent

Replicate	1 ppm	4 ppm	20 ppm
1	0.98	3.90	18.7
2	0.96	4.08	22.7
3	0.97	3.98	22.6
4	1.03	4.38	22.5
5	1.02	4.10	19.4
6	1.01	4.04	19.1
7	0.97	3.86	19.1
8	0.86	4.14	19.4
9	0.87	3.95	21.8
10	0.84	4.33	22.5
11	0.86	4.06	18.8
12	0.90	4.39	19.6
13	0.87	4.06	19.2
14	0.94	4.20	20.3
15	0.96	4.00	23.9
16	0.96	4.27	22.2
17	1.00	3.99	19.2
18	0.97	4.18	18.8
19	0.89	4.21	20.4
20	0.96	4.08	23.1
21	0.94	3.89	21.5
22	0.92	4.04	18.7
23	0.92	3.91	24.8
24	0.93	4.11	17.0
Mean	0.94	4.09	20.6
SD	0.05	0.15	2.02
CV [%]	5.7	3.7	9.8

The coefficient of variation in sample diluent is ranging from 3.7% to 9.8% depending on the concentration.



Validation Report

SENSISpec Pistachio ELISA

(Cat. No. HU0030021)

Version 4

Table 2: Repeatability of the SENSISpec Pistachio ELISA in cookie matrix

Replicate	1 ppm	4 ppm	20 ppm	
1	1.01	4.15	20.6	
2	1.01	4.39	21.6	
3	1.06	4.66	23.6	
4	1.16	4.43	21.2	
5	1.13	4.64	19.7	
6	1.07	4.53	20.3	
7	1.00	4.38	19.2	
8	1.02	4.27	20.3	
9	1.00	4.48	19.7	
10	1.01	4.71	20.3	
11	1.05	4.77	24.0	
12	1.04	4.48	24.3	
13	1.01	4.31	22.7	
14	0.96	4.42	18.0	
15	1.05	4.10	20.6	
16	1.01	4.16	19.5	
17	1.02	4.27	20.6	
18	0.99	4.40	23.6	
19	1.04	4.57	24.9	
20	1.01	4.27	22.3	
21	1.03	4.20	19.0	
22	0.97	4.28	20.4	
23	1.01	4.46	17.6	
24	0.97	4.44	19.4	
Mean	1.03	4.41	21.0	
SD	0.05	0.18	2.00	Mean
CV [%]	4.5	4.1	9.5	6.0

The coefficient of variation in cookie matrix is ranging from 4.1% to 9.5% depending on the concentration.

2.2. Reproducibility

Reproducibility (inter-assay variation) was determined by testing three controls of various concentration levels in four different test runs of the same kit lot, either in sample diluent and cookie matrix.



Validation Report

SENSISpec Pistachio ELISA

(Cat. No. HU0030021)

Version 4

Table 3: Reproducibility of the SENSISpec Pistachio ELISA in sample diluent

Assay No.	1 ppm	4 ppm	20 ppm	
1	0.84	3.79	17.9	
2	0.97	4.20	20.8	
3	0.86	4.16	22.2	
4	1.03	4.05	19.5	
Mean	0.93	4.05	20.1	
SD	0.09	0.18	1.84	Mean
CV [%]	9.7	4.5	9.2	7.8

The coefficient of variation in sample diluent is ranging from 4.5% to 9.7% depending on the concentration.

Table 4: Reproducibility of the SENSISpec Pistachio ELISA in cookie matrix

Assay No.	Level 1	Level 2	Level 3	
1	0.79	3.97	23.3	
2	0.82	3.72	22.0	
3	0.80	3.27	19.3	
4	1.01 ¹⁾	3.63	19.2	
Mean	0.80	3.65	20.9	
SD	0.01	0.29	2.03	Mean
CV [%]	1.7	7.9	9.7	6.4

¹⁾ identified as outlier by 1.5 IQR testing; not considered in the calculation

The coefficient of variation in cookie matrix is ranging from 1.7% to 9.7% depending on the concentration.

2.3. Reproducibility of Extraction

For determining the reproducibility of the extraction process pistachio proficiency test material (DLA 07/2019, sample A) was applied. The homogeneity of the material was assured by the proficiency test provider by analyzing 20 different samples of the material with 3 different test kits, respectively. According PT report the coefficient of variation was ranging from 3.8% to 5.5%.

The material was extracted ten times. The coefficient of variation was calculated by applying the mean results of every extract.

Table 5: Reproducibility of extraction of the SENSISpec Pistachio ELISA

Extraction No.	Mean Result [ppm]
1	36.4
2	38.4
3	39.0
4	38.9



Validation Report

SENSISpec Pistachio ELISA

(Cat. No. HU0030021)

Version 4

5	34.4
6	39.6
7	42.6
8	38.1
9	37.0
10	39.9
Mean	38.4
SD	2.21
CV%	5.8

The coefficient of variation is 5.8% for the inter-extraction reproducibility trial.

3. Analytical Sensitivity

For determination of the analytical sensitivity (LOD), sample diluent and various pistachio free representative sample matrices were assayed in 24fold replicates. After identification of possible outliers the OD mean and standard deviation were calculated. The corresponding concentration of the OD mean + 3x standard deviation was defined as limit of detection. This results in limits of detection according to the following table:

Table 6a: Matrix-dependent and matrix-independent analytical sensitivity of the Immunolab Pistachio ELISA (complete data)

Replicate	Sample diluent [OD]	Cookie matrix [OD]	Sweets matrix [OD]	Ice cream matrix [OD]	Chocolate ¹⁾ matrix [OD]	Spice ¹⁾ matrix [OD]	Food supplement matrix [OD]
1	0.059	0.114	0.063	0.072	0.116	0.199	0.118
2	0.058	0.125	0.064	0.081	0.118	0.220	0.118
3	0.058	0.125	0.065	0.074	0.120	0.240	0.114
4	0.059	0.133	0.065	0.088	0.130	0.220	0.122
5	0.059	0.134	0.067	0.077	0.126	0.254	0.126
6	0.056	0.131	0.064	0.076	0.120	0.233	0.121
7	0.051	0.122	0.074 ²⁾	0.068	0.112	0.221	0.121
8	0.053	0.115	0.066	0.067	0.102	0.205	0.129
9	0.060	0.130	0.066	0.069	0.112	0.250	0.115
10	0.062	0.133	0.060	0.086	0.122	0.205	0.119
11	0.059	0.124	0.068	0.076	0.118	0.222	0.121
12	0.059	0.128	0.070	0.083	0.117	0.226	0.122
13	0.061	0.132	0.066	0.079	0.124	0.247	0.121
14	0.054	0.133	0.065	0.080	0.119	0.247	0.126
15	0.049	0.109	0.062	0.067	0.111	0.230	0.128
16	0.062	0.135	0.067	0.075	0.109	0.224	0.129
17	0.061	0.134	0.069	0.070	0.112	0.203	0.111
18	0.057	0.134	0.063	0.076	0.119	0.204	0.117
19	0.061	0.128	0.064	0.079	0.117	0.201	0.116
20	0.059	0.134	0.068	0.075	0.121	0.219	0.116
21	0.048	0.126	0.063	0.071	0.118	0.228	0.125
22	0.050	0.122	0.065	0.074	0.118	0.233	0.121
23	0.050	0.113	0.062	0.076	0.112	0.214	0.123



Validation Report

SENSISpec Pistachio ELISA

(Cat. No. HU0030021)

Version 4

Replicate	Sample diluent [OD]	Cookie matrix [OD]	Sweets matrix [OD]	Ice cream matrix [OD]	Chocolate ¹⁾ matrix [OD]	Spice ¹⁾ matrix [OD]	Food supplement matrix [OD]
24	0.055	0.121	0.063	0.078	0.111	0.204	0.117
Mean	0.057	0.126	0.065	0.076	0.117	0.223	0.121
SD	0.004	0.008	0.002	0.006	0.006	0.017	0.005
LOD [ppm]	0.07	0.42	0.03	0.16	0.28	0.84	0.27
Corresponding Calibration Curve data [OD]							
0 ppm	0.057	0.058	0.065	0.060	0.065	0.067	0.070
1 ppm	0.277	0.277	0.317	0.270	0.317	0.312	0.313
4 ppm	0.721	0.721	0.825	0.711	0.825	0.791	0.788
10 ppm	1.198	1.198	1.313	1.193	1.313	1.307	1.272
40 ppm	1.836	1.836	2.103	1.794	2.103	2.100	1.929

¹⁾ addition of 1% PVP to sample extraction buffer prior extraction

²⁾ identified as outlier by 1.5 IQR testing; not considered in the calculation

Note: In some cases, various samples of the same kind (e.g. different cookies) were tested. The table above represents the matrix with the highest LOD. Anyway, the table below summarizes all results.

Table 6b: Matrix-dependent and matrix-independent analytical sensitivity of the Immunolab Pistachio ELISA (summary)

Matrix Type	Variant	LOD
Sample diluent	N/A	0.07 ppm
Cookies	"Spoon Biscuit"	0.06 ppm
	"TUC" cookie	0.06 ppm
	"Lemon cookie"	0.42 ppm
Ice-cream	Type Vanilla	0.02 ppm
	Type Lab sample	0.00 ppm
	Type Strawberry	0.16 ppm
Sweets	Vitamin drop	0.04 ppm
	Chewy candy	0.00 ppm
Chocolate ¹⁾	Chocolate bar	0.02 ppm
	Type 1, 59% cocoa	0.15 ppm
	Type 2, 72% cocoa	0.28 ppm
Spices ¹⁾	Curcuma	0.07 ppm
	Paprika	0.84 ppm
	Chili	0.73 ppm
	Cumin	0.20 ppm
Food supplement	"Almased"	0.07 ppm
	"Wellmix Classic"	0.21 ppm
	"Wellmix Vanilla"	0.11 ppm

¹⁾ addition of 1% PVP to sample extraction buffer prior extraction



Validation Report

SENSISpec Pistachio ELISA

(Cat. No. HU0030021)

Version 4

With respect to the sample matrix limits of detection vary from 0.00 to 0.84 ppm. Note that the derived limits of detection are strictly dependent on the coefficient of variation and may thus vary in every individual test.

The lowest positive standard (1 ppm) was defined as limit of quantification to assure that all uncontaminated matrices result in concentrations lower than this value.

4. Recovery

4.1. Food Samples

For recovery experiments different sample matrices were spiked with pistachio extract to obtain various final concentrations after performing all sample pre-treatment steps. The result of the un-spiked matrix was subtracted before evaluation. Tested samples and calculated results were as follows.

Table 7: Recovery of various samples tested with the SENSISpec Pistachio ELISA

Cookie, type "Spoon biscuit", Zero level: 0.00 ppm

Target Value	Actual Concentration	Recovery [%]
1 ppm	0.88	88
4 ppm	3.41	85
10 ppm	8.17	82
Mean		85

Cookie, type "TUC", Zero level: 0.01 ppm

Target Value	Actual Concentration	Recovery [%]
1 ppm	0.95	95
4 ppm	3.84	96
10 ppm	9.43	94
Mean		95

Cookie, type "Lemon Cookie", Zero level: 0.18 ppm

Target Value	Actual Concentration	Recovery [%]
1 ppm	1.00	100
4 ppm	3.66	92
10 ppm	7.90	79
Mean		90



Validation Report

SENSISpec Pistachio ELISA

(Cat. No. HU0030021)

Version 4

Sweets, type "Vitamin drop", Zero level: 0.00 ppm¹⁾

Target Value	Actual Concentration	Recovery [%]
1 ppm	0.85	85
4 ppm	3.40	85
10 ppm	7.61	76
Mean		82

Sweets, type "Chewy candy", Zero level: 0.00 ppm¹⁾

Target Value	Actual Concentration	Recovery [%]
1 ppm	0.76	76
4 ppm	3.09	77
10 ppm	7.01	70
Mean		74

Ice-cream, type Lab sample, Zero level: 0.00 ppm

Target Value	Actual Concentration	Recovery [%]
1 ppm	0.69	69
4 ppm	3.07	77
10 ppm	9.07	91
Mean		79

Ice-cream, type "Vanilla", Zero level: 0.00 ppm

Target Value	Actual Concentration	Recovery [%]
1 ppm	1.05	105
4 ppm	3.42	86
10 ppm	9.65	97
Mean		96

Ice-cream, type "Strawberry", Zero level: 0.00 ppm

Target Value	Actual Concentration	Recovery [%]
1 ppm	0.85	85
4 ppm	3.27	82
10 ppm	8.00	80
Mean		82



Validation Report

SENSISpec Pistachio ELISA

(Cat. No. HU0030021)

Version 4

Chocolate, type Chocolate bar, Zero level: 0.03 ppm¹⁾

Target Value	Actual Concentration	Recovery [%]
1 ppm	0.81	81
4 ppm	3.32	83
10 ppm	7.55	76
Mean		80

Chocolate, type 59% Cocoa, Zero level: 0.12 ppm¹⁾

Target Value	Actual Concentration	Recovery [%]
1 ppm	0.83	83
4 ppm	3.52	88
10 ppm	8.26	83
Mean		84

Chocolate, type 72% Cocoa, Zero level: 0.11 ppm¹⁾

Target Value	Actual Concentration	Recovery [%]
1 ppm	0.79	79
4 ppm	3.31	83
10 ppm	7.01	70
Mean		77

Spices, Curcuma, Zero level: 0.00 ppm¹⁾

Target Value	Actual Concentration	Recovery [%]
1 ppm	0.73	73
4 ppm	3.06	77
20 ppm	6.91	69
Mean		73

Spices, Cumin, Zero level: 0.18 ppm¹⁾

Target Value	Actual Concentration	Recovery [%]
1 ppm	0.70	70
4 ppm	2.81	70
10 ppm	10.69	107
Mean		82



Validation Report

SENSISpec Pistachio ELISA

(Cat. No. HU0030021)

Version 4

Spices, Chili, Zero level: 0.57 ppm¹⁾

Target Value	Actual Concentration	Recovery [%]
1 ppm	0.90	90
4 ppm	3.42	86
10 ppm	7.95	80
Mean		85

Spices, Paprika, Zero level: 0.00 ppm¹⁾

Target Value	Actual Concentration	Recovery [%]
1 ppm	0.86	86
4 ppm	3.44	86
10 ppm	6.77	68
Mean		80

Food supplement, type "Almased", Zero level: 0.17 ppm

Target Value	Actual Concentration	Recovery [%]
1 ppm	0.67	67
4 ppm	3.06	76
10 ppm	7.14	71
Mean		72

Food supplement, type "Wellmix Classic", Zero level: 0.10 ppm

Target Value	Actual Concentration	Recovery [%]
1 ppm	0.66	66
4 ppm	2.52	63
10 ppm	8.28	83
Mean		70

Food supplement, type "Wellmix Vanilla", Zero level: 0.05 ppm

Target Value	Actual Concentration	Recovery [%]
1 ppm	0.71	71
4 ppm	2.72	68
10 ppm	8.32	83
Mean		74

¹⁾ addition of 1% PVP to sample extraction buffer prior extraction

Mean recoveries are ranging from **70%** to **96%** depending on the sample matrix.



Validation Report

SENSISpec Pistachio ELISA

(Cat. No. HU0030021)

Version 4

4.1.1. Heat-treated Matrices

To exemplarily assess changes in analytical sensitivity and recovery for processed food samples, a cookie matrix was spiked with 10 ppm of pistachio and divided in two parts. One part was baked at 150 °C for 20 min. Processed as well as unprocessed extracts were tested undiluted as well as in various dilutions. Additionally, non-spiked zero matrix was tested. All tests were performed in duplicates.

Table 8: Accuracy (heat-treated matrices) of the SENSISpec Pistachio ELISA

Matrix	Level 1 [0 ppm]	Level 4 [1.25 ppm]	Level 5 [2.5ppm]	Level 6 [5 ppm]	Level 7 [10 ppm]
Cookie / unprocessed	0.075	1.42	2.79	5.92	10.7
Cookie / processed	0.052	0.66	1.39	2.70	5.23

Due to the baking process recovery decreases by ~50%.

4.2. Rinse water /CIP

The sensitivity for rinse water / CIP was determined by diluting pistachio extract in different commonly applied water-based rinse solutions to various concentrations. The pH was adjusted to 8.2 ± 0.5 before running the test. Afterwards the samples were diluted 1:5 in pre-diluted sample extraction buffer. Tested samples and calculated results were as follows.

Table 9: Recovery of the SENSISpec Pistachio ELISA for Rinse Solutions

Pure water, Zero level 0.00 mg/L

Target Value	Actual Concentration	Recovery [%]
0.25 mg/L	0.22	88
0.5 mg/L	0.46	93
1 mg/L	0.93	93
2 mg/L	1.99	99
4 mg/L	3.45	86
	Mean	92



Validation Report

SENSISpec Pistachio ELISA

(Cat. No. HU0030021)

Version 4

0.5% Sodium triphosphate, Zero level 0.00 mg/L

Target Value	Actual Concentration	Recovery [%]
0.25 mg/L	0.23	91
0.5 mg/L	0.35	70
1 mg/L	0.76	76
2 mg/L	1.78	89
4 mg/L	4.26	107
Mean		87

2% NaOH, Zero level 0.00 mg/L

Target Value	Actual Concentration	Recovery [%]
0.25 mg/L	0.20	79
0.5 mg/L	0.48	97
1 mg/L	0.92	92
2 mg/L	1.52	76
4 mg/L	2.59	65
Mean		82

0.5% Nitric acid, Zero level 0.00 mg/L

Target Value	Actual Concentration	Recovery [%]
0.25 mg/L	0.23	92
0.5 mg/L	0.43	86
1 mg/L	0.90	90
2 mg/L	1.89	94
4 mg/L	3.84	96
Mean		92

0.5% Citric acid, Zero level 0.0 mg/L

Target Value	Actual Concentration	Recovery [%]
0.25 mg/L	0.26	106
0.5 mg/L	0.48	97
1 mg/L	0.86	86
2 mg/L	2.06	103
4 mg/L	4.69	117
Mean		102

Mean recoveries for rinse water are ranging from **82%** to **102%** depending on the kind of rinse solution.



Validation Report

SENSISpec Pistachio ELISA

(Cat. No. HU0030021)

Version 4

4.3. Swabs

To monitor the swab performance, surfaces (5 cm * 5 cm) were contaminated with defined amounts of a pistachio (in solution). Each amount was applied in duplicate. The surfaces were allowed to dry completely. The swabs were moistened in 1 mL pre-diluted sample extraction buffer. Then the surfaces were swabbed thoroughly. After swabbing the swabs were redissolved in the initial 1 mL of buffer. After redissolving 100 μ L of the solutions were applied in the ELISA. The mean results were recalculated to the amount of pistachio / area.

Table 10: Recovery of the SENSISpec Pistachio ELISA for Swabs

Swabs, Zero level 0.00 ng / cm²

Target Value [ng / cm ²]	Test Result [ng / cm ²]	CV [%]	Recovery [%]
20	15.9	4.9	80
60	42.8	10	71
200	147	39	73
	Mean	18	75

Recovery is ranging from 71-80%. Mean recovery is 75%. Repeated swabbing of the same surface indicates that varying recoveries are an intrinsic feature of the method which can only be improved by multiple repetitions. Recoveries are affected by solubility of proteins, residual liquid and complete swabbing of the surface area.

Note that surfaces were contaminated by protein extracts and not complete commodities. As the latter are more particulate by nature they might be swabbed more easily. Thus, the approach chosen here reflects a worst-case scenario.

5. Linearity

Linearity was determined by spiking various matrices with pistachio extract and testing subsequent dilutions of the resulting extracts. For calculation of the linearity the highest concentration was defined as reference value (100%) and further dilutions were expressed in percent of this reference after consideration of the dilution factor. Regarding the matrices extracted by addition of PVP, further dilution was applied with standard extraction buffer without PVP.



Validation Report

SENSISpec Pistachio ELISA

(Cat. No. HU0030021)

Version 4

Table 11: Matrix dependent linearity of the SENSISpec Pistachio ELISA

Spoon Cookie

Target Value	Concentration [ppm]	Recovery [%]
40 ppm	38.1	100
20 ppm	19.2	101
10 ppm	9.57	100
5 ppm	5.07	107
2.5 ppm	2.66	112
	Mean	105

Vitamin drop

Target Value	Concentration [ppm]	Recovery [%]
40 ppm	33.6	100
20 ppm	19.2	114
10 ppm	10.2	122
5 ppm	5.02	119
2.5 ppm	2.50	119
	Mean	119

Ice-cream, Type Vanilla

Target Value	Concentration [ppm]	Recovery [%]
40 ppm	31.8	100
20 ppm	18.6	117
10 ppm	7.62	96
5 ppm	4.08	103
2.5 ppm	2.29	115
	Mean	108

Chocolate bar¹⁾

Target Value	Concentration [ppm]	Recovery [%]
40 ppm	35.2	100
20 ppm	21.1	120
10 ppm	10.0	114
5 ppm	5.33	121
2.5 ppm	2.48	113
	Mean	117

Spices, Curcuma¹⁾

Target Value	Concentration [ppm]	Recovery [%]
40 ppm	36.0	100
20 ppm	20.0	111
10 ppm	9.21	102
5 ppm	5.54	123
2.5 ppm	2.81	125
	Mean	116



Validation Report

SENSISpec Pistachio ELISA

(Cat. No. HU0030021)

Version 4

Food Supplement “Almased”

Target Value	Concentration [ppm]	Recovery [%]
40 ppm	31.3	100
20 ppm	15.3	98
10 ppm	7.53	96
5 ppm	3.65	93
2.5 ppm	1.99	102
	Mean	97

¹⁾ addition of 1% PVP to sample extraction buffer prior extraction

For different matrices the mean linearity is ranging from 97% to 119%. The linearity is independent of the specific concentration and may only be affected by the precisions stated in chapter 2.

6. Cross-Reactivity

For the following foods no cross-reactivity (results < LOQ) could be detected:

Table 12: Non-cross-reactive food matrices in the SENSISpec Pistachio ELISA

Raw material	c [ppm]	Raw material	c [ppm]
Adzuki bean	0.16	Goat's milk	0.34
Almond	0.43	Guar gum	0.13
Apricot	0.53	Gum arabic	0.86
Barley	0.17	Hazelnut	0.99
Bean, white	0.36	Kidney bean	0.12
Beef	0.09	Kiwi	0.03
Bovine gelatin	0.06	Lamb	0.03
Brazil nut	0.17	Leek	0.41
Buckwheat	0.48	Lentil	0.67
Cabbage, white	0.68	Lupin	0.17
Caraway	0.06	Macadamia	0.95
Cardamom	0.14	Nutmeg	0.10
Carob gum	0.16	Oats	0.10
Carrot	0.08	Onion	0.40
Cayenne	0.19	Paprika	0.18
Celery	0.08	Pea	0.45
Cherry	0.17	Peach	0.03
Chestnut	0.21	Peanut	0.17
Chia	0.63	Pecan	0.53
Chicken	0.53	Pepper, black	0.43
Chickpea	0.16	Pine seed	0.03
Chili	0.65	Pork	0.14
Cinnamon	0.07	Potato	0.15
Clove	0.39	Prawn	0.09
Cocoa	0.03	Pumpkin seed	0.35
Coconut	0.12	Radish	0.57



Validation Report

SENSISpec Pistachio ELISA

(Cat. No. HU0030021)

Version 4

Raw material	c [ppm]	Raw material	c [ppm]
Cod	0.20	Rice	0.49
Corn	0.15	Rye	0.28
Cow' milk	0.14	Saccharose	0.41
Cumin	0.20	Shrimps	0.21
Dill	0.22	Soy flour	0.30
Duck	0.45	Soy lecithin	0.07
Egg, dried	0.14	Split pea	0.97
Fennel	0.84	Sunflower seed	0.30
Fenugreek	0.70	Thyme	0.14
Garden cress	0.48	Tofu	0.23
Garlic, fresh	0.69	Tomato	0.13
Garlic, granulated	0.52	Turkey	0.00
Ginger, fresh	0.10	Turmeric	0.02
Ginger, ground	0.07	Walnut	0.82
Gliadin	0.09	Wheat	0.21

The following cross-reactivities could be determined:

Table 13: Cross-reactive food matrices in the SENSISpec Pistachio ELISA

Raw material	c [ppm]	Cross-reactivity [%]
Poppy seed	1.06	0.0001
Rapeseed	1.71	0.0002
Sesame	1.02	0.0001
Horseradish	1.17	0.0001
Flaxseed	1.03	0.0001
Cashew	33.3	0.0033
Mustard, yellow	1.64	0.0002

7. Robustness

Robustness was determined by variation of different handling parameters as defined in the instruction manual. The results were compared to the results of samples analyzed according to the intended method. An un-spiked cookie sample and a sample spiked with pistachio extract were analyzed, respectively.

7.1. Variation of extraction temperature

The extraction temperature defined as 60°C was changed to 25°C, 40°C and 70°C, respectively using spiked and pure cookie matrix.

Table 14: Variation of extraction temperature in the SENSISpec Pistachio ELISA

Sample	Result 25°C	Result 40°C	Result 60°C	Result 70°C
Cookie 0 ppm	< 1 ppm	< 1 ppm	< 1 ppm	< 1 ppm
Cookie 10 ppm	11.9 ppm	9.06 ppm	10.5 ppm	8.33 ppm



Validation Report

SENSISpec Pistachio ELISA

(Cat. No. HU0030021)

Version 4

Under consideration of the precision data as stated in chapter 2 the results do not differ significantly.

7.2. Variation of extraction time

The extraction time defined as 15 min was changed to 10 min and 20 min respectively.

Table 15: Variation of extraction time in the SENSISpec Pistachio ELISA

Sample	Result 10 min	Result 15 min	Result 20 min
Cookie 0 ppm	< 1 ppm	< 1 ppm	< 1 ppm
Cookie 10 ppm	9.40	8.47	9.06

Under consideration of the precision data as stated in chapter 2 the results do not differ significantly.

7.3. Drift

In contrast to the test procedure as defined in the instruction manual the incubation time of the samples was extended and reduced by 4 minutes compared to the calibrators.

Table 16: Drift in the SENSISpec Pistachio ELISA

Sample	Result 16 min	Result 20 min	Result 24 min
Cookie 0 ppm	< 1 ppm	< 1 ppm	< 1 ppm
Cookie 10 ppm	7.53 ppm	10.1 ppm	11.8 ppm

The results do not differ significantly. Nevertheless, drift in extensive test runs should be avoided by pipetting calibrators once before the samples and once after the samples, using the mean value for calculation.