INGene q PPA

R.11.PPA.K.5TX/Q

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INgene q PPA

R.11.PPA.K.5TX/Q

qPCR kit for the detection of African Swine Fever virus in biological samples. *TaqMan PCR technology.*

KIT FEATURES

APPLICATION

INgene q PPA is a real-time PCR kit suited for the detection of African Swine Fever virus (ASFV) DNA in biological samples from **domestic pigs and wild boars** of clinical interest. Our assay has been optimized and specifically tested for **blood, serum, spleen and liver**.

For veterinary diagnostic use.

TECHNICAL BASE


- INgene q PPA** is an assay for amplification and detection of the **ASFV** genetic material based on the polymerase chain reaction (PCR) and hydrolysis fluorescent probes (*TaqMan* type). The assay consists on a **duplex** real-time PCR in a one-well format.
- This kit contains specific primers and a *TaqMan* probe labelled with FAM fluorochrome for the specific detection of ASFV. In addition, the reaction mix includes an internal control (I.C.) that is amplified with specific primers and detected with a *TaqMan* probe labelled with HEX, allowing identification of false negatives due to PCR inhibition. Furthermore, the master mix contains a passive reference fluorochrome (ROX) to allow signal normalisation in the real-time PCR platforms where this is required.
- The kit has been validated in the following PCR platforms: **Quantstudio™ 5, StepOnePlus™ (Applied Biosystems™), LightCycler® 480 II (Roche)**. It is compatible with other thermocyclers, if they have the appropriate fluorescence channels, but Ct values may vary among them.

RESULTS INTERPRETATION

From each sample, fluorescence data originated from channels FAM and HEX is obtained:

- FAM channel:** pathogen detection. A positive result should show a typical amplification curve with a Ct value <45.
- HEX channel:** Internal Control detection. A positive result during amplification implies a fluorescence curve with a Ct value <45. Normally, values are comprised within a range of $30 < Ct < 45$.
- Both channels.** A negative result implies absence of amplification (undetected) or a Ct ≥ 45 .

Criteria for the kit validation: the assay will be considered as valid when the C+ have a Ct value within the range 32 ± 4 and the C-, a Ct >45 in the FAM channel.

ASSAY VALIDATION

VALIDATION DATA

The previous INgene q PPA kit was validated internally and externally by the European Reference Laboratory for ASFV (CISA-INIA) with a panel of 270 samples. Conservative modifications and improvements were introduced on this assay and both the previous and modified versions were compared in parallel with an additional set of 35 samples provided by the European Reference Laboratory. The improved kit showed 100% correlation with the previous version. Specificity was tested with 161 field samples (65 negatives, 96 positives, one sample of Classical Swine Fever and other tissue samples).

Sensitivity was calculated using 21 reference isolates belonging to 21 of the 22 p72 and 9 tick homogenates (*Ornithodoros Moubata*) obtained from warthogs of east Africa (Kenia y Uganda). The obtained values were: diagnostical sensitivity >98%, analytical sensitivity: detection of ASFV (10 gene copies), diagnostic specificity: >96.8%, repeatability: intra-assay, inter-assay, inter-day, inter-person and inter-batch variabilities were superior to 95%.

After the current registration of the kit, it was re-validated externally by the European Reference Laboratory, using the OIE assay (UPL-PCR) as a reference. Improved performance parameters were obtained in this most recent validation.

- From the analysis of 126 POSITIVE field samples (56 domestic pigs and 70 wild boars) obtained in endemic areas of Europe, the diagnostic sensitivity of the INgene q ASF real-time PCR kit was 98.41% [94.5% - 98.7% CI 95%].
- Based on the analysis of 94 blood samples, the diagnostic specificity of the INgene q PPA real-time PCR kit was 98.9% [94.2% - 99.8% CI95%].
- The value (κ) of 0.97% [CI95%] indicates PERFECT CONCORDANCE between the UPL reference method and the INgene q PPA real-time PCR kit.
- The INgene q PPA assay showed better sensitivity than the UPL assay in sample pools. It is suitable for groups of up to 20 samples, but for subclinical or chronic infections (CT >30) pools of 3-5 samples are recommended.
- The inter-assay and intra-assay variability value was less than 3% and is appropriate in routine tests.

KIT COMPOSITION

- Mixture A: ASFV and IC specific primers and probes
- Mixture B: enzyme mix
- Positive ASFV Amplification control

EXPIRATION: 12 MONTHS

Stored $\leq -18^{\circ}\text{C}$

Assay design follows OIE guidelines (Animal Health Code, 28th Edition).

*Spanish registration nº 3293;
German registration Fli-B 669.*

GOLD STANDARD DIAGNOSTICS MADRID, S.A.

C/Hermanos García Noblejas 39, 8º

28037 MADRID (SPAIN)

Tel: (+34)91 3680501

<http://www.goldstandarddiagnostics.com/>