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CADUCIDAD: 12 meses
Conservado a -20°C

Ingene q PPA R.11. PPA.K.5TX/Q


Ingene q PPA Sondas Detection is a real-time PCR kit suited for the detection of nucleic acid of **African swine fever** in biological samples in a simple way, with high levels of sensitivity and specificity.

The kit contains all necessary reagents and enzymes in **two mixes**. When mixed together, they are at the required concentrations designed to detect the pathogen, simply by adding the nucleic acid of the problem sample. This kit makes use of a UPL® probe labelled with FAM. It also includes another probe, labelled in VIC, as internal control of the reaction, allowing detection of false negatives due to an inhibition of the PCR. In addition, the master mix contains a passive reference fluorochrome, ROX, allowing the normalization of the signal, avoiding possible quantification errors deriving from pipetting.

The kit includes a positive control for detection. A standard positive control with a known number of copies can be supplied under request. This control allows generating a standard curve which links the number of pathogen copies with the cycle threshold (Ct) value.

TECHNICAL BASIS


1. Polymerase chain reaction (PCR) is based on the use of heat-stable DNA polymerase, which assembles a new DNA strand, by using DNA as a template and DNA oligonucleotides (primers) to selectively amplify the target DNA of the sample. The reaction is specific because of the bases complementarities.
2. Mix A1 includes specific primers and probe complementary to the pathogen DNA.
3. Mix B contains all necessary enzymes for specific amplification of a nucleic acid fragment from the pathogen.
4. The probe is labelled with a fluorophore (FAM) in 5' and with a quencher (NFQ) in 3'.
5. During elongation, the enzyme cuts the probe and the fluorophore is separated from the quencher, emitting fluorescence between 465-510 nm. If pathogen is not present in the sample, there is no fluorescence emission.

PREPARATION OF THE PCR

Mix A1 10 µl
Mix B 10 µl
Mix final volumen/sample = 20 µl

MIXES DISTRIBUTION


THERMOCYCLER PROGRAMME

Temperature(°C)	Time	Cycles	Reaction
95	10 min	1	PCR
95	10 s	45	
60*	30 s		

Fluorescence	Positives	Negatives	Inhibited
FAM-PPA 465-510 nm	+	-	-
VIC-CI 533-580 nm	+/-	+	-

AMPLIFICATION PLOT

INTERPRETATION OF RESULTS

A positive result during amplification implies a typical **FAM** fluorescence curve with a Ct value <40. Negative samples will amplify with fluorescence lecture in **VIC** (otherwise they should be analyzed again, diluting the starting template).

Samples with Ct's comprised between **40** and **45** should be regarded as doubtful. In that case we recommend:

- Repeating the analysis.
- Sequencing the amplified product.

PRODUCT MANUFACTURED BY INGENASA

IT-73840 IT-73780 ISO 14001:2004 9191.INGE ISO 9001:2008 9175.ING2

This product has been validated by the European Reference Laboratory for ASF :Centro de Investigación en sanidad Animal (CISA-INIA)

KIT COMPOSITION

Mix A1: (specific primers and probe for **PPA + IC**)
Mix B: Enzymes mixture
Positive control A1

Self life: 12 months
Stored at -20°C