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11.CIR.K2

I. TECHNICAL BASIS

The kit has been designed for the semiquantitative detection of circovirus VP2 protein.

The INGEZIM PCV DAS kit is based on the double antibody sandwich immunoassay technique it described below:

The solid phases are plates coated with a monoclonal antibody (MAb) specific of circovirus VP2. After adding the sample to the well, if it contains virus, they will bind

to the MAb absorbed to the plate. If the sample does not contain virus they will not bind to the MAb. After washing to eliminate the not adsorbed material, a MAb specific to VP2 conjugated with peroxidase is added. If there was virus in the sample, the conjugate will bind to the antigen. This union is detected after adding a substrate able to develop colour in presence of peroxidase.

II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the Kit.
5. Do not use components after expiration dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled.
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. For each utilisation of the Kit, control positive and negative serum must be tested in a systematic way.
10. Stop solution is a strong acid solution that must be used with precaution. In case of accidental contact with skin, wash gently with water.
11. Substrate must be handled with care, it is very sensible to light and contamination.

III. STORAGE OF COMPONENTS

All reagents and plates must be stored between +2°C and +8°C.

IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

V. PREPARATION OF REAGENTS

- **Washing solution:**
Dilute one part of the concentrate washing solution provided in the Kit with 24 parts of distilled or deionised water (i.e. 40 ml of concentrate and 960 ml of water).
- **Diluent :**
Dilute 1 part of 5x concentrated diluent with 4 parts of distilled or deionised water (100 ml of concentrate and 400 ml of water).
- **Preparation of conjugate:**
It is ready to use. Do not dilute.
- **Preparation of positive and negative control.**
They are ready to use. Do not dilute.

VI. PREPARATION OF SAMPLES:

Vaccine samples must be assayed at two fold dilutions from 1/25 in the supplied diluent. It is recommended to make several dilutions because the absorbance obtained depends on the viral charge of the sample.

VII. TEST PROCEDURE

- All reagents must be allowed to come to room temperature before use.
- Sample and controls:**
Add to the respective wells 100 µl of positive and negative control in duplicate. Add 100 µl of dilutions of the samples (prepared as it is described above) in duplicated.
- Seal the plate and incubate **1 hour at room temperature**.
- Wash 5 times following the described procedure.
- Add 100 µl of conjugate to each well. Seal the plate and **incubate 30 minutes at room temperature**.
- Wash 5 times following the described procedure.
- Add 100 µl of substrate to each well and incubate for **15 minutes** at room temperature. (Count the time since the first well had been filled).
- Add 100 µl of stop solution to each well. **ATTENTION:** The stop solution must be dispensed in the same order in which the substrate was added.
- Read the OD of each well at 450 nm within 5 min after the addition of stop solution.

VIII. READING AND RESULT INTERPRETATION

Read the OD at 450 nm

A) Validation of the test:

The test will be valid if:

- The OD obtained with the positive control is higher than 0.8
- The OD obtained with the negative control is lower than 0.3

B) Interpretation of the results:

The viral charge of the sample will be directly proportional to its titer.

The sample titer depends on the user criteria. One general way to define a titer is the last dilution showing an OD higher than negative control OD multiplied by 3

C) Example:

- Negative control OD = 0.092
- Positive control OD = 2.065
- Negative control OD x 3 = 0,092 x 3 = 0,276

	Sample 1	Sample 2
25	3,357	4,99
50	1,843	3,742
100	0,896	2,17
200	0,497	1,121
400	0,317	0,61
800	0,236	0,359
1600	0,152	0,242
3200	0,122	0,143

The titer of Sample 1 will be 1/400 and the sample 2 1/800.

Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tlf: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in _____ by: