## INGEZIM PPA DAS

Prod Ref: 11.PPA.K2

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## I. TECHNICAL BASIS

---

This kit is based on a double antibody sandwich enzymatic immunoassay (DAS or capture Elisa). We make a brief description of the technique below:

The plate is coated with a specific monoclonal antibody (Mab) to VP73 of the ASFV. In this way when a sample is added into the wells, if they contain antigen, the coated Mabs will catch it. After washing a second Mab to VP73 (specific for a different epitope) conjugated with Peroxidase, is added.

If there is antigen captured into the well, the conjugate (Mab, labelled with peroxidase) will bind it and the peroxidase action will be detected after adding the adequate substrate.

Samples to be used in the assay could be any tissue suspected to be infected. For diagnoses purposes, we recommend to use as sample, blood, spleen or lymphatic nodes (tissues in which viral concentration could be high).

## II. PRECAUTIONS AND WARNINGS FOR USERS:

---

1. Read the use instructions carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the kit.
5. Do not use components after expiration dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking where specimens or kit reagents are being handle.
7. Do not pipette by mouth.
8. Use a new tip for each serum sample. Use only distilled water for the preparation of the reagents.
9. For each utilisation of the kit, control positive and negative serum must be tested in a systematic way
10. Substrate and substrate buffer must be handle with care.

## III. STORAGE OF COMPONENTS

---

All components must be stored between +2°C and +8°C. **Once opened**, control sera are stable for one month. In case that they are not going to be used in this period, we recommend to store them at -20°C.

**WARNING! Stop Solution may precipitate, if it happens warm it at +37°C.**

## IV. INFORMATION ABOUT THE WASHING STEPS

---

The washing steps could be done using an automatic washing machine or a multichanel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.

- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use.  
Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

## V. PREPARATION OF REAGENTS

---

- **Washing solution:**  
Dilute one part of the 25x concentrate washing solution provided in the kit with 24 parts of

distilled or deionized water. When ready this solution remains stable at +4°

- **Controls:**  
Before use a 1/50 dilution must be done in diluent. Once opened, controls are stable for 1 month between +2°C and +8°C.  
For long term storage, we recommend to make aliquots of the controls and store it at -20°C till use.
- **Preparation of the conjugate: to make immediately before use.**  
Dilute the needed quantity of conjugate provided in the kit 1/100 into diluent (e.g. 10 µl of conjugate into 1ml of diluent is enough quantity for one strip of 8 wells).  
Shake very well the solution before the use.  
Prepare only the quantity needed for each time because the remainder volume has to be rejected
- **Preparation of the substrate:**  
Ready to use.

## VI. PREPARATION OF SAMPLES:

For diagnosis purposes, samples must be whole blood, spleen or lymphatic nodes.

**For whole blood** (with heparin or EDTA), before use, samples must be frozen and thawed once to avoid false positives results.

Samples are used without dilution (100 µl/well).

**For spleen or lymphatic nodes**, a 1/10 dilution must be done (w/v) using the diluent supplied (i.e. 1g of sample in 10 ml of diluent). Make a good homogenisation of the sample dilution using a stomacher or any other mechanic homogenizator.

## VII. TEST PROCEDURE

1. All reagents (except conjugate) must be allowed to come to room temperature before use.
2. Add 100 µl of positive control to two wells (e.g. A1 and B1), and 100 µl of the negative control (e.g. A2 and B2). Add 100 µl of samples to test on each remainder wells. We recommend the use of two wells per sample. Seal the plate **and incubate for 1 h at 37°C**.
3. **Empty the wells into a receptacle containing 0.1 M NaOH** and wash four times as specified in point IV.
4. Add 100 µl of specific conjugate (prepared following previous instructions) to each well. Seal the plate and **incubate for 1 hour at room temperature (25°C)**.
5. Wash 4 times following the described procedure.
6. Add 100 µl of substrate to each well. Keep the plate **for 10 min at room temperature**.
7. Add 100 µl of stop solution to each well.
8. Read the OD of each well at 405nm within 5 min after the addition of stop solution.

## VIII. READING AND RESULT INTERPRETATION

### Validation of the test

The test could be considered valid when:

- ⇒ The OD of the Positive control is higher that 1.0
- ⇒ The OD of the Negative control is lower that 0,150

### Interpretation:

When you are running duplicates, OD of the sample, will be calculated as the arithmetic mean of OD values in both wells.

- Samples with an OD higher than 0.300 must be considered as **POSITIVES** to ASFV.
- Samples with OD lower that 0,150 must be considered as **NEGATIVES** to ASFV.
- Samples with OD values between both values are considered as **DOUBTFUL**. We recommend re-testing these animals one more time or applying a different technique to check this sample

Desarrollado y fabricado en España por:  
Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA  
APLICADA, S.A.  
C/ Hnos. García Noblejas, 39  
28037- MADRID (SPAIN)  
Tlf: +34 91368.05.01/04  
Fax: +34 91 408.75.98  
E-mail: ingenasa@ingenasa.com  
www.ingenasa.es

Distributed in \_\_\_\_\_ by:

IT-73840  
IT-73780

ISO 14001:2015  
9191.INGE

ISO 9001:2015  
9175.ING2