# INGEZIM PPA CROM ANTICUERPO

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## I. TECHNICAL BASIS

---

The test is based on the technique of immunochromatography, a migration technique that uses purified VP72 protein of ASFV, and is able of detecting specific antibodies against ASFV in porcine serum samples. The diagnostic device consists of two windows:

- Sample window: it contains VP72 protein of ASFV and a protein control, coated to coloured latex particles.
- Result window: it contains a test line (T) formed by the VP72 protein of ASFV and a control line (C) formed by a MAb specific for the protein control.

If the test sample contains antibodies against the VP72 protein of ASFV, they will react with the red particles conjugated to VP72 protein. The latex-protein-antibody complex migrates through the membrane and reacts again with the VP72 protein adsorbed in the test line resulting in the appearance of a red/pink. The appearance of a blue line in the control area (C) indicates that the chromatography has been correctly performed.

## II. STORAGE

---

All the reagents must be kept between +4°C and +25 °C.

## III. PRECAUTIONS

---

1. Read the use instructions carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the kit.
5. Do not use components after expiration date and do not mix components from different batches.
6. Do not eat, drink, or smoke where specimens or kit reagents are being handle.

## IV. SAMPLES

---

The test has been designed for use with serum or blood samples. Serum samples must be fresh, refrigerated up to 7 days at 2-8 °C or stored at - 20 ° C. Blood samples must be fresh or refrigerated up to 4 days at 2-8 °C, and collected with any anticoagulant (EDTA, heparin, etc.). Blood samples without any anticoagulant, could contain micro-clots that can block the device and result in nonspecific reactions.

## V. TEST PROCEDURE

---

- Dilution buffer and device must be at room temperature at the moment of test performance.
- Open the aluminium bag just in the moment to realise the assay. Remove the test device from the bag and place it on a flat surface.

### INSTRUCTIONS FOR WHOLE BLOOD SAMPLES

1. Take 20 µl of whole blood using the provided micropipette and add the sample to the round window. To take accurately the sample touch the drop of sample with the micropipette and it will fill up until the black line. Dispense it quickly to avoid coagulation inside the micropipette.
2. For dispensing the sample on the device, place the micropipette carefully on the round window covering the hole and pressing smoothly.

## INGEZIM PPA CROM *antibody* 11.PPA.K41 WB

3. Wait for **1 minute** and add 3 drops of the running buffer in the same round window (**wait total liquid absorption into the membrane between drop and drop**).
4. Read the result at 10 minutes.

### INSTRUCTIONS FOR SERA SAMPLES

1. Take 10  $\mu$ l of serum sample using the 10  $\mu$ l pipette and the tips.
2. For dispensing the sample on the device, place the micropipette carefully on the round window and press.
3. Add 3 drops of the running buffer in the same round window (**wait total liquid absorption into the membrane between drop and drop**).
4. Read the result at 10 minutes.

## VI. INTERPRETATION OF RESULTS


In order to be sure that chromatographic assay has run properly, a blue line must appear in control area (C).

- **POSITIVE:** Two lines appear: a blue control line and a red/pink line in the result window indicating a positive result. Any readings after 10 minutes are of no diagnostic value.
- **NEGATIVE:** only a blue line appears in the result window.
- **INVALID:** when blue line does not appear in the result window.

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
IT IS VERY IMPORTANT TO FOLLOW THESE INSTRUCTIONS AND RESPECT THE INDICATED  
VOLUMES TO OBTAIN VALID RESULTS

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1. Add 10 µl of serum sample on the round window.

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INMUNOLOGIA Y GENÉTICA APLICADA, S.A.