



SENSIStrip Allergen Lateral Flow Tests

The **SENSIStrip Allergen lateral flow devices** are designed for quick detection of the different allergens that must be labelled on packaged foodstuffs.


The kits are based on the principle of immunoassay. With only a few minutes of incubation time, they provide a very easy-to-use technique for adequate control at different stages of the food manufacturing process, across food matrices, rinse water and equipment.

Discover our new PowerLine portfolio to eliminate false negatives due to the hook effect: excess of allergens in highly contaminated samples.

Easy to use and quick method

The SENSIStrip lateral flow devices are sensitive tools with very low limits of detection.

The tests are performed in a fast and easy way, without the need for a highly qualified end user. The kits are provided with all the necessary and ready-to-use material to perform 20 individual analyses per kit, on food matrices or environmental samples. In addition, an evaluation card is included to enable a semi-quantitative result, based on the intensity of the test line on the strip.



The PowerLine strips include a third Hook (H) line to indicate a very high dose of allergens.



SENSIStrip Lateral Flow Tests

LFD kits for NUTS allergens detection

SENSIStrip Almond	HU0030110
SENSIStrip Brazil nut	HU0030127
SENSIStrip Cashew	HU0030113
SENSIStrip Cashew PowerLine	HU0030207
SENSIStrip Coconut	HU0030114
SENSIStrip Hazelnut	HU0030119
SENSIStrip Hazelnut PowerLine NEW!	HU0030206
SENSIStrip Macadamia nut	HU0030122
SENSIStrip Multi Tree Nut	HU0030136
SENSIStrip Peanut	HU0030128
SENSIStrip Peanut PowerLine	HU0030137
SENSIStrip Pecan nut	HU0030129
SENSIStrip Pistachio	HU0030130
SENSIStrip Pistachio PowerLine NEW!	HU0030255
SENSIStrip Walnut	HU0030133

LFD kits for GLUTEN detection

SENSIStrip Gluten	HU0030118
SENSIStrip Gluten PowerLine	HU0030134

LFD kits for SEAFOOD allergens detection

SENSIStrip Crustaceans (Tropomyosin)	HU0030115
SENSIStrip Crustaceans PowerLine NEW!*	HU0030254
SENSIStrip Fish	HU0030117
SENSIStrip Fish PowerLine	HU0030271

LFD kits for MILK allergens detection

SENSIStrip Beta-Lactoglobulin	HU0030111
SENSIStrip Casein	HU0030112
SENSIStrip Casein PowerLine	HU0030242
SENSIStrip Milk Protein	HU0030123
SENSIStrip Milk Protein PowerLine	HU0030138

LFD kits for EGG allergens detection

SENSIStrip Egg White	HU0030116
SENSIStrip Egg Protein PowerLine	HU0030202
SENSIStrip Lysozyme	HU0030121
SENSIStrip Ovalbumin	HU0030126
SENSIStrip Ovalbumin PowerLine NEW!*	HU0030259

LFD kits for OTHER allergens detection

SENSIStrip Lupin	HU0030120
SENSIStrip Lupin PowerLine NEW!*	HU0030256
SENSIStrip Mustard	HU0030125
SENSIStrip Mustard PowerLine NEW!	HU0030204
SENSIStrip Sesame	HU0030131
SENSIStrip Sesame PowerLine NEW!*	HU0030270
SENSIStrip Total Soy Protein	HU0030132
SENSIStrip Total Soy Protein PowerLine	HU0030135
SENSIStrip Walnut	HU0030133

*Available from July 2026

Lateral Flow Reader for the SENSIStrip Kits: RapidScan ST5-W

No more doubts on interpreting a possible line - get the result in ppm

The RapidScan ST5-W Lateral Flow Assay Reader is an accurate, easy to use and highly flexible imaging platform used for lateral flow assay analysis. The instrument was designed specifically for field and in-process testing applications that require qualitative and quantitative test results.



Article Nr	Products
LFRSCAN002	RapidScan ST5-W Lateral Flow Assay Reader
LFA0030001	LFD Reader Adapter for allergen strips