

Certificate of Product**Product:** BACGene Listeria spp.**Lot:** 26191803**Cat. no.:** 5123222101**Exp. date:** 08-2020

Item	Lot-No.	Store at:
BACGene Listeria spp., PCR plate with pre-dispensed MasterMix	25411803	-20°C
Plasmid DNA pGSE670 (200 copies/μL), 50 μL	25031801	-20°C
Lysisbuffer L, 2.5 mL	25041806	-20°C
Proteinase K for Lysis Buffer L, 500 μL	25081806	-20°C
Lysozym for Lysis Buffer L, 500 μL	25051806	-20°C
Lysis plate, for sample preparation, rippable	19341401	RT
Domed caps, for use with lysis plate, set of 12 strips	19351401	RT
Optical caps, for use with PCR plate, mat of 12 strips, rippable	19361401	RT
PCR support grid	25611601	

Quality Control:

This product has successfully passed quality control and meets our internal specifications. The Batch release was performed with strain Listeria monocytogenes 4b, DSMZ no.: 15675, NCTC no.: 11994.

The production of the product is carried out with best carefulness. We guarantee the above mentioned specification through our QM system. The CoP doesn't release you from performing the standard control upon receipt of incoming goods.