



SENSISwab Swabbing Kit **NEW!**

The Eurofins Technologies SENSISwab kit is a reliable support to monitor surface contamination and to verify the cleanliness of surfaces. The appropriate use of swab samples indicates the presence of allergen proteins.

The kit is specifically designed and validated for the application of SENSISpec allergen ELISA kits and SENSIStrip lateral flow devices. The SENSISwab kit contains consumables and reagents for modules of 5 and 10x5 tests. The items need to be stored at 2–30°C. The shelf life of the kit is 24 months.

Product specifications

Article no.	Description	No. of Swabs
HU0030101	SENSISwab kit, 10 ml Swabbing solution, 5 Swabs, 5 Sample tubes	5
HU0030102	SENSISwab 10x5: multipack of 10 Swab kits	50

Order Our Diagnostic Kits Online! www.eurofins-technologies.com