



## Rapid Quantitative Patulin Test Kit

For juice processors, infant food producers and food laboratories requiring faster turnaround times and quantitative test results for Patulin, a harmful mycotoxin found in fruits, fruit juices and more.

**Gold Standard Diagnostics** Patulin ELISA (enzyme-linked immunosorbent assay) 96 micro-well test kits help identify patulin contamination in juice production on-site and can benefit juice companies by saving time and costs for the quality and fast control of products through the entire juice production/food chain supporting informed decision-making.

### Key Benefits of Patulin ELISA Analysis

- Rapid, laboratory-based analysis - answers in 3-4 hours for up to 41 samples
- Applications for apple juice, apple cider, applesauce and orange juice (contact Gold Standard Diagnostics for additional fruit juice/products)
- Test Kit Range - 7.0 ppb to 300 ppb
- Low per test cost
- Simple sample extraction
- Highest quality antibody with excellent affinity, specificity
- Meets global regulatory action levels
- Improves turnaround by eliminating delays associated with samples shipment for offsite analysis
- Tests operate with standard lab equipment, microplate readers

**Rapid Patulin Kits offer processors brand name protection, compliance, prevention of costly product recalls and ensure consumer confidence.**

### Ordering

Our patulin ELISA plate kits include materials and reagents for sample derivatization and analysis.

Part Number	Product Description
500106	Patulin ELISA 96-test kit
Additional Products (if needed)	
301102	Patulin Standard & Control Set
301103	Patulin Derivatization Set