



VETERINARY DIAGNOSTICS

**ELISA and Lateral Flow Assays
Product Catalogue**





VETERINARY DIAGNOSTICS

In industrial livestock production, operating conditions are optimized to achieve the highest production output while minimizing costs. Given the high density of livestock and the production pressure, the emergence of infectious diseases is a major risk, and therefore prevention and control measures are essential.

Our diagnostic kits support the animal health management by:

- Detecting and monitoring epidemiological risks,
- Supporting vaccines & treatment development,
- Controlling the quality of produced meat,
- Improving animal welfare.

Veterinary Diagnostic Kits

We are engaged in the development, manufacturing and distribution of ELISA, Lateral Flow and PCR diagnostic kits, with an emphasis on infectious diseases for both farm animals and pets. Our easy-to-use and reliable test kits, based on high-quality reagents, are used for the diagnosis of bacterial, viral and parasitic diseases.

Veterinary diagnostic ELISA kits

- Easy and reliable tests, based on high quality reagents (recombinant antigens and monoclonal antibodies)
- Wide range of infectious diseases and animal species covered
- Different presentations available for the best adaptation to the requirements of each client

Veterinary diagnostic Lateral Flow Assays

- Ease of use
- Can be applied directly on the field
- High sensitivity and specificity rates

MULTISPECIES specialities



Product	Product description	No. of tests	Article No.
BRUCELLOSIS			
INgesim Brucella Compac 2.0	Blocking ELISA for the detection of specific antibodies to Brucella spp. in ruminants and porcine serum samples.	2 plates kit	192
		5 plates kit	480
10.BRU.K.3/2			
10.BRU.K.3/5			
CRYPTOSPORIDIOSIS			
INgezim Crypto CROM	Lateral Flow Assay for the detection of Cryptosporidium antigen in biological samples (faeces from all species).	12 individual devices 30 individual devices	12 30
10.CRY.K.42/12			
10.CRY.K.42/30			
GIARDIOSIS			
INgezim Giardia CROM	Lateral Flow Assay for the detection of Giardia antigen in biological samples (faeces from all species).	12 individual devices 30 individual devices	12 30
10.GIA.K.42/12			
10.GIA.K.42/30			
INFLUENZA A			
INgezim Influenza A	Blocking ELISA for the detection of specific antibodies to Influenza type A virus in serum samples from different species.	2 plates kit 5 plates kit	192 480
10.FLU.K.3/2			
10.FLU.K.3/5			
ROTAVIROSIS			
INgezim Rotavirus DAS	Double antibody ELISA for the detection of Rotavirus (Type A) antigen in biological samples (faeces from all species).	1 plate kit	96
10.RT.K.2			
INgezim Rota CROM	Lateral Flow Assay for the detection of Rotavirus (Type A) antigen in biological samples (faeces from all species).	12 individual devices 30 individual devices	12 30
10.RT.K.42/12			
10.RT.K.42/30			
TUBERCULOSIS			
INgezim Tuberculosis DR	Double recognition ELISA for specific antibodies detection to Mycobacterium bovis in serum samples from wild and domestic animals, except cattle.	2 plates kit 5 plates kit	192 480
10.TB.K.0/2			
10.TB.K.0/5			
WEST NILE FEVER			
INgezim West Nile Compac	Blocking ELISA for the detection of specific antibodies to WNV in birds and equine serum samples.	2 plates kit 5 plates kit	192 480
10.WNV.K.3/2			
10.WNV.K.3/5			
SARS-CoV-2			
INgezim Covid Vet	Indirect ELISA for the detection of specific antibodies against S protein of the SARS-CoV-2 virus in serum and plasma samples of minks, ferrets, cats and dogs.	2 plates kit 5 plates kit	192 480
10.CoS.K.1/2			
10.CoS.K.1/5			
BORRELIA			
Borrelia burgdorferi Veterinary IgG ELISA	Indirect ELISA for the detection of specific IgG-antibodies to Borrelia burgdorferi in dog or horse serum. (Virotech)	1 plate kit	96 tests
10.DC122.00			
Borrelia Veterinary plus OspA IgG LINE Immunoblot (anti-dog)	Immunoblot test kit for the quantitative detection of Borrelia burgdorferi sensu lato specific IgG antibodies in dog or horse serum. This test kit can distinguish a wildtype infection from a vaccination. (Virotech)		32 tests
10.DE226G32			

PORCINE specialities



Product	Product description	No. of tests	Article No.
AUJESZKY DISEASE (ADV or PRV)			
INgezim ADV TOTAL	Indirect ELISA for the detection of specific total antibodies to ADV (vaccinated and infected animals) in porcine serum samples.	2 plates kit	192
		5 plates kit	480
11.ADV.K.1/2			
11.ADV.K.1/5			
INgezim ADV gE Plus	Blocking ELISA for the detection of specific antibodies to ADV gE glycoprotein in porcine serum samples.	2 plates kit 5 plates kit	192 480
11.GEP.K.3/2			
11.GEP.K.3/5			
INgezim ADV gB	Blocking ELISA for the detection of specific antibodies to ADV gB glycoprotein in porcine serum samples.	2 plates kit 5 plates kit	192 480
11.GB.K.3/2			
11.GB.K.3/5			
CLASSICAL SWINE FEVER (CSFV)			
INgezim CSF Compac	Blocking ELISA for the detection of specific antibodies to Classical Swine Fever Virus in porcine serum samples.	2 plates kit 5 plates kit	192 480
11.CSF.K.3/2			
11.CSF.K.3/5			
AFRICAN SWINE FEVER (ASFV)			
INgezim PPA DAS 2.0	Double antibody ELISA for the detection of ASFV antigen in biological samples.	1 plate kit 5 plates kit	96 480
11.PPA.K.2			
11.PPA.K.2/5			
INgezim ASF CROM Ag	Lateral Flow Assay for the detection of ASFV antigen in porcine blood samples.	30 individual devices 100 individual devices	30 100
11.ASF.K.42/30			
11.ASF.K.42/100			
INgezim ASFV-R	Immunoenzymatic assay based on indirect ELISA technique for the detection of specific antibodies to cp312 and p30 proteins of ASFV in serum, blood (fresh or on paper) and spleen exudate samples, from swine and wild boar.	2 plates kit 5 plates kit	192 480
11.ASF.K1/2			
11.ASF.K1/5			
INgezim PPA CROM Anticuerpo	Lateral Flow Assay for the detection of specific antibodies to ASFV in porcine serum and whole blood samples.	30 individual devices 100 individual devices	30 100
11.PPA.K.41/30			
11.PPA.K.41/100			
INgezim PPA Compac	Blocking ELISA for the detection of specific antibodies to ASFV in porcine serum samples.	2 plates kit 5 plates kit	192 480
11.PPA.K.3/2			
11.PPA.K.3/5			
INgezim ASFV/CSFV CROM Ab	Specific antibodies against to classical swine fever and African swine fever virus	25 individual devices	25
11.SFV.K41			
PORCINE BRUCELLOSIS			
INgezim Brucella Porcina	Indirect ELISA for the detection of specific antibodies to Brucella suis in porcine serum samples.	2 plates kit 5 plates kit	192 480
11.BP.K.1/2			
11.BP.K.1/5			
ERYSIPELAS			
INgezim Mal Rojo	Indirect ELISA for the detection of specific antibodies to E. Rhusiopathiae in porcine serum samples.	2 plates kit 5 plates kit	192 480
11.MR.K.1/2			
11.MR.K.1/5			
GLÄSSER DISEASE			
INgezim Haemophilus	Indirect ELISA for the detection of specific antibodies to Haemophilus parasuis in porcine serum samples.	2+2 plates kit	192
11.HPS.K.1/4			
ENZOOTIC PNEUMONIA			
INgezim M. hyo Compac	Blocking ELISA for the detection of specific antibodies to Mycoplasma hyopneumoniae in porcine serum samples.	2 plates kit 5 plates kit	192 480
11.MHY.K.3/2			
11.MHY.K.3/5			
PORCINE INFLUENZA			
INgezim Influenza Porcina	Indirect ELISA for the detection of specific antibodies to SIV in porcine serum samples.	2 plates kit 5 plates kit	192 480
11.FLU.K.1/2			
11.FLU.K.1/5			
PORCINE CIRCOVIRUS			
INgezim PCV DAS	Double antibody ELISA for the detection of PCV2 antigen in Virus culture samples.	1 plate kit	96
11.CIR.K.2			
INgezim Circovirus IgM/IgG	Capture ELISA for the detection of specific IgMs and IgGs to PCV2 in porcine serum samples.	2+2 plates kit	192
11.PCV.K.2			
INgezim Circo IgG	Indirect ELISA for the detection of specific antibodies to PCV2 in porcine serum samples.	2 plates kit 5 plates kit	192 480
11.PCV.K.1/2			
11.PCV.K.1/5			
PORCINE PARVOVIRUS (PPV)			
INgezim PPV	Indirect ELISA for the detection of specific antibodies to PPV in porcine serum samples.	2 plates kit 5 plates kit	192 480
11.PPV.K.1/2			
11.PPV.K.1/5			
INgezim PPV Compac	Blocking ELISA for the detection of specific antibodies to PPV in porcine serum samples.	2 plates kit 5 plates kit	192 480
11.PPV.K.3/2			
11.PPV.K.3/5			
INgezim PPV DAS	Double antibody ELISA for the detection of PPV antigen in biological samples including virus cultures.	1 plate kit	96
11.PPV.K.2			
PORCINE CORONAVIROSI (PRCV) AND TRANSMISSIBLE GASTROENTERITIS (TGEV)			
INgezim Corona Diferencial 2.0	Blocking ELISA for the detection and differentiation of antibodies to TGEV and PRCV in porcine serum samples.	2 plates kit 5 plates kit	192 480
11.DIF.K.3/2			
11.DIF.K.3/5			
INgezim TGEV 2.0	Blocking ELISA for the detection of specific antibodies to TGEV in porcine serum samples.	2 plates kit 5 plates kit	192 480
11.TGE.K.3/2			
11.TGE.K.3/5			
INgezim Coronavirus Porcino 2.0	Blocking ELISA for the detection of specific antibodies to PRCV/TGEV in porcine serum samples.	2 plates kit 5 plates kit	192 480
11.PRC.K.3/2			
11.PRC.K.3/5			
INgezim PEDV	Indirect ELISA for detection of specific antibodies to PEDV in porcine serum samples.	2 plates kit 5 plates kit	192 480
11.PED.K.1/2			
11.PED.K.1/5			

Product	Product description	No. of tests	Article No.	
PORCINE PLEUROPNEUMONIA				
INgezim APP MIX	Indirect ELISA for the detection of antibodies against Actinobacillus pleuropneumoniae serotypes 1-2-9-11, 3-6-8-15 and 4-5-7 in porcine serum.	6 plates kit	576	11.APP.K1/6
INgezim APP MIX 10,12	Indirect ELISA for the detection of antibodies against Actinobacillus pleuropneumoniae serotypes 10 and 12 in porcine serum.	2 plates kit	192	11.APX.K1/2
PORCINE REPRODUCTIVE AND RESPIRATORY SYNDROME (PRRSV)				
INgezim PRRS DR	Double recognition ELISA for the detection of specific antibodies to PRRSV in porcine serum samples.	2 plates kit 5 plates kit	192 480	11.PRS.K.0/2 11.PRS.K.0/5
INgezim PRRS Universal	Indirect ELISA for the detection of specific antibodies to PRRSV (EU and US isolates) in porcine serum samples.	2 plates kit 5 plates kit	192 480	11.PRU.K.1/2 11.PRU.K.1/5
INgezim PRRS 2.0	Indirect ELISA for the detection of specific antibodies to PRRSV (EU and US strains) in porcine serum samples.	2 plates kit 5 plates kit	192 480	11.PR2.K.1/2 11.PR2.K.1/5
PORCINE ROTAVIROSI				
INgezim Rotavirus Porcino	Indirect ELISA for the detection of specific antibodies to Rotavirus (Type A) in porcine serum samples.	2 plates kit	192	11.RTP.K.1/2
PORCINE TUBERCULOSIS				
INgezim TB porcine	Indirect ELISA for the detection of specific antibodies to Mycobacterium bovis in wild boar and pig serum.	2 plates kit 5 plates kit	192 480	11.TBP.K.1/2 11.TBP.K.1/5
INgezim TB CROM Ab	Double recognition Lateral Flow Assay for specific detection of Mycobacterium bovis antibodies in swine and wild boar blood and serum.	30 individual devices 100 individual devices	30 100	11.TBP.K41/30 11.TBP.K41/100

RUMINANTS specialities



Product	Product description	No. of tests	Article No.	
BOVINE BESNOITIOSIS				
INgezim Besnoitia	Indirect ELISA for the detection of specific antibodies to Besnoitia besnoiti in bovine serum samples.	2 plates 5 plates	192 480	12.BES.K.1/2 12.BES.K.1/5
BOVINE LEUKOSIS (BLV)				
INgezim BLV Compac 2.0	Blocking ELISA for the detection of specific antibodies to BLV in bovine serum samples and milk samples.	2 plates kit 5 plates kit 10 plates kit	192 480 960	12.BLV.K.3/2 12.BLV.K.3/5 12.BLV.K.3/10
INgezim BLV Confirmación	Indirect ELISA for the confirmation of specific antibodies to BLV in bovine serum samples.	5+5 plates kit	480	12.BLV.K.1/5
BOVINE RESPIRATORY SYNCYTIAL (BRSV)				
INgezim BRSV Compac	Blocking ELISA for the detection of specific antibodies to BRSV in bovine serum samples.	2 plates kit 5 plates kit	192 480	12.BRS.K.3/2 12.BRS.K.3/5
BRUCELOSIS				
INgezim Brucella Bovina Milk	Indirect ELISA for the detection of specific antibodies to Brucella in bovine milk samples (individual or tank samples)	6 plates kit	576	12.BB.K.1/6
INgezim Brucella Bovina 2.0	Indirect ELISA for the detection of specific antibodies to Brucella in bovine serum samples (individual or pooled samples).	2 plates kit 5 plates kit 10 plates kit	192 480 960	12.BB2.K.1/2 12.BB2.K.1/5 12.BB2.K.1/10
INgezim Brucella ovis	Indirect ELISA for the detection of specific antibodies to Brucella ovis in small ruminants serum samples.	2 plates kit 5 plates kit	192 480	13.BO.K.1/2 13.BO.K.1/5
INgezim Brucella Small Ruminants	Indirect ELISA for the detection of specific antibodies to Brucella in ruminants serum samples.	2 plates kit 5 plates kit	192 480	13.BM.K.1/2 13.BM.K.1/5
INgezim AGID	Agar Gel immuno-diffusion (AGID) assay for the detection of specific antibodies (infection Ab and no vaccination Ab) present in ruminant serum samples infected with Brucella abortus and melitensis.	6 Agar plates with 14 wells per plate	84	12.BB2.K.7/84

Product	Product description	No. of tests	Article No.	
BOVINE ROTAVIROSI				
INgezim Rotavirus Bovino	Indirect ELISA for the detection of specific antibodies to RTB (Type A) in ruminants serum samples.	2 plates kit	192	12.RT.K.1/2
BLUE TONGUE (BTV)				
INgezim BTV DR	Double recognition ELISA for the detection of specific antibodies to BTV in ruminants serum samples.	2 plates kit 5 plates kit	192 480	12.BTV.K.0/2 12.BTV.K.0/5
INgezim BTV DAS	Double antibody ELISA for the detection of BTV antigen in BTV culture samples (to evaluate the performance of the culturing).	1 plate kit	96	12.BTV.K.2
INgezim BTV CROM	Indirect Lateral Flow Assay for the detection of specific antibodies to BTV in ruminants serum samples.	30 individual devices	30	12.BTV.K.41/30
INgezim BTV Compac 2.0	Blocking ELISA for the detection of specific antibodies to BTV in bovine serum samples (individual samples).	2 plates kit 5 plates kit	192 480	12.BTV.K.3/2 12.BTV.K.3/5
MAEDI VISNA - CAPRINE ARTHRITIS ENCEPHALITIS				
INgezim Maedi Screening	Indirect ELISA for the detection of specific antibodies to small ruminants Lentivirus. in small ruminants serum samples.	2 plates kit 5 plates kit	192 480	13.SLS.K.1/2 13.SLS.K.1/5
INgezim Maedi Confirmación	Indirect ELISA for the detection of specific antibodies to small ruminants Lentivirus in small ruminants serum samples.	2+2 plates kit	192	13.SLC.K.1
BOVINE VIRAL DIARRHEA (BVDV)				
INgezim BVD DAS	Double Antibody Sandwich ELISA, for the detection of BVDV antigen in blood, sera, plasma, leukocytes, ear tissue and cell culture in bovine.	1 plate kit 5 plate kit	96 480	12.BVD.K.2/1 12.BVD.K.2/5
INFECTIOUS BOVINE RHINOTRACHEITIS (IBRV)				
INgezim IBR gE Compac	Blocking ELISA for the detection of specific antibodies to gE protein from Infectious bovine Rhinotracheitis Virus (IBRV) in bovine serum samples.	2 plates kit 5 plates kit	192 480	12.IBE.K.3/2 12.IBE.K.3/5
INgezim IBR Compac 2.0.	Blocking ELISA for the detection of specific antibodies to IBRV gB glycoprotein (total) in bovine serum samples.	2 plates kit 5 plates kit	192 480	12.BHV.K.3/2 12.BHV.K.3/5
INgezim IBR 2.0.	Indirect ELISA for the detection of specific antibodies to IBRV (total) in bovine serum, plasma and milk samples.	2 plates kit 5 plates kit	192 480	12.BHV.K.1/2 12.BHV.K.1/5
RIFT VALLEY FEVER (RVFV)				
INgezim FVR Compac	Blocking ELISA for the detection of specific antibodies to RVFV in ruminants serum samples.	2 plates kit 5 plates kit	192 480	13.FVR.K.3/2 13.FVR.K.3/5
INgezim FVR IgM	Capture enzymatic immunoassay for the detection of specific antibodies (IgM) to RVFV in ruminants serum samples.	2 plates kit 5 plates kit	192 480	13.FVR.K.2/2 13.FVR.K2/5
BOVINE NEOSPOROSIS				
INgezim Neospora 3.0.	Indirect ELISA for the detection of specific antibodies to Neospora caninum in bovine serum samples.	2 plates kit 5 plates kit	192 480	12.NC.K.1/2 12.NC.K.1/5
SCHMALLEBERG DISEASE (SBV)				
INgezim Schmallenberg Compac	Blocking ELISA for the detection of specific antibodies to Schmallenberg virus in ruminants serum samples.	2 plates kit 5 plates kit	192 480	13.SBV.K.3/2 13.SBV.K.3/5
PESTE DES PETITS RUMINANTS (PPRV)				
INgezim PPR Compac	Detection of specific antibodies to pestiviruses (BVDV y BDV) in serum and milk samples of bovine, ovine and goat.	2 plates kit 5 plates kit 10 plates kit	192 480 960	13.PPR.K.3/2 13.PPR.K.3/5 13.PPR.K3/10
PESTIVIRUS				
INgezim Pestivirus Compac	Indirect ELISA for the detection of specific antibodies to Pestivirus in serum and milk ruminant samples.	2 plates kit 5 plates kit	192 480	12.BVD.K.3/2 12.BVD.K.3/5
PARATUBERCULOSIS				
INgezim Paratuberculosis	Indirect ELISA for the detection of specific antibodies to Mycobacterium avium spp paratuberculosis in serum and milk ruminant samples.	2 plates kit 5 plates kit	192 480	12.PTB.K.1/5 12.PTB.K.1/2

CANINE specialities



Product	Product description	No. of tests	Article No.
CANINE CORONAVIRUS DISEASE (CCV)			
INgezim Corona Canino	Indirect ELISA for the detection of specific antibodies to CCV in canine serum samples.	1 plate kit	96 15.CCV.K.1
CANINE DISTEMPER (CDV)			
INgezim Moquillo IgG	Indirect ELISA for the detection of specific antibodies to CDV (IgGs) in canine serum samples.	1 plate kit	96 15.CDG.K.1
INgezim Moquillo IgM	Capture antibody ELISA for the detection of specific antibodies to CDV (IgMs) in canine serum samples.	1 plate kit	96 15.CDM.K.2
CANINE DISTEMPER (CDV)			
INgezim CDV-IC	Lateral Flow Assay for detection of Canine Distemper Virus in biological samples (tears, conjunctiva samples, nasal discharges, etc.).	12 & 30 individual devices	12 15.CDV.K.42/12 30 15.CDV.K.42/30
CANINE PARVOVIRUS DISEASE (CPV)			
INgezim Parvovirus IgM	Capture antibody ELISA for the detection of specific antibodies to CPV (IgMs) in canine serum samples.	1 plate kit	96 15.CPM.K.2
INgezim CPV DAS	Double antibody ELISA for the detection of CPV antigens in biological samples (dog's faeces).	1 plate kit	96 15.CPV.K.2
INgezim Parvo Canino	Indirect ELISA for the detection of specific antibodies to CPV in canine serum samples.	1 plate kit	96 15.CPV.K.1
CANINE PARVOVIRUS DISEASE (CPV)			
INgezim Parvo CROM	Lateral Flow Assay for the detection of CPV antigen in canine serum samples.	12 individual devices 50 individual devices	12 15.CP2.K.42/12 50 15.CP2.K.42/50
CANINE PARVOVIRUS DISEASE (CPV) & GIARDIA			
INgezim CPV/Giardia IC	Lateral Flow Assay for detection of canine parvovirus and Giardia lamblia antigen in biological samples (faeces and intestinal content).	12 individual devices	12 15.CGD.K.42/12
CANINE PARVOVIRUS DISEASE (CPV) & CORONAVIRUS			
INgezim CPV/CCV IC	Lateral Flow Assay for detection of Canine Parvovirus and Canine Coronavirus antigen in biological samples (faeces and intestinal content).	12 individual devices	12 15.CCD.K.42/12
CANINE PARVOVIRUS DISEASE (CPV), GIARDIA & CORONAVIRUS			
INgezim CPV/GIA/CCV IC	Lateral Flow Assay for the determination of parvovirus, canine coronavirus and Giardia lamblia in faeces samples.	12 individual devices	12 15.PGC.K.42/12
CANINE EHRLICHIOSIS			
INgezim Ehrlichia	Indirect ELISA for the detection of specific antibodies to Ehrlichia canis in canine serum samples.	1 plate kit 10 plates kit	96 15.EHR.K.1 960 15.EHR.K.1/10
INgezim Ehrlichia Vet	Indirect ELISA for the detection of specific antibodies to Ehrlichia canis in canine serum samples.	32 wells kit	32 15.EHR.K.8/32
CANINE EHRLICHIOSIS			
INgezim Ehrli CROM	Specific antibodies against to Ehrlichia canis	12 individual devices 30 individual devices	12 15.EHR.K.4/12 30 15.EHR.K.4/30

Product	Product description	No. of tests	Article No.
CANINE ERLICHIOSIS AND ANAPLASMOSIS			
INgezim Ehr/ANA IC	Lateral Flow Assay for the detection of anti-Ehrlichia canis and anti-Anaplasma antibodies.	12 individual devices	12 15.EHA.K.41/12
CANINE LEISHMANIOSIS			
INgezim Leishmania	Indirect ELISA for the detection of specific antibodies to Leishmania spp infantum in canine serum samples.	1 plate kit 10 plates kit	96 15.LSH.K.1 960 15.LSH.K.1/10
CANINE LEISHMANIOSIS			
INgezim Leishma CROM	Lateral Flow Assay for detection of antibodies to LSH in canine serum samples.	12 individual devices 50 individual devices	12 15.LSH.K.41/12 50 15.LSH.K.41/50
INgezim Leishmania Vet	Indirect ELISA for the detection of specific antibodies to Leishmania infantum in canine serum samples.	32 wells kit	32 15.LSH.K.8/32
RABIES			
E.I.A. Rabia	Double antibody enzymeimmunoassay for the quantification of rabies virus glycoprotein in infected cell culture samples.	2 plates kit	192 40.RAB.K.2/2
LYME DISEASE - BORRELIA BURGDORFERI			
INgezim Lyme CROM	Indirect Lateral Flow Assay for the detection of specific antibodies against to B. burgdogferi	12 individual devices	12 15.LY.K.4/12
DIROFILARIA IMMITIS - HEARTWORM			
INgezim Heart CROM	Indirect Lateral Flow Assay for the detection of dirofilaria immitis antigen	12 Individual devices 30 Individual devices	12 15.HW.K.4/12 30 15.HW.K.4/30
CANINE LEPTOSPIRA			
INgezim Leptospira IgM IC	Lateral flow test for the detection of IgM antibodies to Leptospira spp. in canine serum, plasma and whole blood.	12 individual devices	12 15.LEPTO.K.41/12
BRUCELLA CANIS			
INgezim Brucella IC	Lateral Flow Assay for the detection of antibodies to Brucella canis in canine serum, plasma and whole blood.	12 individual devices	12 15.BC.K.41/12

RABBIT specialities



Product	Product description	No. of tests	Article No.
MYXOMATOSIS			
INgezim Myxomatosis	Indirect ELISA for the detection of specific antibodies to Myxomatosis virus in rabbit serum samples.	2 plates kit	192 17.MIX.K.1/2
RABBIT HAEMORRHAGIC DISEASE (RHDV)			
INgezim RHDV	Indirect ELISA for the detection of specific antibodies to RHDV in rabbit serum samples.	2 plates kit	192 17.RHD.K.1/2
INgezim RHDV DAS	Double antibody ELISA for the detection of RHDV antigen in rabbit liver and other biological samples.	1 plate kit	96 17.RHD.K.2
INgezim® RHDV1/2 DIF CROM	Duplex double antibody sandwich lateral flow assay for the detection of rabbit haemorrhagic disease virus antigen, differentiating between RHDV and RHDV2 in rabbit liver and liver exudate samples.	12 individual devices 30 individual devices	12 17.RHD.K.42/12 30 17.RHD.K.42/30

FELINE specialities



Product	Product description	No. of tests		Article No.
FELINE CORONAVIRUS (FCV)				
INgezim Corona Felino	Indirect ELISA for the detection of specific antibodies to FCV in feline serum samples.	1 plate kit	96	16.FC.V.K.1
INgezim FCoV IC	Lateral Flow Assay for detection of specific antibodies to FCV in serum and blood samples.	12 individual devices	12	16.FC.V.K.41/12
FELINE IMMUNODEFICIENCY (FIV)				
INgezim FIV	Indirect ELISA for the detection of specific antibodies to FIV in feline serum samples.	1 plate kit	96	16.FIV.K.1
INgezim FIV Vet	Indirect ELISA for the detection of specific antibodies to FIV in feline serum samples.	32 wells kit	32	16.FIV.K.8/32
INgezim FIV Crom	Lateral Flow Assay for the detection of specific antibodies to FIV in feline serum and blood samples.	12 individual devices	12	16.FI2.K.41/12
		30 individual devices	30	16.FI2.K.41/30
FELINE LEUKAEMIA (FeLV)				
INgezim FeLV	Indirect ELISA for the detection of specific antibodies to FeLV gp70 in feline serum samples.	1 plate kit	96	16.FLV.K.1
INgezim FeLV CROM	Lateral Flow Assay for detection of FeLV (p27) antigen in feline serum samples.	12 individual devices	12	16.FLV.K.42/12
		30 individual devices	30	16.FLV.K.42/30
INgezim FeLV DAS	Double antibody ELISA for the detection of FeLV (p27) antigen in feline serum samples.	1 plate kit	96	16.FLV.K.2
INgezim FeLV Vet	Double antibody ELISA for the detection of FeLV (p27) antigen in feline serum samples.	32 wells kit	32	16.FLV.K.9/32
INgezim FeLV/FIV IC	Lateral Flow Assay for detection of FeLV antigen and specific antibodies to FIV in feline blood, serum or plasma samples.	12 individual devices	12	16.FFD.K.40/12
FELINE PANLEUCOPENIA				
INgezim FPV IC	Lateral Flow Assay for detection of Feline Panleucopenia Virus antigen in biological samples (feces and intestinal content)	12 individual devices	12	16.FPV.K.42/12
FELINE TOXOPLASMA				
INgezim Feline toxoplasma IC	Lateral Flow Assay for the detection of IgG and IgM antibodies to Toxoplasma gondii in feline serum, plasma and whole blood.	12 individual devices	12	16.TOX.K.41/12

EQUINE specialities



Product	Product description	No. of tests		Article No.
AFRICAN HORSE SICKNESS (AHSV)				
INgezim AHSV Compac Plus	Blocking ELISA for the detection of specific antibodies to AHSV in equine serum samples.	2 plates kit	192	14.AHS.K.3/2
		5 plates kit	480	14.AHS.K.3/5
INgezim PEA DAS	Double antibody ELISA for the detection of specific AHSV antigen in equine biological samples (spleen and tissues).	1 plate kit	96	14.PEA.K.2
INgezim AHSV CROM Anticuerpo	Lateral Flow Assay for the detection of specific antibodies to AHSV in equine serum samples.	30 individual devices	30	14.AHS.K.41/30
EQUINE INFECTIOUS ANEMIA (EIAV)				
INgezim Anemia DR	Double recognition ELISA for the detection of specific antibodies to EIAV in equine serum samples.	2 plates kit	192	14.AIE.K.0/2
INgezim Anemia CROM	Lateral Flow Assay for the detection of specific antibodies to EIAV in equine serum samples.	30 individual devices	30	14.AIE.K.41/30
EQUINE ARTERITIS (EAV)				
INgezim Arteritis 2.0	Indirect ELISA for the detection of specific antibodies to EAV in equine serum samples.	2+2 plates kit	192	14.EA2.K.1/4
EQUINE RHINOPNEUMONITIS				
INgezim Rinoneumonitis	ELISA for the detection of specific antibodies to Equine herpesvirus (EHV) subtypes 1 & 4 in equine serum samples	2 plates kit	192	14.HVE.K.1/2
		5 plates kit	480	14.HVE.K.1/5
WEST NILE FEVER (WNV)				
INgezim West Nile IgM	Capture antibody ELISA for the detection of specific IgMs to WNV in equine serum samples.	1+1 plates kit	96	14.WNV.K.2

VETLINE ASSAYS

VetLine is the veterinary product line of our strategic partner, Novatec Immundiagnostica (part of Gold Standard Diagnostics Europe), providing a wide range of ELISA kits for antibody detection across a range of species including canine, feline, bovine, porcine, caprine, and poultry.

- Ready-to-use and colour coded reagents
- Standardized procedures
- ELISAs with cut off control
- Easy to use on automated systems
- Snap-off microtiter strip
- Long shelf-life



CANINE specialities

Product	Product description	No. of tests		Article No.
VetLine Anaplasma	Indirect ELISA for the detection of antibodies to Anaplasma in canine serum and plasma samples.	1 plate kit	96	15.ANAV0850
VetLine Borrelia	Indirect ELISA for the detection of antibodies to Borrelia in canine serum and plasma samples.	1 plate kit	96	15.BORVT0040
VetLine Canine Distemper Virus (CDV)	Indirect ELISA for the detection of antibodies to CDV in canine serum and plasma samples.	1 plate kit	96	15.MEAVT0330
VetLine Canine Parvovirus (CPV)	Indirect ELISA for the detection of antibodies to CPV in canine serum and plasma samples.	1 plate kit	96	15.PARVT0370
VetLine Dirofilaria Antigen	Double antibody ELISA for the detection of Dirofilaria antigen in canine serum and plasma samples.	1 plate kit	96	15.DIRVT4760
VetLine Ehrlichia	Indirect ELISA for the detection of antibodies to Ehrlichia in canine serum and plasma samples.	1 plate kit	96	15.EHRVT0930
VetLine Leishmania	Indirect ELISA for the detection of antibodies to Leishmania in canine serum and plasma samples.	1 plate kit	96	15.LEIVT0310
VetLine Leptospira	Indirect ELISA for the detection of antibodies to Leptospira in canine serum and plasma samples.	1 plate kit	96	15.LEPVT0660
VetLine Babesia*	Indirect ELISA for the detection of antibodies to Babesia in canine serum and plasma samples.	1 plate kit	96	15.BABVT0890
VetLine Canine C-Reactive* Protein (CCRP) - coming soon	Indirect ELISA for the detection of antibodies to canine C-Reactive protein in canine serum and plasma samples.	1 plate kit	96	15.CCRPV4050
VetLine Echinococcus*	Indirect ELISA for the detection of antibodies to Echinococcus in canine serum and plasma samples.	1 plate kit	96	15.ECHVT0130
VetLine Leptospira IgM*	Indirect ELISA for the detection of IgM antibodies to Leptospira in canine serum and plasma samples.	1 plate kit	96	15.LEPVM0660
VetLine Toxocara*	Indirect ELISA for the detection of antibodies to Toxocara in canine serum and plasma samples.	1 plate kit	96	15.TOCVT0450

FELINE specialities



Product	Product description	No. of tests		Article No.
VetLine Feline Immunodeficiency Virus (FIV)	Indirect ELISA for the detection of antibodies to FIV in feline serum and plasma samples.	1 plate kit	96	16.FIVVT0750
VetLine Feline Leukemia Virus (FeLV) Antigen	Double antibody ELISA for the detection of FeLV antigen in feline serum and plasma samples.	1 plate kit	96	16.FELVT4800
VetLine Bartonella*	Indirect ELISA for the detection of antibodies to Bartonella in feline serum and plasma samples.	1 plate kit	96	16.BARVT0900

MULTISPECIES specialities



Product	Product description	No. of tests		Article No.
VetLine Brucella	Indirect ELISA for the detection of antibodies to Brucella in bovine and porcine serum, pooled serum and milk samples.	1 plate kit	96	10.BRUVT0050
VetLine TBE/FSME	Indirect ELISA for the detection of antibodies to TBE/FSME in Multispecies serum and plasma samples.	1 plate kit	96	10.TICVT0440
VetLine Toxoplasma	Indirect ELISA for the detection of antibodies to Toxoplasma in canine/feline serum and plasma samples.	1 plate kit	96	10.TOXVT0460
VetLine Toxoplasma IgM	Indirect ELISA for the detection of IgM antibodies to Toxoplasma in canine/feline serum and plasma samples.	1 plate kit	96	10.TOXVM0460
VetLine Giardia Antigen*	Double antibody ELISA for the detection of Giardia antigen in stool samples.	1 plate kit	96	10.GIAVT4160
VetLine TBE/FSME IgM*	Indirect ELISA for the detection of IgM antibodies to TBE/FSME in Multispecies serum and plasma samples.	1 plate kit	96	10.TICVM0440
VetLine Sarcoptes	Indirect ELISA for the detection of antibodies to Sarcoptes scabiei in serum samples.	1 plate kit	96	10.SARVT0980

PORCINE specialities



Product	Product description	No. of tests		Article No.
VetLine Trichinella	Indirect ELISA for the detection of antibodies to Trichinella in porcine serum and plasma samples.	1 plate kit	96	11.TRIVT0480
VetLine Pasteurella multocida Toxin Antibody *	Indirect ELISA for the detection of antibodies to Pasteurella multocida Toxin in porcine serum and plasma samples.	1 plate kit	96	11.PASVT0960

AVIAN specialities



Product	Product description	No. of tests		Article No.
VetLine Riemerella	Indirect ELISA for the detection of antibodies to Riemerella in poultry serum and plasma samples.	1 plate kit	96	18.RIEMVT0880

RUMINANT specialities



Product	Product description	No. of tests		Article No.
VetLine Bovine Haptoglobin*	Indirect ELISA for the detection of antibodies to Bovine Haptoglobin in bovine serum and plasma samples.	1 plate kit	96	12.BHAPV4040
VetLine Coxiella	Indirect ELISA for the detection of antibodies to Coxiella in plasma and milk.	1 plate kit	96	12.COXVT0600
VetLine Coxiella Phase 1*	Indirect ELISA for the detection of antibodies to Coxiella Phase 1 in serum, plasma and milk.	1 plate kit	96	12.COX1VT0600
VetLine Coxiella Phase 2*	Indirect ELISA for the detection of antibodies to Coxiella Phase 2 in serum plasma and milk.	1 plate kit	96	12.COX2VT0600

*For Research Use Only



RapidScan ST5-W Lateral Flow Assay Reader

Companies and laboratories that use lateral flow devices as a rapid test method need to obtain quick and reliable results for contamination testing.

The RapidScan ST5-W Lateral Flow Assay Reader is an accurate, easy to use and highly flexible imaging platform used for lateral flow assay analysis. The instrument was designed specifically for field and in-process testing applications that require qualitative and quantitative test results.



Key Benefits

- **Ease-of-use:** the reader is easy to set up and operate, featuring an intuitive user interface with a colour touchscreen display.
- **Quantitative results** are obtained within a few minutes. Data can be stored on the reader then exported through a USB port.
- **Applications:** the INgezim ASF CROM Ag, Lateral Flow Assay for the detection of ASFV antigen in porcine blood samples, is validated on the RapidScan reader.

RapidScan ST5-W

Specifications

Product Name	RapidScan ST5-W Lateral Flow Assay Reader
Article No.	LFRSCAN002
Supplied Materials	RapidScan ST5-W Reader Power Cable and International Adapters QR Code Adapter Check Chart User Guide
Item Dimensions	11 x 14 x 21cm
Item Weight	930 g
Display	5" LCD touch screen
Platform	Lateral Flow
Power	DC 12V 1.5A
Interface	USB 2.0 port
Application Adapters	Reader price includes one application adapter of choice.



ELISA Analyzers

Comprehensive diagnostic solutions improve lab efficiency and minimize overall costs. We offer precise, compact and cost-efficient instruments supported with outstanding customer service.

These innovative, industry-leading instruments - ThunderBolt® and The BOLT™ - are the standard platforms for Eurofins Technologies ELISA-based analytical test kits including allergens, food pathogens, mycotoxins, veterinary drug residues, and animal health kits. We are delighted to provide one of the most complete, integrated, and innovative test offerings in the market.

ThunderBolt®



The ThunderBolt® is an innovative 2-plate, fully automated, open platform that is compact and cost efficient. The user friendly and flexible software is capable of programming a broad range of protocols. The ThunderBolt® is available for both ELISA and Chemiluminescent (CLIA) applications.

Key Benefits

- Fully-automated: load and walk away
- High capacity: 192 samples
- Open architecture: program any ELISA or Chemiluminescent (CLIA) protocol* (CLIA protocols require ELISA+CLIA reader option)
- Adaptable: capable of running many different protocols in a single batch
- Built-in reader: ultra-compact and fully-automated reader options
- Compact: smallest and lightest instrument in its class
- User-friendly and flexible software
- Extensive configuration – run up to 8 protocols in a single batch
- On-board camera: monitor internal operations in real-time

The Bolt™



The Bolt™, a one-plate ELISA+CLIA processor, shares many features of the ThunderBolt®, with a more customizable, lower throughput and cost-effective design. Our interactive user-friendly software provides a simple laboratory experience for instrument operation.

Key Benefits

- Space saving design – fits standard 60 cm lab bench
- 96 sample capacity
- Program any EIA and CLIA protocol
- High precision syringes
- Low consumables costs (no disposable tips)
- Exterior status indicator light
- Connectivity to a laboratory information system (LIS)

Modular Design with Optional:

- Convention incubator
- Linear shaker
- ELISA or ELISA+ CLIA Reader*

* CLIA protocols require ELISA + CLIA reader option



ThunderBolt®

Specifications	
Intelligent Racks	Automatic sample location management
System Architecture	Open, fully customizable
High Precision Micro-Syringe	Aspirate 1uL with ≤ 3% CV
On-board Reader(s)	Spectrophotometer or Spectrophotometer + Fluorescence Combo Reader
Forced Convection Incubator	Evenly heats using forced air to eliminate "Edge Effect"
On-board Camera	Probe-mounted camera allows remote troubleshooting
Sliding Sample Tray	Easy loading and unloading of samples
Built-in Barcode Reader	Streamline loading and decrease sample tracking errors
Orbital MTP Shaker	Precise adherence to assay protocols
Small Footprint	Fits on a 60cm deep standard laboratory bench; weight: 30 kg (65 lbs)
Sample Capacity	192 patient sample positions
User Interface	MS Windows Graphical User Interface (Windows 7 or above)
CE Marked	Yes
Article No.	00300



The Bolt™

Specifications	
System Architecture	Open, fully customizable
High Precision Micro-Syringe	Aspirate 1uL with ≤ 3% CV
On-board Reader(s)	Spectrophotometer or Spectrophotometer + Fluorescence Combo Reader
Forced Convection Incubator	Evenly heats to eliminate "Edge Effect"
Linear Shaker	Precise adherence to assay protocols / No-spill design
Small Footprint	Fits on a 60cm deep standard laboratory bench; weight: 27 kg (59.5 lbs)
Sample Capacity	96 patient sample positions
User Interface	MS Windows Graphical User Interface (Windows 7 or above)
CE Marked	Yes
Article No.	00500



