

*SENSISpec Spike Solution Gluten**Cat.No. M.30.GLU.A04003m*

<b>Product:</b>	SENSISpec Spike Solution Gluten
<b>Catalog Number:</b>	<b>M.30.GLU.A04003m</b>
<b>Lot:</b>	09-2021
<b>Quantity:</b>	3 mL, 25 ppm Gluten, Uncertainty $\pm$ 20%
<b>Physical state:</b>	liquid
<b>Storage:</b>	2-8°C
<b>Expiry:</b>	2023.03.

**Intended Use:**

The SENSISpec INgezim Spike Solution Gluten is intended for the execution of spiking experiments in combination with the SENSISpec INgezim GLUTEN R5 ELISA kits (30.GLU.K2, 30.GL2.K2 and 30.GLH.K2)

**Application example:**

In order to get a sample containing 10 ppm of GLUTEN, add 100  $\mu$ L of SENSISpec Spike Solution Gluten to 0.25 g of sample. Extract the sample according to the SENSISpec INgezim GLUTEN R5 ELISA kit instruction. After centrifugation, the supernatant must be diluted 1:25 with the prediluted diluent included in the kit.

Date: 2022.01.17



---

Control Manager