



# Microgen

## Legionella Latex test

A rapid latex agglutination test providing results in two minutes directly from a selective medium for the confirmatory identification of:

- *Legionella pneumophila* Serogroup 1 and Serogroup 2-15
- Non-*L. pneumophila* species of clinical significance grown on selective media



### Fast & Easy

- Results in 2 minutes
- No equipment or highly trained personnel required



### Inclusive

- Detects *Legionella pneumophila* Serogroup 1 and 2-15, together with 10 of the most clinically relevant species



### Efficient

- Quickly integrates into laboratory workflow improving efficiency
- 100% specificity and sensitivity compared to traditional culture methods

Microgen *Legionella pneumophila* identifies a range of non-*L. pneumophila* species of clinical significance with high specificity and no cross-reaction with *L. pneumophila*. The latexes used in Microgen *Legionella* are coated with immunoglobulins which have been specially processed to minimize issues with potentially cross reacting organisms thereby eliminating the need for a control latex.

The test is **compliant with ISO 11731 Water Quality** – Detection and enumeration of *Legionella*.

Detects the following species:

- *L. pneumophila* Serogroup 1 and 2-15
- *L. micdadei*
- *L. bozemanii* I & II
- *L. dumoffii*
- *L. longbeachae* I & II
- *L. jordanis*
- *L. gormanii*
- *L. anisa*
- *L. feelei*



### Ordering Information

Code	Description	Size	Storage
M45	Microgen <i>Legionella</i> latex test	50 tests	2 - 8°C



**Gold Standard Diagnostics**  
 contact@eu.goldstandarddiagnostics.com • www.goldstandarddiagnostics.com