

SENSISpec Allergen Spike Solutions

Our SENSISpec Allergen Spiking Solutions can be used to produce positive control samples suitable for validation of Allergen ELISA kits, such as determination of recovery rate in a spike/recovery assay and matrix interference in a linearity assay. The procedures for spiking samples are available upon request.

Article No.	Product Name	Quantity	Concentration	Shelf life
HU0030049	SENSISpec Spike Solution Almond	3 mL	1000 ppm	24
HU0030050	SENSISpec Spike Solution Beta-Lactoglobulin	3 mL	10 ppm	13
HU0030051	SENSISpec Spike Solution Casein	3 mL	1000 ppm	13
HU0030052	SENSISpec Spike Solution Cashew	3 mL	1000 ppm	13
HU0030054	SENSISpec Spike Solution Coconut	3 mL	1000 ppm	13
HU0030055	SENSISpec Spike Solution Tropomyosin (Crustacean)	3 mL	10 ppm	24
HU0030056	SENSISpec Spike Solution Egg White	3 mL	1000 ppm	24
HU0030057	SENSISpec Spike Solution Fish (Cod)	3 mL	1000 ppm	24
HU0030058	SENSISpec Spike Solution Gliadin	3 mL	1000 ppm	18
HU0030059	SENSISpec Spike Solution Hazelnut	3 mL	1000 ppm	18
HU0030060	SENSISpec Spike Solution Lupin	3 mL	1000 ppm	13
HU0030061	SENSISpec Spike Solution Lysozyme	3 mL	10 ppm	24
HU0030062	SENSISpec Spike Solution Macadamia nut	3 mL	1000 ppm	13
HU0030063	SENSISpec Spike Solution Milk	3 mL	1000 ppm	13
HU0030064	SENSISpec Spike Solution Tropomyosin (Molluscs)	3 mL	10 ppm	13
HU0030065	SENSISpec Spike Solution Mustard	3 mL	1000 ppm	24
HU0030066	SENSISpec Spike Solution Ovalbumin	3 mL	10 ppm	24
HU0030067	SENSISpec Spike Solution Brazil nut	3 mL	1000 ppm	13
HU0030068	SENSISpec Spike Solution Peanut	3 mL	1000 ppm	13
HU0030069	SENSISpec Spike Solution Pecan nut	3 mL	1000 ppm	13
HU0030070	SENSISpec Spike Solution Pistachio	3 mL	1000 ppm	13
HU0030071	SENSISpec Spike Solution Sesame	3 mL	1000 ppm	18
HU0030072	SENSISpec Spike Solution Soy Trypsin Inhibitor	3 mL	10 ppm	24
HU0030073	SENSISpec Spike Solution Walnut	3 mL	1000 ppm	13

